




## 8359 Genu Direxa Stable

## 8368 Genu Direxa Stable wraparound

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	9
<b>FR</b> Instructions d'utilisation .....	14
<b>IT</b> Istruzioni per l'uso .....	19
<b>ES</b> Instrucciones de uso .....	24
<b>NL</b> Gebruiksaanwijzing .....	29
<b>SV</b> Bruksanvisning .....	34
<b>DA</b> Brugsanvisning .....	39
<b>NO</b> Bruksanvisning .....	43
<b>FI</b> Käyttöohje .....	48
<b>PL</b> Instrukcja użytkowania .....	53
<b>SK</b> Návod na používanie .....	57
<b>BG</b> Инструкция за употреба .....	62
<b>TR</b> Kullanma talimatı .....	67
<b>ET</b> Kasutusjuhhis .....	72
<b>RU</b> Руководство по применению .....	77
<b>ZH</b> 使用说明书 .....	82

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#### INFORMATION

Last update: 2017-07-05

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.

- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 8359 Genu Direxa Stable and 8368 Genu Direxa Stable wraparound knee joint braces.

## 2 Components/design

Item	Quantity	Designation	Article number	Figures
	1	Genu Direxa Stable, Genu Direxa Stable wraparound	8359, 8368	
		<b>Consisting of:</b>		
1		Donning loops		(see fig. 2, see fig. 7)
2		Upper wide hook-and-loop closure, front		
3		Upper narrow hook-and-loop closure, front		
4		Patella cutout		
5		Lower narrow hook-and-loop closure, front		
6		Lower wide hook-and-loop closure, front		
7		Lower wide hook-and-loop closure, back		
8		Lower narrow hook-and-loop closure, back		
9		Upper narrow hook-and-loop closure, back		
10		Upper wide hook-and-loop closure, back		
		Joint bar		(see fig. 9)
		Extension stop		(see fig. 10)

## 3 Intended use

### 3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 3.2 Indications

Knee complaints with severe and/or complex ligament laxities and/or moderate to severe feelings of instability where temporary or permanent extension and/or flexion limitation of the knee joint is indicated, for example:

- Knee pain with pre-operative or chronic instability, functional conservative treatment of knee ligament injuries and deformities
  - Stretched/ruptured ACL (anterior cruciate ligament)
  - Stretched/ruptured PCL (posterior cruciate ligament)
  - Stretched/ruptured collateral ligaments
  - Genu recurvatum
- Severe instability due to patellofemoral pain syndrome
  - Chondromalacia
  - Chondropathy patellae
  - Osteoarthritis of the patellofemoral joint
- Stabilization of severe and/or complex post-traumatic and post-operative instabilities
  - Ligament and/or meniscus reconstruction
  - Use during rehabilitative measures
- Severe and/or complex degenerative instabilities
  - Osteoarthritis of the knee

- Rheumatoid arthritis

The indication must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

None known.

#### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

### 3.4 Mechanism of Action

The brace improves neuromuscular function, proprioception and circulation. It protects against hyperextension, limits the adjustable range of motion in extension and flexion, supports sensorimotor functions and circulation, promotes the resorption of edemas and can relieve pain.

## 4 Safety

### 4.1 Explanation of warning symbols

 Warning regarding possible risks of accident or injury.

 Warning regarding possible technical damage.

### 4.2 General safety instructions

#### Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

#### Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

#### Contact with ointments, lotions, or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to ointments, lotions, or other products that contain oils or acids.

## 5 Handling

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 5.1 Selecting the Size

- 1) Measure the circumference of the thigh 15 cm above the center of the patella.
- 2) Measure the circumference of the lower leg 15 cm below the center of the patella.
- 3) Determine the size of the brace (see size chart).

### 5.2 Fitting and Application

#### NOTICE

#### Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

#### 5.2.1 Adaptation

The joint bars are fitted individually by trained qualified personnel.

### INFORMATION

The joint stops (medial and lateral) must have identical degree settings.

Available extension stops (square grip):	0°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

- 1) Pull the black strap upwards to remove the bracket over the joint bar (see fig. 8).
- 2) Pull the joint bar out of the brace (see fig. 9).
- 3) Insert the desired extension and flexion stops (see fig. 10).  
→ The stops will engage with slight resistance.
- 4) Move the brace joints and check the position of the stops.
- 5) Slide the joint bar into the brace (see fig. 9).
- 6) Pull the black strap first upwards and then downwards to position the bracket over the joint bar (see fig. 8).

#### 5.2.2 Application

### CAUTION

#### Use with heat-sensitive skin

Skin irritation due to overheating

- ▶ Do not use the product in case of a known heat allergy.
- ▶ Do not continue to use the product if skin irritation occurs.

## **CAUTION**

### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- Ensure that the brace is applied properly and fits correctly.

## **INFORMATION**

After correct adjustment, leave the rear hook-and-loop closures fastened.

### **8359 Genu Direxa Stable**

- > Open all hook-and-loop closures on the brace.
  - > The patient sits on the front edge of a chair.
- 1) Hold the brace on both sides with the donning loops (item 1) and pull it over the knee until the patella cutout (item 4) is centered on the patella (see fig. 1 and see fig. 2).
  - 2) Fasten the wide hook-and-loop closure at the top (see fig. 2, item 2).
  - 3) Fasten the wide hook-and-loop closure at the bottom (see fig. 2, item 6).
  - 4) The two wide hook-and-loop closures at the rear may only be adjusted and fastened by qualified personnel (see fig. 2, item 7 and item 10).
  - 5) Fasten the lower front hook-and-loop closure (see fig. 2, item 5).
  - 6) Fasten the upper front hook-and-loop closure (see fig. 2, item 3).
  - 7) The two hook-and-loop closures at the rear may only be adjusted and fastened by qualified personnel (see fig. 2, item 8 and item 9, and see fig. 3). **Optional:** For the corresponding indication, fasten the hook-and-loop closures crosswise (see fig. 4).

### **8368 Genu Direxa Stable wraparound**

- > Open all hook-and-loop closures on the brace.
  - > The patient sits on the front edge of a chair.
- 1) Wrap the brace around the knee from behind (see fig. 5).
  - 2) Center the patella cutout over the patella (see fig. 6).
  - 3) Fasten the hook-and-loop closure below the patella cutout while applying tension. Make sure that the hook-and-loop sections are placed over one another (see fig. 6).
  - 4) Fasten the hook-and-loop closure above the patella cutout while applying tension. Make sure that the hook-and-loop sections are placed over one another (see fig. 6).
  - 5) Fasten the wide hook-and-loop closure at the top (see fig. 7, item 2).
  - 6) Fasten the wide hook-and-loop closure at the bottom (see fig. 7, item 6).
  - 7) The two wide hook-and-loop closures at the rear may only be adjusted and fastened by qualified personnel (see fig. 7, item 7 and item 10).
  - 8) Fasten the lower front hook-and-loop closure (see fig. 7, item 5).
  - 9) Fasten the upper front hook-and-loop closure (see fig. 7, item 3).
  - 10) The two hook-and-loop closures at the rear may only be adjusted and fastened by qualified personnel (see fig. 7, item 8 and item 9, and see fig. 3). **Optional:** For the corresponding indication, fasten the hook-and-loop closures crosswise (see fig. 4).

## **5.3 Cleaning**

### **NOTICE**

#### **Use of improper cleaning agents**

Damage to the product due to use of improper cleaning agents

- Only clean the product with the approved cleaning agents.

- 1) Pull the black strap upwards to remove the bracket over the joint bar (see fig. 8).
- 2) Pull the joint bars out of the brace (see fig. 9).

- 3) Fasten all hook-and-loop closures.
- 4) Recommendation: Use a laundry bag or net.
- 5) Wash the brace in warm water at **40 °C (104 °F)** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 6) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 7) Reinsert the joint bars.  
**INFORMATION: The labelling on the joint bars must be visible from the outside.**
- 8) Pull the black strap first upwards and then downwards to position the bracket over the joint bar.

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

### 7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

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Otto Bock HealthCare GmbH