
17K29*, 17K32*, 17K33*


DE	Gebrauchsanweisung	5
EN	Instructions for use	11
FR	Instructions d'utilisation	17
IT	Istruzioni per l'uso	23
ES	Instrucciones de uso	29
PT	Manual de utilização	36
NL	Gebruiksaanwijzing	42
SV	Bruksanvisning	48
DA	Brugsanvisning	54
NO	Bruksanvisning	60
PL	Instrukcja użytkowania	66
HU	Használati utasítás	72
CS	Návod k použití	78
TR	Kullanma talimatı	84
EL	Οδηγίες χρήσης	90
RU	Руководство по применению	97

1

17K29*


2

17K32*

17K33*

3

INFORMATION

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INFORMATION

Date of last update: 2021-05-11

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17K29*, 17K32* and 17K33* orthotic knee joints.

2 Product description

2.1 Available sizes

Knee joint bar			
Reference number	Size	Width in mm	Material
17K29*, 17K32*, 17K33*	6, 5, 4	12, 14, 16	Stainless steel

2.2 Components/design

Scope of delivery for 17K29* (see fig. 1)			
Item	Quantity [pcs.]	Designation	Article number
1	1	Upper bar	
2	1	Joint screw	501S32=M4x10x9.5 (size 6, 5), 501S32=M6 (size 4)
3	1	Ring lock	

Scope of delivery for 17K29* (see fig. 1)			
Item	Quantity [pcs.]	Designation	Article number
			17Y13=5 (size 6, 5), 17Y13=4 (size 4)
4	1	Lower bar	
5	1	Clamping sleeve	506S1=3x14
6	1	Bearing nut	17Y93=6x6.7xM4 (size 6, 5), 17Y93=9x7.2xM6 (size 4)

Not included in the scope of delivery (see fig. 1)			
Item	Quantity [pcs.]	Designation	Article number
7	1	Spring-loaded thrust piece	501D1

Scope of delivery for 17K32*, 17K33* (see fig. 2)			
Item	Quantity [pcs.]	Designation	Article number
1	1	Upper bar	
2	1	Joint screw	501S32=M4x10x9.5 (size 6, 5), 501S32=M6 (size 4)
3	1	Lower bar	
4	1	Bearing nut	17Y93=6x6.7xM4 (size 6, 5), 17Y93=9x7.2xM6 (size 4)

3 Intended use

3.1 Indications for use

The orthotic knee joints are to be used **exclusively** in pairs to provide orthotic treatment for a patient's lower limbs.

3.2 Indications

- Partial or total paralysis of the leg muscles
 - Orthopaedic diseases of the lower limbs
- Indications must be determined by the physician.

3.3 Lifetime

Given intended use and proper installation, the product is designed for a lifetime of **1 year**.

3.4 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

 **CAUTION**

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

 **CAUTION**

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **CAUTION**

Improper use of the lock

Risk of injury due to unexpected deactivation of the orthotic joint lock

- ▶ Fully extend the orthotic joint to activate the lock. Check whether the lock is engaged.
- ▶ Only deactivate the lock when you are able to control the flexion of the orthotic joint.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

5 Preparing the product for use

⚠ CAUTION

Repeated bending at the same position

Risk of injury due to breakage of the bar, change in or loss of functionality

- ▶ Avoid repeated bending at the same position.

NOTICE

Orthotic joints not aligned in parallel

Malfunction of the locking mechanism

- ▶ Use the 743R6 orthotic joint alignment fixture.
- ▶ Make sure the orthotic joints lock simultaneously in the extension stop.

NOTICE

Shaping under the effect of heat

Damage due to material changes

- ▶ Do not apply heat to shape the product.

INFORMATION

Increased wear and tear

The orthotic knee joints are intended only for use in combination with free-motion or foot lifting orthotic ankle joints. Wear occurs more rapidly in combinations involving orthotic ankle joints with a dorsal stop.

Shaping the bars

- > **Required tools:** 711S5 bending irons, 711S3 horizontal bending irons, 711S2 bending pusher set.
- 1) Use the 711S5 bending irons for bending.
- 2) Use the 711S3 horizontal bending irons or the 711S2 bending pusher set for horizontal bending.
- 3) Use a bending radius over 10 mm (see fig. 3).
- 4) Remove grooves and burrs by grinding.

Using protective devices

- ▶ Use joint protectors to protect the product against soiling.

Improving corrosion protection

Polish or sinter the product for improved corrosion protection. Ottobock recommends the 618T40* sintering powder for surface coating.

- ▶ Do not heat the material longer than **5 minutes** at max. **150 °C** for surface coating.

6 Cleaning

After contact with water containing salt, chlorine or soap, or if they get dirty, the orthotic joints must be promptly cleaned.

- 1) Rinse the orthotic joints with clean fresh water.

- 2) Dry the orthotic joints with a cloth or allow them to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Spare parts are listed in the section "Components/design".

Only use 633F7 special lubricant.

17K29: Replacing the ring lock

- 1) Loosen the joint screw and remove the bearing nut.
- 2) Disassemble the orthotic joint.
- 3) Remove the clamping sleeve and disassemble the ring lock.
- 4) Replace the ring lock and assemble the clamping sleeve.
- 5) Assemble the orthotic joint with the screw and the bearing nut.
- 6) For final assembly, secure the joint screw with Loctite 241.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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