1C50 Taleo, 1C53 Taleo LP


DE	Gebrauchsanweisung	4
EN	Instructions for use	13
FR	Instructions d'utilisation	22
IT	Istruzioni per l'uso	32
ES	Instrucciones de uso	42
PT	Manual de utilização	52
NL	Gebruiksaanwijzing	62
SV	Bruksanvisning	71
DA	Brugsanvisning	80
NO	Bruksanvisning	89
FI	Käyttöohje	97
PL	Instrukcja użytkowania	107
HU	Használati utasítás	117
CS	Návod k použití	126
RO	Instrucțiuni de utilizare	135
HR	Upute za uporabu	145
SL	Navodila za uporabo	154
SK	Návod na používanie	163
BG	Инструкция за употреба	172
TR	Kullanma talimatı	182
EL	Οδηγίες χρήσης	191
RU	Руководство по применению	201
JA	取扱説明書	211
ZH	使用说明书	220
KO	사용 설명서	227

1

2

3

4

INFORMATION

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INFORMATION

Last update: 2019-09-19

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 1C50 Taleo and 1C53 Taleo LP prosthetic feet are suitable for walking on various surfaces and for a wide range of speeds. The heel strike is absorbed by the large, replaceable heel wedge.

Carbon and polymer spring elements permit perceptible plantar flexion at heel strike, a natural rollover movement and high energy return.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

The table that follows shows the suitable spring stiffness of the prosthetic foot, matching the patient's body weight and activity.

Spring stiffness relative to body weight and activity level		
Body weight [kg]	Normal activity level	High activity level
up to 51	1	2
52 to 58	2	3
59 to 67	3	4
68 to 77	4	5
78 to 88	5	6
89 to 100	6	7
101 to 115	7	8
116 to 130	8	9
131 to 150	9	–

2.3 Environmental conditions

Allowable environmental conditions
Temperature range: –10 °C to +45 °C (14 °F to 113 °F)
Chemicals/liquids: fresh water, salt water, perspiration, urine, soapsuds, chlorine water
Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand

Allowable environmental conditions

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see Page 20).

Unallowable environmental conditions

Solids: Highly hygroscopic particles (e.g. talcum), continuous contact with sand

Chemicals/liquids: Acids, continuous use in liquid media

Storage and transport

Temperature range $-20\text{ }^{\circ}\text{C}$ to $+60\text{ }^{\circ}\text{C}$ ($-4\text{ }^{\circ}\text{F}$ to $+140\text{ }^{\circ}\text{F}$), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

2.4 Service life

Prosthetic foot

Depending on the patient's level of activity, the service life of the product is 2 to 3 years.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **CAUTION**

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 14).

⚠ CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.

- Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the carbon spring are indications of loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Protective sock	SL=Spectra-Sock-7
1	Heel wedge kit	2F50*

Spare parts/accessories (not included in the scope of delivery)

Designation	Reference number
Footshell	2C15*
Connection cover	2C19*, 2C20*
Bolt cover (1C50 only)	2F51*

5 Preparation for use

CAUTION

Incorrect alignment, assembly or adjustment

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

NOTICE

Grinding the prosthetic foot or footshell

Premature wear resulting from damage to the product

- Do not grind the prosthetic foot or footshell.

5.1 Alignment

INFORMATION

There is a plastic adapter cover on the proximal connection of the prosthetic foot. This protects the connecting section from scratches during the alignment and trial fitting.

- ▶ Remove the adapter cover before the patient leaves the workshop/fitting area.

5.1.1 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.
 - ▶ Always use the prosthetic foot with the footshell.
- ▶ Apply or remove the footshell as described in the footshell instructions for use.

5.1.2 Bench Alignment

INFORMATION

Only use the prosthetic foot with the heel wedge.

TT bench alignment

Bench alignment process

Required materials: 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)

Assemble and align the prosthetic components in the alignment apparatus according to the following specifications:

①	Heel height: effective heel height (shoe heel height - sole thickness of forefoot) + 5 mm
②	Sagittal alignment: alignment reference line on a-p markings of the footshell (see fig. 2)
③	Frontal alignment: alignment reference line on m-l markings of the footshell (see fig. 2)
④	Connect the prosthetic foot and prosthetic socket using the selected adapters. Follow the adapter instructions for use.
⑤	Sagittal alignment: Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: individual residual limb flexion + 5°
⑥	Observe the abduction or adduction position.

TF bench alignment

- ▶ Observe the information in the prosthetic knee joint instructions for use.

5.1.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.1.4 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- Remove the adapter cover from the prosthetic foot after completion of the dynamic trial fitting and gait training exercises.

5.1.4.1 Optimising the heel characteristics

The behaviour of the prosthetic foot at heel strike and during heel contact in the mid-stance phase can be adapted by replacing the heel wedge. The prosthetic foot is delivered with 3 heel wedges.

Hardness of the heel wedges: The heel wedges are colour-coded: (transparent=soft, grey=medium, black=hard). Ottobock recommends starting with the grey heel wedge.

- 1) Pull the prosthetic foot apart slightly and remove the existing heel wedge.
- 2) Align the other heel wedge so the Ottobock lettering is the right way up and the tip points to anterior.
- 3) Insert the heel wedge into the prosthetic foot (see fig. 3).

5.2 Optional: Installing the foam cover

The foam cover sits between the prosthetic socket and prosthetic foot. It is cut longer in order to compensate for the movements of the prosthetic foot and prosthetic knee joint. During flexion of the prosthetic knee joint, the foam cover undergoes posterior compression and anterior elongation. The foam

cover should be stretched as little as possible in order to increase its service life. There is a connecting element (such as a connection plate, connection cap or connection cover) on the prosthetic foot.

> **Required materials:** degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive

1) Measure the length of the foam cover on the prosthesis and add the length allowance.

TT prostheses: Distal allowance for movement of the prosthetic foot.

TF prostheses: Allowance proximal of the knee rotation point for flexion of the prosthetic knee joint and distal allowance for movement of the prosthetic foot.

2) Cut the pre-shaped foam cover to length and fit it in the proximal area on the prosthetic socket.

3) Pull the foam cover over the prosthesis.

4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or rests on the foot adapter.

5) Install the prosthetic foot on the prosthesis.

6) Mark the outer contour of the connecting element on the distal face of the foam cover.

7) Disassemble the prosthetic foot and remove the connecting element.

8) Clean the connecting element using a degreasing cleaner.

9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.

10) Let the glue dry (approx. **10 minutes**).

11) Install the prosthetic foot and adapt the exterior cosmetic shape. Take into account compression caused by cosmetic stockings or SuperSkin.

6 Cleaning

> **Allowable cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)

1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**

Clean the product with clear fresh water and a pH neutral soap.

2) Rinse the soap away with clear fresh water. In doing so, rinse the footshell until all dirt has been removed.

3) Dry the product with a soft cloth.

4) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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