9X25

DE Gebrauchsanweisung (Fachpersonal)	2
EN Instructions for Use (Qualified Personnel)	6
FR Instructions d'utilisation (Personnel spécialisé)	9
IT Istruzioni per l'uso (Personale tecnico specializzato)	13
ES Instrucciones de uso (Personal técnico especializado)	17
PT Manual de utilização (Pessoal técnico)	21
NL Gebruiksaanwijzing (Vakmensen)	25
SV Bruksanvisning (Fackpersonal)	29
DA Brugsanvisning (Faguddannet personale)	32
NO Bruksanvisning (Fagpersonell)	36
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	40
JA 取扱説明書 (有資格担当者)	44

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1 Timeliness of the Instructions for Use

INFORMATION

Date of the last update: 2013-01-22

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

2 Product Description


2.1 Description

The 9X25 Rocker Switch is suitable for incorporation in the definitive socket. The 9X25 Rocker Switch is suitable for patients who cannot generate enough muscle signal to control a myoelectric prosthesis but who want to take advantage of the high rehabilitation value of myoelectric prostheses.

Should the rocker switch be mounted in the provisional socket, the patient can practise the control or switching as early as the first stage of the rehabilitation. This will ensure a smooth transition to the use of the definitive prosthesis.

The rocker switch is used to control the prosthetic components with digital control and to switch between the prosthetic components. It can be used in combination with MYOBOCK electrodes, linear transducer and 4-step transducer.

2.2 Function

The rocker switch features 2 control points that are independently controlled. A function for control of the prosthetic components is assigned to each control point.

By pressing one side of the switch rocker, the switch will reach control point S1 and activate the control signal A. Control signal A is active.

By pressing the other side of the switch rocker, the switch will reach control point S2 and activate the control signal B. Now the control switch B is solely active.

3 Application

3.1 Indication

The 9X25 Rocker Switch is to be used **solely** in exoprosthetic fittings of the upper limbs for myoelectric control of arm prostheses.

3.2 Application

The 9X25 Rocker Switch is suitable for both right or left arm amputees as well as bilateral amputees.

3.3 Conditions of Use

The 9X25 Rocker Switch is suitable for indoor and outdoor use.

3.4 Qualification

The fitting of patients with the 9X25 Rocker Switch may only be carried out by a prosthetist.

4 Safety Instructions

Failure to follow the safety instructions given below can lead to faulty control or malfunction of the 9X25 Rocker Switch and result in a risk of injury for the patient.

5 General/Patient-Specific Safety Instructions

5.1 General Safety Instructions

- Check the rocker switch housing regularly for damage.
- When mounted into a provisional socket, the rocker switch may only be worn by an amputee within an orthopaedic workshop or as an inpatient in a rehabilitation facility.
- When installing the connection cable, make sure there are no tight angles or kinks.
- The rocker switch may not be used for bathing prostheses.
- Instruct the patient in the proper handling of the rocker switch according to the section "Patient-Specific Safety Instructions".

5.2 Patient-Specific Safety Instructions

- Handle the rocker switch with care and check the housing and pivoting mount for damage (e.g. cracks or breakage) before every use. Should the rocker switch be damaged or not function properly, please contact your prosthetist immediately.
- Driving a vehicle of any kind or operating machinery is strictly prohibited while wearing a provisional socket fitted with a rocker switch.
- Do not allow debris or liquids to get into the rocker switch.

- The rocker switch should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- The rocker switch and any damaged components may only be opened or repaired by certified Ottobock Myo-Service technicians.
- Regularly clean the rocker switch with a cloth slightly dampened with a soap solution. Do not allow moisture to enter into the rocker switch or into the contacts.

6 Delivery

6.1 Scope of Delivery

1 pc.	9X25 Rocker Switch
2 pc.	501S46=M2x8 Oval head screw
2 pc.	501S46=M2x5 Oval head screw
1 pc.	647G402 Instructions for Use

7 Preparation for Use

7.1 Connecting the Connection Cable

The rocker switch is delivered without connecting cable. It must be ordered separately. Appropriate cables are 13E50=*, 13E97=*, 13E98=1200 and 13E99=1200 Connecting Cables. Abundantly apply 633F11 Silicone Grease before connecting the contacts of the connecting cable. The grease will protect the connectors from sweat and humidity.

8 Disposal

This product must not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

9 Legal Information

9.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination

with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

9.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

9.3 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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