





8E38=9, 8E39=9, 8E41=9

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	19
<b>FR</b> Instructions d'utilisation .....	37
<b>IT</b> Istruzioni per l'uso .....	55
<b>ES</b> Instrucciones de uso .....	73
<b>PT</b> Manual de utilização .....	91
<b>NL</b> Gebruiksaanwijzing .....	107
<b>SV</b> Bruksanvisning .....	125
<b>DA</b> Brugsanvisning .....	141
<b>NO</b> Bruksanvisning .....	157
<b>HU</b> Használati utasítás .....	172
<b>TR</b> Kullanma talimatı .....	188
<b>RU</b> Руководство по применению .....	204
<b>EL</b> Οδηγίες χρήσης .....	222
<b>JA</b> 取扱説明書 .....	238
<b>ZH</b> 使用说明书 .....	253

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English

Date of the last update: 2021-03-30

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

## 1 Intended use

### 1.1 Medical Purpose

The MyoHand VariPlus Speed (8E38=9, 8E39=9, and 8E41=9) are to be used **solely** for the exoprosthetic fitting of the upper extremities.

### 1.2 Application

The MyoHand VariPlus Speed (8E38=9, 8E39=9, and 8E41=9) can be used for unilateral or bilateral amputees.

### 1.3 Conditions of use

The MyoHand VariPlus Speed (8E38=9, 8E39=9, and 8E41=9) have been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

The MyoHand VariPlus Speed (8E38=9\*, 8E39=9\* and 8E41=9\*) is designed exclusively for use on one patient. Use of the product by other persons is not approved by the manufacturer.

### 1.4 Prosthetist qualifications

The fitting of a patient with the MyoHand VariPlus Speed (8E38=9, 8E39=9, and 8E41=9) may only be carried out by prosthetists who have been trained and certified by Ottobock for the use of MyoBock components.

## 1.4 Contraindications

All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

## 2 Safety instructions

### Non-observance of safety instructions

Non-observance of the following safety instructions can lead to product malfunction or defects and result in a risk of injury to the patient.

- Do not use silicone spray when putting on the cosmetic glove; silicone spray does not allow for a proper fit. We recommend using Procomfort Gel (633S2) as a mounting aid.
- In order to reduce the malfunctioning risk of the MyoHand VariPlus Speed and the connected components due to strong electromagnetic radiation (e.g. when near high-tension power lines, transmitters or security systems for goods in department stores), set the electrode gains as low as possible.
- Before interrupting or establishing electrical connections (such as when removing the Hand from the prosthesis), the battery must be removed from the battery receptacle to disconnect the system from the power supply.
- If you do not achieve the expected results by adjustment and by selection of the appropriate control mode, please contact an Ottobock Myo-Service facility.
- Use the "Patient Information" section of this pamphlet (Section 3) to inform your patient about the function and use of the MyoHand VariPlus Speed.

## 3 Patient information

- Do not allow foreign particles or liquids to get into the MyoHand VariPlus Speed mechanism.
- The MyoHand VariPlus Speed should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Avoid staying near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores), as this can lead to malfunction of the MyoHand VariPlus Speed.
- The MyoHand VariPlus Speed has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Myo-Service.
- The MyoHand VariPlus Speed and MyoHand VariPlus Speed components may be opened or repaired only by certified Ottobock Myo-Service technicians.
- If the MyoHand VariPlus Speed will not be used for an extended period of time, it should be left in the open position to protect the mechanism.
- Please provide your patient with the 646D165 Patient Information on the MyoHand VariPlus Speed.
- If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WLAN devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour. It is recommended to maintain the following minimum distances:
  - Mobile phone GSM 850/GSM 900: 0.99m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7m

- DECT cordless phones incl. base station: 0.35m
- WLAN (routers, access points,...): 0.22m
- Bluetooth devices (third-party products not approved by Ottobock): 0.22m
- When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

### 3.1 Driving vehicles

An upper extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test center. For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork). It is imperative that the driver is able to operate the vehicle without any risk with the MyoHand VariPlus Speed turned off. Driving with the MyoHand VariPlus Speed turned ON may present a risk if the MyoHand VariPlus Speed inadvertently moves due to unintentional muscle contraction or other causes.

#### Caution:

Before a MyoHand VariPlus Speed with a quick-disconnect wrist can be used to drive a vehicle, it must be positioned in such a way that a slight turn of the Hand, which can occur during steering, will not disconnect the Hand from the prosthesis.

## 4 Scope of delivery

1 pc. of 8E\*\*=9 MyoHand VariPlus Speed

1 pc. Instructions for Use

## 5 Description and function

The MyoHand VariPlus Speed is a myo-electrically controlled prosthetic hand that is characterised by a high grip speed in connection with an innovative and sensitive control concept. The MyoHand VariPlus Speed is equipped with a programming function.

Different proportional and digital control modes with 1 or 2 electrodes, linear transducers or switches allow individual adaptation to the respective patient.

The MyoHand VariPlus Speed is based on the Ottobock DMC System (DMC = Dynamic Mode Control): This system, which has been developed by Ottobock, uses two independent measurement and control systems in order to optimally adapt the grip speed and grip force to the patient's muscle signal.

The proportional DMC-Control (Dynamic Mode Control) allows the amputee to control grip speed and grip force proportionally to his muscle signal. As soon as the strength of the muscle signal changes, grip speed and grip force immediately adapt to the new muscle signal.

The coaxial bushing of the MyoHand VariPlus Speed (8E38=9) is marked with a green ring; the connection cable of the MyoHand VariPlus Speed (8E39=9 or 8E41=9) is marked with a green bushing.

## 5.1 Programming function

For optimal adaptation to the patient, six different control modes are available. Connect the MyoSelect (757T13) to select the control mode.

## 5.2 Switching the MyoHand VariPlus Speed ON and OFF (Fig. 1, Fig. 2)

By pressing the cosmetic glove you can actuate the ON/OFF switch integrated into the control electronics.

By pressing the back of the hand:                      Function ON                      (Fig. 1)

By pressing the thumb:                                      Function OFF                      (Fig. 2)

## 6 Adjusting the electrode(s)

The best possible function of the MyoHand VariPlus Speed is achieved through optimal positioning and adjustment of the electrodes. This is easy with the MyoBoy (757M10) or with the MyoBoy (757M11) with USB Port.

### Attention!

Allow the patient to rest during the adjustment of the electrodes to avoid fatigue and thus incorrect electrode settings (i.e. excess sensitivity). The electrodes should be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interferences from electric devices, the position of the electrodes should be verified and, if need be, repositioned. If the interference cannot be eliminated, please contact the Ottobock Myo-Service.

## 7 Adjusting the control modes (Fig. 3 to Fig. 6)

For optimal adaptation to the patient there are 6 control modes, which can be selected and adjusted using the MyoSelect (757T13). The MyoHand VariPlus Speed is delivered with control mode 1, DMC plus, activated.

### Control modes 1, 3, and 4:

Adjust the electrodes so that the patient can maintain the respective electrode signal for approximately two seconds at a point above the value HIGH (Fig. 3).

### Control mode 2:

Adjust the CLOSING electrode so that the patient can maintain the electrode signal for approximately two seconds at a point above the value ON (Fig. 2). Adjust the OPENING electrode so that the patient can maintain the electrode signal for approximately two seconds at a point above the value LOW (Fig. 4).

### Control mode 5:

Adjust the electrodes so that the patient can maintain the respective electrode signal for approximately two seconds at a point above the value ON (Fig. 5).

### Control mode 6:

The upper threshold (see control mode description) can be adjusted using the Ottobock MyoSelect (757T13).

In the PAULA program, select "Digital Twin Hand without function plug". Adjust the slider on the scale to 3.5. The MyoHand VariPlus Speed is delivered with the **adjustment regulator** of the electronic control already set to 3.5. Now adjust the electrode so that the electrode signal reaches a value just below the threshold – marked in PAULA with a broken line at 3.5 (Fig. 6): The MyoHand

VariPlus Speed closes. Exceed the threshold with a quick and strong electrode signal to open the MyoHand VariPlus Speed.

If the electrode signal is too slow or too weak, reduce the threshold in PAULA (slider). Then set the adjustment regulator C of the electronic control to the same value using MyoSelect.

## 8 Overview of control modes for the MyoHand VariPlus Speed

Control mode														
		No.												
Control mode														
Name														
		1	2		3	4	5		6					
		DMC plus	AutoControl LowInput		VarioControl	VarioDual	Digital Control		Double Channel Control					
Grip force		Proportional to signal decline				●								
		Time proportional		●	●	●		●	●	●				
		Proportional to signal strength	●					●						
Speed		CLOSING		Digital (sustained above threshold)		●	●	●						
				Proportional to signal decline				●						
				Proportional to signal strength	●				●					
		OPENING		Digital (sustained above threshold)				●		●	●	●	●	
				Proportional to signal increase				●						
				Proportional to signal strength	●	●	●			●				
Electrode signal strength (according to MyoBoy scale)		CLOSING		ON		●				●				
				LOW									●	
				HIGH	●				●		●			
		OPENING		ON						●	●			
				LOW		●	●							●
				HIGH	●				●		●			
Signal sources		Switch				●	●			●	●			
		Linear transducer						●						
		One electrode				●		●			●		●	
		Two electrodes		●	●					●	●			

The terms ON, LOW and HIGH correspond to the respective scale values of MyoBoy and PAULA.

## 9 Description of control modes

### 9.1 Control mode 1: **DMC plus** Control with 2 electrodes

The grip speed and grip force are determined by the strength of the electrode signal (resulting from the muscle contraction). With another, stronger electrode signal, the amputee can always increase the grip force up to the maximum of approx. 100 N.

**Opening:** Proportionally controlled by the OPENING electrode.


**Closing:** Proportionally controlled by the CLOSING electrode.

Example 1: A gentle electrode signal applies the smallest grip force to the object.

Example 2: A stronger electrode signal creates a more powerful grip force. Another, stronger electrode signal increases the grip force up to the maximum of approx. 100 N.

Control mode 1	OPEN	CLOSE	Indication
DMC plus	Sustained electrode signal  <b>Speed:</b> proportional Adjustment regulator A	Sustained electrode signal <b>Grip force:</b> proportional After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed with unwanted electrode signals is prevented. <b>Speed:</b> proportional Adjustment regulator A	For patients with two strong electrode signals

### Display of the MyoSelect (757T13)

Adjustment regulator A is used to adjust the proportional speed of the “OPENING /CLOSING” function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

**9.2 Control mode 2:**

**AutoControl - LowInput  
Control with 2 electrodes,  
1 electrode and 1 switch or  
1 switch**

The MyoHand VariPlus Speed closes at constant speed, and the grip force increases proportionally to the duration of gripping. The closing speed and the speed of grip force build-up can be adjusted with adjustment regulator B.


**9.2.1 Control with 2 electrodes**

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** At constant speed by muscle signal of any level above the ON-threshold to the CLOSING electrode.

Control mode 2	OPEN	CLOSE	Indication
AutoControl LowInput	<p>Sustained electrode signal</p> <p>Reduced range of proportionality: Maximum speed once the LOW-threshold has been reached</p> <p><b>Speed:</b> proportional Adjustment regulator A</p>	<p>Sustained electrode signal</p> <p><b>Grip force:</b> Time proportional After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed with unwanted electrode signals is prevented.</p> <p><b>Speed:</b> constant Adjustment regulator B</p>	<p>For amputees with two weak electrode signals</p>

**Display of the MyoSelect (757T13)**

Adjustment regulator A is used to adjust the proportional speed of the "OPENING" function.

Adjustment regulator B is used to adjust the digital speed of the "CLOSING" function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.


### 9.2.2 Control with 1 electrode and 1 switch

**Opening:** Proportionally controlled by the OPENING electrode.

**Closing:** At constant speed by operation of the switch.

Control mode 2	OPEN	CLOSE	Indication
AutoControl LowInput	<p>Sustained electrode signal</p> <p>Reduced range of proportionality: Maximum speed once the LOW-threshold has been reached</p> <p><b>Speed:</b> proportional Adjustment regulator A</p>	<p>Signal from the switch</p> <p><b>Grip force:</b> Time proportional After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed with unwanted electrode signals is prevented.</p> <p><b>Speed:</b> constant Adjustment regulator B</p>	<p>For amputees with only one muscle and weak electrode signal</p>

#### Display of the MyoSelect (757T13)

Adjustment regulator A is used to adjust the proportional speed of the “OPENING” function.

Adjustment regulator B is used to adjust the constant speed of the “CLOSING” function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

### 9.2.3 Control with 1 switch


This control mode can be used with any MyoBock switch.

**Opening:** At constant speed, as long as the OPEN-side of the switch is operated.  
The MyoHand VariPlus Speed remains open.

**Closing:** At constant speed, as long as the CLOSE-side of the switch is operated.

Control mode 2	OPEN	CLOSE	Indication
AutoControl LowInput	MyoHand VariPlus Speed opens as long as the OPEN-side of the switch is operated.  <b>Speed:</b> constant Adjustment regulator A	MyoHand VariPlus Speed closes as long as the CLOSE-side of the switch is operated.  <b>Grip force:</b> Time proportional <b>Speed:</b> constant Adjustment regulator B	For amputees with weak or no electrode signal

### Display of the MyoSelect (757T13)

Adjustment regulator A is used to adjust the constant speed of the “OPENING” function.

Adjustment regulator B is used to adjust the constant speed of the “CLOSING” function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

### 9.3 Control mode 3: **VarioControl** Control with 1 electrode or 1 linear transducer

With this control mode, the opening speed is determined by the strength and speed of the muscle contraction. The closing speed and grip force depend on the reduction of muscle tension.

#### 9.3.1 Control with 1 electrode

**Opening:** Proportional. The opening speed is determined by the speed and strength of the muscle contraction.

**Closing:** Proportional. The closing speed is determined by the speed and strength of the muscle relaxation.


**Stopping:** By very slow muscle relaxation to the electrode; the MyoHand VariPlus Speed stops in an opened position.

Example 1: After opening the amputee slowly relaxes the muscle. Closing takes place proportional to the decline of the electrode signal. The object is grasped with low force.

Example 2: After opening the amputee relaxes the muscle at maximum speed. The MyoHand VariPlus Speed closes at maximum speed and begins to grasp an object with maximum grip force (approx. 100 N).

Control mode 3	OPEN	CLOSE	Indication
VarioControl	<p>Increasing electrode signal through muscle contraction</p> <p>Speed and strength of muscle contraction to the electrode</p> <p><b>Speed:</b> proportional Adjustment regulator A</p>	<p>Declining electrode signal through muscle relaxation</p> <p><b>Grip force:</b> Proportional to the decline of the electrode signal</p> <p>After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed with unwanted electrode signals is prevented.</p> <p><b>Speed:</b> proportional Adjustment regulator A</p>	<p>For amputees with one strong electrode signal or tendency to co-contraction</p>

### Display of the MyoSelect (757T13)

Adjustment regulator A is used to adjust the proportional speed of the “OPENING / CLOSING” function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

### 9.3.2 Control with 1 linear transducer

**Opening:** Proportional. The speed and strength of the pull on the linear transducer determine the opening speed of the MyoHand VariPlus Speed.

**Closing:** Proportional. The speed of the release of the linear transducer determines the closing speed.

**Stopping:** By very slow release of the linear transducer the MyoHand VariPlus Speed stops in an opened position.


Example 1: After opening the amputee slowly releases the pull on the linear transducer. The object is grasped with low force.

Example 2: After opening the amputee releases the pull on the linear transducer at maximum speed. The MyoHand VariPlus Speed closes at maximum speed and grasps an object with a grip force of approx. 100 N.

Control mode 3	OPEN	CLOSE	Indication

VarioControl	<p>Speed and strength of pull on the linear transducer</p> <p><b>Speed:</b> proportional Adjustment regulator A</p>	<p>Speed of release of pull on the linear transducer</p> <p><b>Grip force:</b> Proportional to the release of the pull on the linear transducer</p> <p>After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the Myo-Hand VariPlus Speed with unwanted electrode signals is prevented.</p> <p><b>Speed:</b> proportional Adjustment regulator A</p>	For amputees with weak or no electrode signal
--------------	---	---	---

### Display of the MyoSelect (757T13)

Adjustment regulator A is used to adjust the proportional speed of the “OPENING/CLOSING” function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

### 9.4 Control mode 4: VarioDual Control with 2 electrodes

With this control mode, the opening speed is determined by the strength and speed of the muscle contraction. The closing speed depends on the speed of the muscle relaxation. The following or simultaneous electrode signal to the second electrode determines the grip force. Another, stronger electrode signal to the second electrode increases the grip force.

#### Electrode 1:

**Opening:** Proportional. The opening speed is determined by the speed of the muscle contraction.

**Closing:** Proportional. The closing speed is determined by the speed of the muscle relaxation.

**Stopping:** By very slow muscle relaxation to the electrode; the MyoHand VariPlus Speed stops in an opened position.

#### Electrode 2:

**Grasping:** The level of the electrode signal to the second electrode determines the build-up of the grip force. The maximal grip force is approx. 100 N.


Example 1: After opening, the amputee relaxes the muscle at any speed. The MyoHand VariPlus Speed closes proportionally to the speed of the muscle relaxation. The object is grasped with minimal force.

Example 2: After grasping as described in example 1, the MyoHand VariPlus Speed is to grasp the object with a higher grip force. For this purpose, the amputee must create an

electrode signal to the second electrode. Proportionally, a grip force of up to approx. 100 N can be built up.

Control mode 4	OPEN	CLOSE	Indication
VarioDual	<p>Increasing electrode signal through muscle contraction to the first electrode</p> <p>Speed and strength of muscle contraction to the electrode</p> <p><b>Speed:</b> Proportional to the increase of the electrode signal Adjustment regulator A</p>	<p>Declining electrode signal through muscle relaxation to the first electrode or: Sustained electrode signal to the second electrode</p> <p><b>Grip force:</b> Proportional to the signal strength to the second electrode. If only the first electrode is used, the MyoHand VariPlus Speed closes up to the minimal grip force.</p> <p><b>Grip force increase:</b> By another, stronger electrode signal to the second electrode. After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed with unwanted electrode signals is prevented.</p> <p><b>Speed:</b> Proportional to the decline of the electrode signal to the first electrode / proportional to the signal strength to the second electrode Adjustment regulator A</p>	<p>For amputees with two strong electrode signals</p>

### Display of the MyoSelect (757T13)

Adjustment regulator A is used to adjust the proportional speed of the "OPENING/CLOSING" function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

### 9.5 Control mode 5:

#### Digital Control Control with 2 electrodes, 1 electrode and 1 switch or 1 switch

The MyoHand VariPlus Speed closes at constant speed, and the grip force increases proportionally to the duration of gripping. The OPENING/CLOSING speed and the speed of grip force build-up can be adjusted with adjustment regulator B.


### 9.5.1 Control with 2 electrodes

**Opening:** Digital, signal to the OPENING electrode opens the MyoHand VariPlus Speed at constant speed.

**Closing:** Digital, signal to the CLOSING electrode closes the MyoHand VariPlus Speed at constant speed.

Control mode 5	OPEN	CLOSE	Indication
Digital Control	Sustained electrode signal  Adjustment regulator B	Sustained electrode signal  <b>Grip force:</b> Duration of the signal Adjustment regulator B	For amputees with two weak electrode signals

#### Display of the MyoSelect (757T13)

Adjustment regulator B is used to adjust the constant speed of the "OPENING/CLOSING" function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.


### 9.5.2 Control with 1 electrode and 1 switch

**Opening:** Electrode signal to the OPENING electrode.

**Closing:** At constant speed, as long as the switch is operated.

Control mode 5	OPEN	CLOSE	Indication
Digital Control	Sustained electrode signal  Adjustment regulator B	Signal from the switch  <b>Grip force:</b> Duration of the signal Adjustment regulator B	For amputees with only one muscle and weak electrode signal

#### Display of the MyoSelect (757T13)

Adjustment regulator B is used to adjust the speed of the "OPENING/CLOSING" function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

### 9.5.3 Control with 1 switch

This control mode can be used with any MyoBock switch.


**Opening:** At constant speed, as long as the OPEN-side of the switch is operated. The MyoHand VariPlus Speed remains open.

**Closing:** At constant speed, as long as the CLOSE-side of the switch is operated.

Control mode 5	OPEN	CLOSE	Indication
Digital Control	MyoHand VariPlus Speed opens as long as the OPEN-side of the switch is operated.  Adjustment regulator B	MyoHand VariPlus Speed closes as long as the CLOSE-side of the switch is operated.  <b>Grip force:</b> Duration of the signal Adjustment regulator B	For amputees with weak or no electrode signals

**Display of the MyoSelect (757T13)**



Adjustment regulator B is used to adjust the constant speed of the “OPENING/CLOSING” function and the speed of grip force build-up.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

**9.6 Control mode 6: Double Channel Control Control with 1 electrode**


With this control mode, a quick and strong signal opens the MyoHand VariPlus Speed, and a slow, gently signal closes it.

Control mode 6	OPEN	CLOSE	Indication

<p>Double Channel Control</p>	<p>Quick, strong electrode signal that reaches the upper threshold within 80 ms and is maintained above the upper threshold for at least 30 ms</p> <p>Adjustment regulator B</p> 	<p>Slow, gentle electrode signal that does not reach the upper threshold within 80 ms</p> <p><b>Grip force:</b> Duration of the signal</p> <p>Adjustment regulator B</p> 	<p>For amputees with one strong electrode signal</p>
-------------------------------	--	--	--

**Notice:** The upper threshold (B) can be adjusted using the Ottobock MyoSelect (757T13) (adjustment regulator C).

### Display of the MyoSelect (757T13)

Adjustment regulator B is used to adjust the speed of the “OPENING/CLOSING” function and the speed of grip force build-up.

Adjustment regulator C is used to adjust the upper threshold.

**Notice:** Always check the adjustments on the prosthesis and adapt them to the amputee.

## 10 Handling the Batteries

As soon as the MyoHand VariPlus Speed recognizes a fully charged EnergyPack (757B20, 757B21) or MyoEnergy Integral (757B25=\* / 757B35=\*), it automatically switches to Li-Ion battery technology. Once this has happened, the system should only be operated with this battery technology. If a 757B15 X-ChangePack is employed after the control has switched to Li-Ion battery technology, the battery can not be used to its full capacity. The 757T13 Ottobock MyoSelect can be used, however, to reset the MyoHand VariPlus Speed control for use with X-ChangePack. For operational security and reliability, only 757B20, 757B21, 757B25=\* / 757B35=\* or 757B15 Batteries are to be used.

The default setting is optimized for NiMH Batteries.

## 11 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

The grace period is no more than one month before or three months after maintenance is due.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components must always be sent in for maintenance and repairs:

The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

## 12 Disposal

This product may not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

## 13 Technical data

Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Static current	1 mA
Opening width	100 mm
Proportional speed	15–300 mm/s
Proportional grip force	15–300 mm/s

Power supply	0 – approx. 100 N 757B20 EnergyPack (7.2 V)
Service life of the gripping component	5 years
Battery service life	2 years

***With limitations:***

757B21 EnergyPack	(7.2 V)
757B25=* / 757B35=* MyoEnergy Integral	(7.4 V)
X-ChangePack (757B15), Interchangeable Battery NiMH	(6.0 V)

## 14 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer assumes no liability for damage caused by component combinations (use of third-party products) and applications which it did not authorise.

The product may only be opened and repaired by authorized Ottobock technicians.

## 15 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective-applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 16 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 17 Appendices

### ***17.1 Operating states/error signals***

#### **Battery charge level**

An intelligent "battery manager" alerts the patient to low voltage conditions – the hand will become slower and the gripping force will decrease. At the same time the battery is protected against harmful deep discharge.

A second charged battery should be kept on hand as a replacement.

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