


ottobock.


MovolinoArm Friction 12K12

DE	Gebrauchsanweisung	5
EN	Instructions for use	13
FR	Instructions d'utilisation	20
IT	Istruzioni per l'uso	28
ES	Instrucciones de uso	36
PT	Manual de utilização.....	44
NL	Gebruiksaanwijzing	52
SV	Bruksanvisning	60
DA	Brugsanvisning.....	67
NO	Bruksanvisning	74
FI	Käyttöohje.....	81
PL	Instrukcja użytkowania.....	89
HU	Használati utasítás	97
CS	Návod k použití.....	104
RO	Instrucțiuni de utilizare.....	111
HR	Upute za uporabu	119
TR	Kullanma talimatı.....	126
EL	Οδηγίες χρήσης	134
RU	Руководство по применению	142
JA	取扱説明書.....	151
ZH	使用说明书.....	157

INFORMATION

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INFORMATION

- ▶ Date of last update: 2020-12-07
- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

Explanation of symbols

English

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

1 Single components and accessories

1.1 Single components

For available single components refer to the catalogue.

1.2 Accessories

For myoelectric fittings	For cable-controlled prosthesis	For passive prosthesis
13E129=G100/G300 Electrode cables with straight plug and 13E121 plug connection		
11D12=34 Lamination ring for Electric Hand 2000	10V18=34 or 10V36=34 Ottobock Wrist Unit	10A40 Hand Adapter Wood

For myoelectric fittings	For cable-controlled prosthesis	For passive prosthesis
757B35=1 MyoEnergy Integral		
9E420=* 7in1 Controller		

2 Description

2.1 Indications for use

The Ottobock MovolinoArm Friction 12K12 **may only be used** for the exoprosthetic fitting of the upper extremities in children in combination with Ottobock components.

2.2 Field of application

The Ottobock MovolinoArm Friction 12K12 may **only** be used for passive, cable-controlled and myoelectrically controlled arm prostheses for children.

2.3 Safety instructions

Please make your patients and their parents / guardians aware of the following safety instructions:

CAUTION

Risk of injury as a result of excessive strain. The Ottobock prosthetic components described here have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety.

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

NOTICE

Risk of corrosion. Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids.

Using this medical device in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

Risk of damage due to incorrect environmental conditions.

Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthesis and its components. Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

NOTICE

Risk of damage caused by improper cleaning agents. The use of aggressive cleaning agents must be avoided. They may cause damage to bearings, seals and plastic parts. The product may only be cleaned with a damp cloth and mild soap.

2.4 Function

The MovolinoArm Friction is available in one size. The elbow component can be used for both right and left-side fittings. It is suitable for children ages 3-5.

The MovolinoArm Friction has one friction setting for humeral rotation and one for flexion / extension of the forearm.

Extension/flexion starts at 0–5° pre-flexion. The forearm can be flexed by a maximum of 140°. The terminal flexion/ extension stops are dampened. The friction can be adjusted via the adjustment screw using the Allen key (wrench size 3 mm) included in the scope of delivery. Adjustment should be performed at maximum flexion (fig. 1) or maximum extension (fig. 2) of the forearm. During adjustment, please consider the individual strengths, moving abilities and needs of the child.

Humeral rotation is limited to a total angle of 160°. The friction can be adjusted via the adjustment screw (fig. 3) using the Allen key (wrench

size 3 mm) included in the scope of delivery. During adjustment, please consider the individual strengths, moving abilities and needs of the child.

3 Technical data

Article number	12K12
Weight [g]	180
Upper arm connection Ø [mm]	43.5
Wrist connection Ø [mm]	34
Length by which the forearm component can be shortened relative to the elbow pivot point, from - to [mm]	165 – 115

4 Handling

4.1 Installation preparation

The two-socket procedure with outer and inner sockets is recommended for fittings using the 12K12. During the two-socket procedure, the lamination ring is laminated into the outer socket. This allows for easy assembly / disassembly of the elbow component.

4.2 Assembly

For a unilateral fitting, the end of the thumb and the pivot point of the elbow joint are relevant to determine the length of the components. In case of bilateral fittings, priority is given to the prosthesis' functionality. When adapting the length, make sure that the prosthesis will offer the patient a high degree of freedom of movement.

Changes in length can be realized by shortening the elbow joint on the forearm side (see table chapter 3) and by changing the upper arm socket.

The lamination ring must be separated from the component before lamination. The screw is removed from the elbow joint (fig. 4 arrow) and the lamination ring with the cover and screw is pulled off. The cover and screw are not required for lamination and must be removed from the lamination ring.

Then assemble the lamination ring using the two lamination dummies

included in the kit (fig. 5). Assemble the black lamination dummy first. This ensures proper orientation. Figure 6 intentionally illustrates incorrect installation. If the marked lamination dummy does not protrude evenly above the lamination ring, then it can be aligned by turning the dummy (fig. 6 arrows). The lamination dummy locks into place audibly, providing acoustic feedback of the correct position. The markings allow for correct alignment of the lamination ring when the outer socket is made. Dummy marking A indicates the Anterior side, marking P indicates the Posterior side. Proper alignment of the lamination ring is essential for the physiologically optimal adjustment of the humeral rotation of the elbow joint. Subsequently, assemble the white lamination dummy. To do this, the dummy is to be pressed into the lamination ring until a recess develops. This recess constitutes the force-fit seal of the lamination ring. When orienting the MovolinoArm Friction for humeral rotation, maximum benefit of the prosthesis for the user ought to be provided. For close-to-body activities it is advisable to allow more internal rotation in proportion to external rotation.

Care must be taken to ensure that the upper arm socket is fabricated to be as circular and slender as possible in order to attain the maximum possible flexion angle.

The included o-ring is installed over the still-soft laminate below the dummy after lamination. It provides for a good transition between the lamination resin and the lamination ring (fig. 8).

After the lamination has hardened, the lamination dummies are removed from the lamination ring. Then the MovolinoArm Friction 12K12 can be mounted.

NOTICE

No disassembly of the forearm and elbow ball:

The forearm and the elbow ball are connected by a cable and inseparable. Disassembly of the forearm is not possible!

4.3 Adjustments and final assembly

4.3.1 Myoelectric prosthetic fitting

To install the myoelectrical components in the forearm, the assembly cover must be removed from the forearm. To do so, loosen the screw using the

Torx screwdriver 710Y12=T10. After carefully lifting the assembly cover, it can be removed in the direction of the elbow ball (fig. 9).

The following procedure is recommended for installation of the myoelectrical components.

Installation of MyoEnergy Integral 757B35=1

The charging receptacle contour inside the forearm component is recessed to allow for installation of the MyoEnergy Integral charging receptacle. This contours must be ground free. The recesses which are intended as guides for drilling can be utilized with an 8 mm (5/16") twist drill. Then the charging receptacle can be mounted inside the crook of the arm according to instructions. Do not bend the flat ribbon cable which runs from the battery to the charging receptacle (fig. 11)!

Roughen the inside of the forearm socket with 634A3 acetone in order to glue the lamination ring 11D12=34. Allow the acetone to evaporate prior to gluing. Use 636K18 Sealing Resin Compact Glue.

After the lamination ring 11D12=* bonding has hardened, install the 7in1 Controller 9E420 contact plate. Press the contact plate firmly into the lamination ring.

Installation of the electric cables in the forearm.

The electrode cables are connected in the proximal section of the elbow joint (fig. 12). The electrode signals are transmitted to the forearm socket. Lift the circuit board out of the forearm (fig. 10). Plug the 13E129=G100/G300 electrode cable into the circuit board while observing the contact markings. Carefully reinsert the circuit board in its initial position into the forearm.

This electrode cable transmits the electrode signals to the 9E420=* 7-in-1 Controller.

Connect all cables to the 7in1 Controller 9E420 as specified.

Insert the assembly cover carefully and install. Make sure the cables are placed properly. Cables may not be kinked or pinched.

To assemble the elbow component, pull the electrode cables through the outer socket of the upper arm and the plug hole in the cover. Then the cover and the countersunk head screw (fig. 4) are mounted. Both parts are routed through the inner part of the outer socket and connected to the lower part of the elbow component. The countersunk head screw must be tightened with a torque of 5 Nm.

NOTICE

Risk of damage due to improper handling.

The MovolinoArm Friction 12K12 must not be flexed during charging and use of the MyolinoLink 60X6.

4.3.2 Passive Prosthetic Fittings

The 10A40 wooden adapter is used in the fabrication of a passive prosthesis. The adapter is to be milled down to a diameter of 34mm. Roughen the inside of the forearm socket with 634A3 acetone. Allow the acetone to evaporate prior to gluing the wooden adapter into the forearm with 634K18 Sealing Resin Compact Adhesive.

Alignment of the hand should always be performed in the neutral-zero position. However, the individual conditions and needs of the patients must be considered in all cases.

4.3.3 Body-Powered Prosthetic Fitting

For a body-powered prosthetic fitting, 10V18=34 or 10V36=34 Wrist Unit is being used. After separating the wrist unit, the lamination ring is glued into the roughened inside of the forearm using 636K18 Sealing Resin Compact Adhesive.

After performing the prosthetic fitting (chapter 4.3.2 to 4.3.3), the cover and the countersunk head screw (fig. 12) are mounted. Both parts are routed through the inside of the outer socket and connected to the lower part of the elbow component. The countersunk head screw is tightened with a torque of 5 Nm.

5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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