## AxonMaster 13E500

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	2
<b>EN</b> Instructions for use (qualified personnel) .....	12
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	23
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	34
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	44
<b>PT</b> Manual de utilização (Pessoal técnico) .....	54
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	64
<b>SV</b> Bruksanvisning (Fackpersonal) .....	74
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	83
<b>NO</b> Bruksanvisning (Fagpersonell) .....	93
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	103

## INFORMATION

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### INFORMATION

Last update: 2015-03-09

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product Description

### 2.1 Function

The 13E500 AxonMaster is the central control unit in the Axon-Bus prosthesis system. The AxonMaster samples the patient's control signals and routes them to the corresponding prosthesis components via the Axon-Bus (joint control and power supply connection).

This allows the user to control the prosthesis and switch between the prosthesis components. The AxonMaster also controls the data communication process of the Axon-Bus system.

Depending on the Axon-Bus prosthesis components currently used, specific standard programmes are available. With the 560X500 AxonSoft adjustment software, these standard programmes can be adapted/modified for the specific needs of each respective patient.

### 2.2 Combination Possibilities

A maximum of 3 of the following control devices can be connected to the AxonMaster:

- 13E200=\* MyoBock Electrode
- 13E202=\* Suction Socket Electrode
- 9X50/9X52 Linear Transducer
- Cable Pull Switch

## 3 Application

### 3.1 Indications for Use

The 13E500 AxonMaster is intended **exclusively** for exoprosthetic fittings of the upper limbs and serves to control and coordinate all active Ottobock Axon-Bus prosthetic components within the system.

### 3.2 Use/Field of Application

The 13E500 AxonMaster may only be used in conjunction with components of the Ottobock Axon-Bus prosthetic system.

### 3.3 Conditions of Use

The 13E500 AxonMaster was developed for everyday activities only.

### 3.4 User/Qualifications

Fitting a user with the 13E500 AxonMaster may only be carried out by prosthetists who have been authorised by Ottobock after completion of a corresponding training course.

## 4 Safety

### 4.1 Explanation of Warning Symbols

#### **WARNING**

Warnings regarding possible risks of severe accident or injury.

#### **CAUTION**

Warnings regarding possible risks of accident or injury.

#### **NOTICE**

Warnings regarding possible technical damage.

### 4.2 General Safety Instructions

#### **CAUTION**

#### **Incorrect electrode settings**

Injury due to faulty control or malfunction of the product.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact Ottobock Myo-Service.
- ▶ Set the electrode gains as low as possible to reduce the risk of malfunctioning due to strong electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ Allow the patient to rest during the adjustment of the electrodes; otherwise, muscle fatigue will lead to inconsistent results. As a result, the therapist will tend to set the electrode gains too high.

**⚠ CAUTION****Incorrect assembly**

Injury due to faulty control or malfunction of the product.

- ▶ During assembly, pay special attention to the section "Preparation for Use" in these Instructions for Use.

**⚠ CAUTION****Use of unsuitable prosthesis components**

Injury due to faulty control or malfunction of the product.

- ▶ Only combine the product with prosthesis components approved by Ottobock (see the Instructions for Use and catalogues).

**⚠ CAUTION****Penetration of liquids**

Injury due to faulty control or malfunction of the product.

- ▶ Make sure that no liquids can get into the product.
- ▶ Do not use the product for bathing prostheses.

**⚠ CAUTION****Incorrect laying of the connecting cables**

Injury due to faulty control or malfunction of the product as a result of damage to the connecting cables due to incorrect laying.

- ▶ When installing the cables, make sure there are no tight angles or kinks.

**NOTICE****Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Product malfunction due to interference during data transfer.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.50 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.35 m
  - DECT cordless phones incl. base station: 0.18 m
  - WiFi (routers, access points,...): 0.11 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.11 m

**NOTICE****Improper care**

Damage to the casing through the use of acetone, petrol or similar solvents.

- ▶ Only use a damp cloth and mild soap to clean the product.

**INFORMATION**

Prior to initial operation of the product, charge the supplied battery.

**INFORMATION**

The product serves to adapt the prosthesis to individual user requirements. Please be aware that only a carefully adjusted prosthesis will improve the user's quality of life. If you do not

achieve the expected results by adjustment and by selection of the appropriate programme, please contact your prosthetist or Ottobock Myo-Service.

## **5 Delivery/Storage**

### **5.1 Scope of Delivery**

- 1 pc. 13E500 AxonMaster
- 1 pc. Lamination dummy
- 1 pc. Casting template
- 1 pc. 647G590 Instructions for Use (Qualified Personnel)

### **5.2 Storage**

Observe the specifications in the chapter "Technical Data".

## **6 Preparation for Use**

### **6.1 Installation in the forearm/upper arm**

The 13E500 AxonMaster can be positioned in the forearm socket or in the upper arm socket in the case of a fitting with elbow joints 12K500 AxonArm Hybrid/12K501 AxonArm Ergo, for example.

The rounded design allows it to be secured in the outer socket or between the outer and inner socket, depending on space availability. If there is sufficient space between the outer and inner socket, the 13E500 AxonMaster can be secured using a hook-and-loop closure strap. When little space is available, the lamination dummy of the 13E500 AxonMaster is laminated into the outer socket.

### **6.2 Power supply**

Power is supplied to the AxonMaster and the prosthesis components by the 757B500/757B501 AxonEnergy Integral.

## **7 Handling**

### **7.1 Establishing a connection**

Adjustments can be made to the prosthesis components via Bluetooth data transfer and the 560X500 AxonSoft software. For this purpose, a wireless connection must be established between the 13E500 AxonMaster and the PC with the help of the 60X5 BionicLink PC.

### **7.2 Selecting a program**

For optimum adaptation to the patient, various standard programmes are provided as well as the option of making refined patient-specific settings. These standard programmes can be selected using the 560X500 AxonSoft adjustment software and include predefined control options and switching methods.

Two 13E200/202=\* MyoBock Electrodes and/or a cable pull switch can be connected to the 13E500 AxonMaster. The electrode for closing the hand is connected to input 1, the electrode for opening the hand to input 2. The cable pull switch for a pull cable support, for example, connects to input 3. The remaining inputs cannot be used at this time and therefore remain closed with cover caps.

### **7.3 Feedback Messages/Operating States**

An intelligent battery manager of the 13E500 AxonMaster alerts the user to low voltage conditions of the battery – the hand will become slower and the gripping force will decrease. Finally, when there is very little battery capacity remaining, the system will switch off. This protects the battery against harmful deep discharge.

## 8 Control Options

### 8.1 2-Multi

Proportional control with 2 electrodes.

The neutral position of the Michelangelo Hand is the initial position for opposition mode and lateral mode.

The gripping speed and gripping force are determined by the strength of the electrode signal (resulting from the muscle contraction). The gripping force can be increased up to the maximum gripping force at any time through a stronger electrode signal.

Starting from an open hand position, the Michelangelo Hand automatically switches to the neutral position when the hand is not controlled via the electrode signals. While closing to the neutral position, the system builds up the minimum gripping force (approx. 15 N).

#### **Opening:**

Proportionally controlled by the OPENING electrode (e.g. input 2)

#### **Closing:**

Proportionally controlled by the CLOSING electrode (e.g. input 1)

#### **Example 1:**

The hand is opened by an OPEN signal in opposition mode. After the OPEN muscle has relaxed, the hand slowly moves to the neutral position.

#### **Example 2:**

The hand is in the neutral position. Upon receiving a CLOSE signal, the hand closes in the last used gripping mode (lateral or opposition).

#### **Example 3:**

When a low electrode signal is applied, the system builds up the lowest gripping force in lateral mode or opposition mode to grip an object.

#### **Example 4:**

When there is a higher electrode signal, a greater gripping force is generated; through secondary gripping with a stronger electrode signal, the gripping force is increased up to the maximum gripping force.

OPEN	CLOSE
Sustained electrode signal	Sustained electrode signal
	<b>Gripping force:</b> Proportional
<b>Speed:</b> Proportional	<b>Speed:</b> Proportional

### 8.2 2-ELE

Proportional control with 2 electrodes.

The open position of the Michelangelo Hand is the initial position for lateral mode, opposition mode and the neutral position. **It is only possible to exit the neutral position by fully opening the Michelangelo Hand.**

The gripping speed and gripping force are determined by the strength of the electrode signal (resulting from the muscle contraction). The gripping force can be increased up to the maximum gripping force at any time through a stronger electrode signal.

Only when the Michelangelo Hand is in the fully opened hand position does it change to the neutral position if the hand is not controlled via the electrode signals. While closing to the neutral position, the system builds up the minimum gripping force (approx. 15 N).

#### **Opening:**

Proportionally controlled by the OPENING electrode (e.g. input 2)

#### **Closing:**

Proportionally controlled by the CLOSING electrode (e.g. input 1)

**Example 1:**

When a low electrode signal is applied, the system builds up the lowest gripping force to grip an object.

**Example 2:**

When a higher electrode signal is applied, a greater gripping force is generated. Secondary gripping with a stronger electrode signal increases the gripping force up to the maximum.

**Example 3:**

The hand is fully open. After the OPEN muscle has relaxed, the hand slowly moves to the neutral position. An OPEN signal allows the hand to move out of the neutral position again.

OPEN	CLOSE
Sustained electrode signal	Sustained electrode signal
	Only possible from the fully opened position
	<b>Gripping force:</b> Proportional
<b>Speed:</b> Proportional	<b>Speed:</b> Proportional

**8.3 1-Double**

Proportional control with 1 electrode.

The open position of the Michelangelo Hand is the initial position for lateral mode, opposition mode and the neutral position. **It is only possible to exit the neutral position by fully opening the Michelangelo Hand.**

The Michelangelo Hand is opened with a quick and strong signal and closed with a slow, gentle signal.

The gripping speed and gripping force are determined by the strength of the electrode signal (resulting from the muscle contraction). The gripping force can be increased up to the maximum gripping force at any time through a stronger electrode signal.

Only when the Michelangelo Hand is in the fully opened hand position does it change to the neutral position if the hand is not controlled via the electrode signals. While closing to the neutral position, the system builds up the minimum gripping force (approx. 15 N).

**Opening:**

At a proportional speed by means of a fast and strong muscle signal.

**Closing:**

At a proportional speed by means of a slow and gentle muscle signal.

**Example 1:**

With a fast and strong electrode signal, the hand opens at a proportional speed.


**Example 2:**

The hand is fully open. After the OPEN muscle has relaxed, the hand slowly moves to the neutral position. An OPEN signal allows the hand to move out of the neutral position again.

**Example 3:**

With a slow and gentle electrode signal, the hand closes at a proportional speed and gripping force.

OPEN	CLOSE
Quick, strong electrode signal that reaches the upper threshold within 80 ms and is maintained above the upper threshold for at least 40 ms.	Slow, gentle electrode signal that does not reach the upper threshold within 80 ms.
	<b>Gripping force:</b> Proportional

### 8.4 1-Vario

Proportional control with 1 electrode.

The open position of the Michelangelo Hand is the initial position for lateral mode, opposition mode and the neutral position. **It is only possible to exit the neutral position by fully opening the Michelangelo Hand.**

Only when the Michelangelo Hand is in the fully opened hand position does it change to the neutral position if the hand is not controlled via the electrode signals. While closing to the neutral position, the system builds up the minimum gripping force (approx. 15 N).

#### Opening:

Proportional. The OPENING speed is determined by the level and speed of the muscle contraction.

#### Closing:

Proportional. The CLOSING speed and gripping force is determined by the speed of the muscle relaxation.

#### Stopping:

Through very slow muscle relaxation via the electrode. The Michelangelo Hand remains in the current position.

#### Example 1:

After opening, the user relaxes the muscle at low speed. Closing takes place proportionally to the decline of the electrode signal. The object is grasped with low force.

#### Example 2:

After opening, the user relaxes the muscle at the fastest speed. The Michelangelo Hand closes at the fastest speed and grasps using the maximum gripping force.

OPEN	CLOSE
Increasing electrode signal through muscle contraction	Declining electrode signal through muscle relaxation
<b>Speed:</b> Proportional	<b>Gripping force:</b> Proportional to the decline of the electrode signal
<b>Speed:</b> Proportional	<b>Speed:</b> Proportional

## 9 Switching Methods

### 9.1 Co-Contraction

**Control:** 2 electrodes

**Indication:** For patients with co-contraction ability

**Execution:** Both electrode signals must exceed the co-contraction thresholds within 80 ms and then decline. This time can be configured using the adjustment software.

**Short co-contraction:** The signal declines quickly (default within 500 ms)

**Long co-contraction:** The signal is held longer (default more than 500 ms) and then declines

## 9.2 Four-Channel Control

**Control:** 2 electrodes

**Indication:** For patients with 2 fast signals

**Execution, slowly rising signal:** By slowly beginning muscle contraction. It takes more than 80 ms for the electrode signal to reach the set upper switching threshold after exceeding the ON-threshold.

**Execution, quickly rising signal:** By quickly and strongly beginning muscle contraction. The electrode signal reaches the set upper switching threshold within 80 ms after exceeding the ON-threshold.

## 9.3 Pulse

**Control:** Electrode or switch for input 1 and/or input 2

**Execution:** The electrode signal or switch signal must exceed the set pulse threshold within 80 ms and then decline.

## 9.4 Long Open

**Control:** 1 or 2 electrodes

**Indication:** For patients who can maintain the open hand position with a long muscle signal.

**Execution:** The electrode signal must be held at the value for the upper threshold of the operating range (default 100 or HI) for the duration of the switching time when the hand is already fully open.

The duration of the switching time can be configured in the adjustment software.


## 9.5 Automatic switching after neutral position is reached

**Control:** By means of the neutral position

**Indication:** For patients who have a preferred type of grip.

**Execution:** Switching is triggered after the neutral position has been reached and the delay time has passed. The duration of the delay time can be configured in the adjustment software.

## 10 Open the Axon-Bus Gripping Component in an Emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open. A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

## 11 Disposal

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

## 12 Legal Information

### 12.1 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 12.2 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

### 12.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

This product meets the requirements of the European Directive 1999/5/EC for radio equipment and telecommunications terminal equipment. The conformity assessment was drawn up by the manufacturer in accordance with Annex IV of the directive.

### 12.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

—Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution: Exposure to Radio Frequency Radiation.**

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP  
3820 West Great Lakes Drive  
Salt Lake City, Utah 84120-7205 USA  
Phone + 1-801-956-2400  
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road  
L7L 5N5 Burlington, Ontario  
Canada  
Phone + 1-800-665-3327

---

**Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.**

## 13 Appendices

### 13.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the European Directive for Medical Devices 93/42/EEC and Radio and Telecommunications Terminal Equipment 1999/5/EC

LOT:PPPP YYYY WW Lot number

Legal manufacturer

Compliance with the requirements according to "FCC Part 15" (USA)

Non-ionising radiation

Compliance with the requirements according to the "Radiocommunications Act" (AUS)

### 13.2 Technical data

Environmental conditions	
Storage and transport in original packaging	-20 °C/-4 °F to +40 °C/+104 °F
Storage and transport without packaging	-20 °C/-4 °F to +40 °C/+104 °F Max. 80% relative humidity, non-condensing
Operation	0 °C/+32 °F to +60 °C/+104 °F Max. 80% relative humidity, non-condensing

General Information	
Reference number	13E500
Dimensions	53 x 25 x 9 mm/2.08 x 0.98 x 0.35 inch
Weight	15 g/0.53 oz
Operating voltage	11.1 V $\equiv$
Power supply	757B500 AxonEnergy Integral or 757B501 AxonEnergy Integral
Service life	5 years

Data transfer	
Wireless technology	Bluetooth 2.0
Range	approx. 5 m/16.4 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK
Data rate (over the air)	704 kbps
Maximum output power (EIRP)	-9.4 dBm (0.116 mW)

### 13.3 Glossary

The term "Axon" stands for **A**daptive **e**xchange **o**f **n**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

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### **Caution: Exposure to Radio Frequency Radiation.**

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Responsible party:  
Otto Bock Health Care, LP  
3820 West Great Lakes Drive  
Salt Lake City, Utah 84120-7205 USA  
Phone + 1-801-956-2400  
Fax + 1-801-956-2401

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<http://www.hc-sc.gc.ca/rpb>.

Responsible party:  
Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road  
L7L 5N5 Burlington, Ontario  
Canada  
Phone + 1-800-665-3327

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