















757B35=0

DE	Gebrauchsanweisung (Fachpersonal)
EN	Instructions for Use (Qualified Personnel)
FR	Instructions d'utilisation (Personnel spécialisé)
IT	Istruzioni per l'uso (Personale tecnico specializzato)
ES	Instrucciones de uso (Personal técnico especializado)
PT	Manual de utilização (Pessoal técnico)
NL	Gebruiksaanwijzing (Vakmensen)
SV	Bruksanvisning (Fackpersonal)
DA	Brugsanvisning (Faguddannet personale)
NO	Bruksanvisning (Fagpersonell)
TR	Kullanma talimatı (Uzman personel)
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό)
RU	Руководство по применению (Квалифицированный персонал)

1 Vorwort

INFORMATION

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INFORMATION

Last update: 2014-01-14

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions and the precautions specified in this accompanying document.
- ▶ If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

2 Product Description

2.1 Function

The 757B35=0 MyoEnergy Integral is a complex power supply system. It consists of a battery, charging receptacle, 2-pole supply cable, 3-pole communication cable and electronics. The components are inseparable and ready for mounting. Use the 757L35 MyoCharge Integral to charge the battery.

2.2 Components

2.2.1 Rechargeable Battery

The battery (Fig. 1, item 1) consists of two cells.

2.2.2 Charging Receptacle

The charging receptacle (Fig. 1, item 2) with integrated button, LED, and beeper has several functions.

- Contacts for battery charging
- LED for displaying the current battery charge level
- Button for turning the prosthesis on and off and for activating the battery charge level display
- Beeper for providing feedback on operating states

2.2.3 Supply Cable

The 2-pole supply cable (Fig. 1, item 3) with receptacle establishes the connection between the battery and the respective prosthesis component. The 757VAK1 Battery Adapter Cable for VASI System is required for operation with VASI components. This adapter cable is not included in the scope of delivery and can be ordered from Ottobock.

2.2.4 Communication Cable

The 3-pole communication cable (Fig. 1, item 4) with receptacle is used for the exchange of data. For functionality, please refer to the instructions for use of the respective prosthetic components.

2.2.5 Electronics

The integrated electronics (Fig. 1, item 5) protect the battery against short circuits, overvoltage, undervoltage and charging outside the allowable temperature range. They also communicate with the battery charger or the MyoLino Link.

3 Application

3.1 Intended Use

The 757B35=0 MyoEnergy Integral is to be used exclusively as a power supply for the MyoBock prosthetic system and VASI prosthetic components. The 757B35=0 MyoEnergy Integral is to be used for unilateral or bilateral exoprosthetic fittings of the upper limbs.

3.2 Area of Application

The 757B35=0 MyoEnergy Integral can be used for upper arm and forearm prosthetic components.

3.3 Conditions of Use

Using the 757B35=0 MyoEnergy Integral with the 4in1 Controller is prohibited.

3.4 Required Qualifications

The fitting of a patient with the 757B35=0 MyoEnergy Integral must be carried out exclusively by prosthetists.

- ▶ **INFORMATION: The 757B35=0 MyoEnergy Integral will be referred to as the MyoEnergy Integral in the remainder of this text.**

4 Safety

4.1 Explanation of Warning Symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

5 Safety Instructions

5.1 Product-Specific Safety Instructions

⚠ CAUTION

Non-observance of safety instructions

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- ▶ Destruction of the product.
- Follow the safety instructions and the precautions specified in this accompanying document.

⚠ CAUTION

Disconnecting or establishing electrical connections when switched on

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Before replacing prosthetic components (e.g. gripping components), be sure to switch off the product as specified in the instructions for use.

⚠ CAUTION

Penetration of liquids

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Ensure that no liquid penetrates into the battery.
- Do not use the battery for bathing prostheses.

⚠ CAUTION

Improper use of the product

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Instruct the patient in the proper handling of the product according to the section "Patient-Specific Safety Instructions".

INFORMATION

Prior to initial operation of the product, charge the supplied battery.

5.2 Patient-Specific Safety Instructions

WARNING

Simultaneous use of the product with active, implanted systems

- ▶ Injuries due to temporary influences of active implantable systems (e.g. pacemakers, defibrillators, etc.) due to electromagnetic interference of the product.
- Do not bring the product into the immediate proximity of active, implantable systems.

CAUTION

Penetration of dirt and humidity

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Ensure that neither solid particles nor liquids can penetrate into the product.

CAUTION

Mechanical overload caused by mechanical influences or external loads

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Do not subject the product to mechanical vibrations or impacts.
- Check the product for visible damage before each use.

CAUTION

Wearing in extreme ambient temperatures

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Avoid wearing the prosthesis in areas with temperatures outside of the permissible range. See the section "Technical Data" in the appendix to these instructions for use.

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.99 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
 - DECT cordless phones incl. base station: 0.35 m
 - WiFi (routers, access points,...): 0.22 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

CAUTION

Electromagnetic interference

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Avoid areas where there are sources of strong electromagnetic interference (e.g. EAS systems in department stores).

⚠ CAUTION**Independent manipulation of the product**

- ▶ Injury due to faulty control or malfunction of the prosthetic system.
- Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- Do not damage the battery or separate the connection cables between the battery packs.
- The product and any damaged components may only be opened and/or repaired by certified Ottobock Myo-Service technicians.

⚠ CAUTION**Dirtying of the electrical contacts**

- ▶ Risk of injuries due to faulty control or malfunction of the prosthetic system caused by an insufficient supply of power to the prosthetic components.
- Ensure that the contacts of the charging receptacle are always clean and free of grease.
- Clean the electrical contacts of the battery charger and charging receptacle regularly using cotton swabs and a mild soap solution.
- Avoid damaging the contact surface coating with pointed or sharp objects.

6 Delivery/Storage

6.1 Scope of Delivery

Item shown in Fig. 2	Quantity	Product description
Item 1	1	757B35=0 MyoEnergy Integral
Item 2	2	Battery lamination dummy
Item 3	1	Electronics lamination dummy
Item 4	1	Charging receptacle lamination dummy
Item 5	1	Drilling template for charging receptacle
Not illustrated	1	647G885 Instructions for Use

6.2 Storage

Please see the specifications in the "Technical Data" section.

7 Recommended Laminating Instructions for the MyoEnergy Integral

The product can be positioned in the upper arm or forearm socket. The battery is secured between the outer and inner socket according to space availability.

INFORMATION

The length of the flexible flat cable between the electronics and the charging receptacle can be modified by carefully bending or folding it. Only light pressure must be applied when pinching and sliding down the folded edge. If the flexible flat cable is bent in the opposite direction at a spot where it has already been folded, serious damage can occur.

7.1 Laminating Instructions for the MyoEnergy Integral in the Upper Arm Socket

Under consideration of the physiological patient indications, as little weight as possible should be placed in the forearm, so the most proximal positioning of the battery is desirable. Therefore, use

the 743A23 Alignment Aid for ErgoArm Elbow Components for the socket alignment of a 12K44=* or 12K50=* ErgoArm. Use adhesive tape to fasten the battery to the bars of the alignment aid. Please observe the information in the 647H447 Instructions for Use, which are included with the 743A23 Alignment Aid for ErgoArm Elbow Components.

7.2 Laminating Instructions for the MyoEnergy Integral in the Forearm Socket

NOTICE

Gluing the charging receptacle with an unsuitable adhesive

- ▶ Damage to the outer socket during removal of the charging receptacle due to the use of unsuitable adhesives (e.g. Cyamet quick-drying adhesive).
- 617H46 Bonding Agent must be used to glue the charging receptacle because the battery and charging receptacle must be able to be removed from the outer socket in case of a defect.

INFORMATION

- Before handling and processing the 617H12=* Pedilen Rigid Foam 200 and the 617P21=* Pedilen Hardener, please read the corresponding material safety data sheet.
- Before handling and processing the 617H55=* C-Orthocryl Lamination Resin or the 617H19=* Orthocryl Lamination Resin and the 617P37=* Hardening Powder, please read the corresponding material safety data sheet.
- Ensure sufficient ventilation of rooms and workstations.
- For your own protection, wear an appropriate respirator, protective gloves and goggles when sanding the socket.

INFORMATION

Note that for the following lamination technique, the lamination dummies for the battery and charging receptacle may not be fixed to the inner socket with double-sided tape, screws or similar. They may only be held in place by the 99B25 Tube Sock. Modelling with Plasticine is, therefore, not required.

- 1) Fabricate the inner socket. Fix the lamination dummies for the charging receptacle, batteries, electronics and electrodes to the inner socket. Cover the inner socket with standard plastic wrap.
- 2) Pull the 99B25 Tube Sock down so it is free of wrinkles and tie off proximally. Verify the position of the lamination dummies. Warm the 99B71=* PVC Bag and pull it over the residual limb. Tie the bag off and evacuate (see Fig. 4).
- 3) Fabricate the foam positive using 617H12=* Pedilen Rigid Foam 200 and 617P21=* Pedilen Hardener using known methods (see Fig. 5).
- 4) Sand the foam positive (see Fig. 9).
- 5) After sanding, cover the foam positive with PVA foil and sufficient 623T3=* Perlon Stockinette and tie off.
- 6) In order to increase the stability of the socket, the laminate can be strengthened with a layer of 616G12 Carbon Fibre Cloth between the first and second layers of 623T3=* Perlon Stockinette. Ensure that the carbon fibre cloth does not overlap the lamination ring for reasons of electrical safety.
- 7) Laminate the outer socket using known methods and allow it to harden.
- 8) Place the two bore holes of the charging receptacle template on the nubs of the outer socket and mark the contour on the outer socket using a thin felt-tip pen or scriber (see Fig. 10).
- 9) Burn the nubs out from inside. The resulting holes indicate the centre of the bore holes to be drilled (see Fig. 11).

- 10) Drill both holes that have been burned out up to the marked contour using a 726W9=* Tapered Cutter (see Fig. 12).
- 11) Cut the contour free. Check the contour using the charging receptacle template. The template should fit the cutout without play (see Fig. 13).
- 12) Perform a final check of the cutout using the charging receptacle. For this purpose, push the charging receptacle out through the cutout from the inside of the socket (see Fig. 14).
- 13) Insert the charging receptacle into the cutout. The charging receptacle must fit flush to the outside of the socket without pressure being applied (see Fig. 15).
- 14) If the outer socket has a wall thickness of more than 2 mm, it will not be possible to push the charging receptacle in. The inside of the socket must be sanded in the area of the flat ribbon cable of the charging receptacle (see Fig. 16).

INFORMATION: For demonstration purposes, the socket is displayed cut open.

- 15) Example of a charging receptacle that has been positioned correctly in the outer socket (see Fig. 17).

INFORMATION: For demonstration purposes, the socket is displayed cut open.

- 16) Apply 617H46 Bonding Agent to the small rim of the charging receptacle. The bonding agent will make the charging receptacle stick to the outer socket (see Fig. 18).

- 17) Insert the charging receptacle into the outer socket. Fix the charging receptacle to the socket with adhesive tape and let the bonding agent cure (see Fig. 19).

INFORMATION: Please observe the cure time of the bonding agent.

- 18) Fix the batteries and electronics to the inside of the outer socket (e.g. with double-sided tape) and insert the inner socket. The inner socket will fix the batteries and electronics in the trough of the outer socket (see Fig. 20).

INFORMATION: For demonstration purposes, the socket is displayed cut open.

8 Use/Operation

INFORMATION

Before switching on the MyoEnergy Integral for the first time, fully charge it using the 757L35 MyoCharge Integral.

8.1 Charging the Battery

Connect the charging plug to the charging receptacle. The beeper will emit one long tone. The prosthesis is turned off and charging starts automatically.

INFORMATION

If you want to use the prosthesis, you have to remove the charging plug and turn the prosthesis on again.

8.2 Switching On and Off

The prosthesis is switched on or off by pushing the button on the charging receptacle (see Fig. 3, arrow). This requires pressing the button for more than 1 second.

Two short beeps are heard and the LED is briefly illuminated when the system powers up. Different colours indicate the current battery charge level. See also the section "Beeper Functions".

Switching off is confirmed by 1 long beep. The LED does not light up.

8.3 Query Battery Charge Level

The battery charge level can be displayed at any time whilst the prosthesis is turned on. In order to do so, press the charging receptacle button (see Fig. 3, arrow) for less than 1 second. The LED is illuminated. The colour of the LED (see table) indicates the current battery charge level.

Battery empty
 Battery 50% charged
 Battery fully charged

Red LED is illuminated
 Yellow LED is illuminated
 Green LED is illuminated

8.4 Safety Shutoff

The purpose of the safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature during charging
- Short circuit
- Overvoltage or undervoltage

INFORMATION

After a safety shutdown, the battery charger must be connected to reactivate the electronics.

8.5 Beeper Functions

The following beeper functions are set by default:

	One long beep	Two short beeps
Switch on the prosthesis: Press button for more than 1 second.		•
Switch off the prosthesis: Press button for more than 1 second.	•	
Connect charging plug: Charging begins.	•	
Disconnect charging plug: Charging ends.	•	

The following beeper functions exist only in conjunction with the 9E420=* 7in1 Controller:

	One long beep	Two short beeps	Continuous pulses (as long as the button is being pressed)
Parental access switch			•
Connect the MyoLino Link		•	
Disconnect the MyoLino Link: Prosthesis was previously switched on		•	
Disconnect the MyoLino Link: Prosthesis was previously switched off	•		

9 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Air dry to remove residual moisture.

10 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that does not comply with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

11 Legal Information

11.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

11.2 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11.3 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

12 Appendices

12.1 Symbols on the Battery

Symbol	Meaning
	Declaration of conformity according to the European Directive for Medical Devices 93/42/EEC.
	This product must not be disposed of with regular domestic waste. Disposal that does not comply with the regulations of your country may have a detrimental impact on health and the environment. Please note the information provided by the responsible authorities in your country regarding return and collection processes.
LOT:PPPP YYYY WW	Lot number
	Legal manufacturer

12.2 Accessories

757VAK1 Battery Adapter Cable for the VASI System

12.3 Technical Data

12.3.1 Specific Technical Data

Designation	Value
Capacity	300 mAh
Type	Lithium polymer
Dimensions	approx. 35 x 20 x 20 mm
Weight	approx. 25 g (battery pack without accessories)

12.3.2 General Technical Data

Designation	Value
Output voltage	approx. 7.4 V
Charging voltage	approx. 8.2 V
Charging time	approx. 2.0 h
Operating temperature	0 °C to 60 °C (32 °F to 140 °F)
Charging temperature	0 °C to 45 °C (32 °F to 113 °F)
Storage and transport temperature (with and without packaging)	-20 °C to 50 °C (-4 °F to 122 °F)
Operating time of prosthesis	approx. 1000 grip cycles with MyoBock Adult System approx. 2500 grip cycles with MyoBock Child System
Service life of the MyoEnergy Integral	2 years
Battery service life	2 years

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