# 12K110N=50-S

<b>DE</b> Gebrauchsanweisung .....	2
<b>EN</b> Instructions for use .....	11
<b>FR</b> Instructions d'utilisation .....	20
<b>IT</b> Istruzioni per l'uso .....	29
<b>ES</b> Instrucciones de uso.....	38
<b>PT</b> Manual de utilização .....	47
<b>NL</b> Gebruiksaanwijzing .....	56
<b>SV</b> Bruksanvisning.....	65
<b>DA</b> Brugsanvisning .....	74
<b>NO</b> Bruksanvisning.....	83
<b>EL</b> Οδηγίες χρήσης .....	92

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English

Date of last update: 2021-04-15

- Please read this document carefully before using the product and observe the safety notices.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

**Before fitting the patient with the 12K110N=50-S DynamicArm Plus , please first read the 647G526=all\_INT Instructions for Use !**

## 1 Intended use

### 1.1 Medical purpose

The 12K110N=50-S DynamicArm Plus is to be used **solely** for the duration of the service and for trial fittings for exoprosthetic fittings of the upper limbs.

### 1.2 Application

The 12K110N=50-S DynamicArm Plus is suitable for unilateral and bilateral above-elbow amputees. The 12K110N=50-S DynamicArm Plus is to be used **exclusively** for the fitting of adults.

### 1.3 Application conditions

The 12K110N=50-S DynamicArm Plus has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

### 1.4 Prosthetist qualification

The fitting of a patient with a 12K110N=50-S DynamicArm Plus may only be carried out by prosthetists who have been authorized by Ottobock based on completion of a corresponding training course.

## 2 Safety instructions

Please learn the safety instructions as well as patient information from the 647G526=all\_INT Instructions for Use for the 12K110N DynamicArm Plus .

Failure to follow the safety instructions mentioned below can lead to faulty control or malfunction of the 12K110N=50-S DynamicArm Plus and result in risk of injuries for the patient.


## 3 Serviceinformation

The 12K110N=50-S DynamicArm Plus must be sent for examination to the authorized Ottobock Service after every patient fitting. In the service the complete DynamicArm will be examined and possibly readjusted and any wearing parts will be exchanged if necessary.

## 4 Adjustment, Fittings, Application

Before fitting the patient with the 12K110N=50-S DynamicArm Plus , please first read the 647G526=all\_INT Instructions for Use for the 12K100N DynamicArm . They include important additional information about the application of the 12K110N=50-S DynamicArm Plus .

## 5 Scope

### 5.1 System components

- Item 1 Adapter long
- Item 2 Adapter short
- Item 3 Threaded tube, short, with lamination ring and counter nut
- Item 4 Threaded tube, long, with lamination ring and counter nut
- Item 5 9E169 Coaxial Plug
- Item 6 10S17 Electric Wrist Rotator
- Item 7 Retaining ring
- Item 8 12K110N=50-S DynamicArm Plus Service
- Item 9 Disposable foam cover
- Item 10 Red assembly clamp

### 5.2 Accessories and tools

The following items have been supplied additionally to the system components:

- Item 11 Hook wrench
- Item 12 Strap wrench
- Item 13 Suitcase key
- Item 14 Euro mains cable
- Item 15 US mains cable
- Item 16 Battery charger
- Item 17 Red assembly clamp, band clamp

### 5.3 Manuals

All manuals necessary for the fitting can be found on the inner side of the suitcase under the foam cover.

- no illustr. Instructions for Use for DynamicArm
- no illustr. Alignment and Assembly Instructions for DynamicArm
- no illustr. Note tag for reusable package
- no illustr. Instructions for Use for Electric Wrist Rotator
- no illustr. Instructions for Use for Battery Charger
- no illustr. Patient Information for DynamicArm

## 6 Disassemble the patient's prosthesis

Before adapting the 12K110N=50-S DynamicArm Plus to the patient, the patient's elbow harness must be removed from the socket.

**Note !** The 12K110N=50-S DynamicArm Plus is compatible with upper arm sockets of the 12K110N DynamicArm Plus.

Example: Prosthetic socket with 12K100N DynamicArm and system electric hand.

Rotate the prosthetic component (system electric hand or system electric greifer) until the quick disconnect wrist unlocks. Carefully pull off the prosthetic components.

**Caution!** The DynamicArm must be switched off before removing the prosthetic components.

Turn the finger wheel counter-clockwise to set minimum compensation. This will prevent the band from being pulled in.

**Note:**

A DynamicArm with a retracted band must only be repaired by an authorized Ottobock Myo-Service provider!

Extend the elbow and lock it with the Bowden cable. Switch off the DynamicArm .

Remove the band clamp.

**Attention:** The lifter band should not be removed from the elbow ball. Fixing the lifter strap in place (e.g. with a haemostat) is no longer necessary.


Lift the band clamp out.


Insert red assembly clamp.

Tighten the Allen head screw.

Turn the friction screw until the socket can be screwed off the elbow joint.  
Unscrew the socket from the DynamicArm .

Determine the length. Measure the distance between the front edge of the socket and the rotation center of the joint. Document value 'X'.

## 7 Length ranges

### 7.1 Length range 189 mm – 211 mm

Required accessories (see also section 5):

- Coaxial plug (item 5)
- Threaded tube, short, with lamination ring and counter nut (item 3)
- DynamicArm (item 8)

### 7.2 Length range 214 mm – 256 mm

Required accessories (see also section 5):

- Retaining ring (item 7)
- Electric wrist rotator (item 6) / coaxial plug (item 5)
- Adapter short (item 2)
- Threaded tube, long, with lamination ring and counter nut (item 4)
- DynamicArm (item 8)

### 7.3 Length range 259 mm – 305 mm

Required accessories (see also section 5):

- Retaining ring (item 7)
- Electric wrist rotator (item 6) / coaxial plug (item 5)
- Adapter long (item 1)
- Threaded tube, long, with lamination ring and counter nut (item 4)
- DynamicArm (item 8)



## 8 Length adjustment of the DynamicArm

1. On the basis of the measured length, you assemble the corresponding version.

2. Screw the threaded tube in or out as required and adjust the measurement 'X' (see section 6).

3.  **Attention:** The threaded tube must be screwed in at least to the point where the marker bore holes (arrow) are covered by the counter nut. If this is not possible, the next longer version will have to be chosen.

4. Tighten the adapter and counter nut of the threaded tube clockwise (hand-tight).

5. Clamp the correctly positioned strap wrench into the vise. Fix the DynamicArm in the strap wrench.

First apply the hook wrench to the adapter (optional) and finally to the counter nut of the threaded tube. Pull the nut tight (clockwise) until the torque moment is reached.

- 5.1 The torque moment of 16 Nm will be reached as soon as the handles touch each other at the hook wrench.

## 9 Completion of the DynamicArm

Slip the disposable foam cover over the DynamicArm .

Mark the corresponding length according to the determined measurement 'X' (see section 6, item 7).


Use scissors to cut-off the disposable foam cover.

Connect the cables with the electric wrist rotator / coaxial plug.

Next mount the electric rotator / coaxial plug into the lamination ring of the threaded tube. See the 647G526=all\_INT Instructions for Use.

Pull the disposable foam cover previously cut to length over the DynamicArm .

Mount the prosthetic components.

Pull the cosmetic glove over the system electric hand.

In case of a fitting with the electric wrist rotator, make sure that the rotational movement is working smoothly.

## 10 Assemble the DynamicArm

After adaptation of the 12K110N=50-S DynamicArm Plus , you can fit the prosthesis with the prosthetic socket.

1. Dismantle the lamination ring of the 12K100N=50-S as well as the prosthetic socket of the 12K100N DynamicArm (see section 6, items 4 – 7).
2. Screw the socket onto the elbow joint. Or turn the friction screw until the socket can be screwed on the elbow joint (see 647G526=all\_INT Instructions for Use).
3. Remove the Allen head screw. Then replace the red assembly clamp with the band clamp and reinstall the Allen head screws (see 647G526=all\_INT Instructions for Use).
4. Turn the friction screw to set the friction of humeral rotation (see 647G526=all\_INT Instructions for Use).
5. Adjust compensation to the patient by turning the finger wheel (see 647G526=all\_INT Instructions for Use).

## 11 Technical data

Power supply	Li-Ion battery
Charging temperature:	> 0 °C
Weight (dependent on forearm length)	approx. 1000 g
Maximum lifting force:	50 N
Maximum load	approx. 60 N
Flexion angle	approx. 15° to 145°
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F, Max. 93% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F, Max. 93% relative humidity, non-condensing
Charging temperature	+5 °C/+41 °F to +40 °C/+104 °F
Operation	+5 °C/+41 °F to +45 °C/+113 °F, Max. 93% relative humidity, non-condensing

### 11.1 Disposal

This product may not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

## 12 Liability

Otto Bock Healthcare Products GmbH, hereinafter called manufacturer, assumes liability only if the user observes the product's specified fabrication and use instructions, adheres to its care instructions, and follows its service schedule. This product may only be used in combination with components that have been authorized by the manufacturer (please refer to the instructions for use and catalogs). The manufacturer assumes no liability for damage caused by component combinations (use of third-party products) and applications which it did not authorise.

This product may only be opened and repaired by authorized Ottobock staff.

## 13 Trademarks

All product names mentioned in this accompanying document are subject without restriction to the respective applicable trademark laws and the property of the respective owners. All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 14 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the 2014/53/EU directive.

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### FCC-Statement:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

### Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care LP  
Two Carlson Parkway North, Suite 100  
55447 Minneapolis, Minnesota, USA  
Phone + 1-763-553-9464

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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### Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website <http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road  
L7L 5N5 Burlington, Ontario  
Canada

Otto Bock Healthcare Products GmbH  
Brehmstraße 16 · 1110 Wien · Austria  
T +43-1 523 37 86 · F +43-1 523 22 64  
[info.austria@ottobock.com](mailto:info.austria@ottobock.com) · [www.ottobock.com](http://www.ottobock.com)

Ottobock has a certified Quality Management System in accordance with ISO 13485.