## DynamicArm 12K100N

## DynamicArm Plus 12K110N

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	2
<b>EN</b> Instructions for use (user) .....	14
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	26
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	39
<b>ES</b> Instrucciones de uso (Usuario) .....	51
<b>PT</b> Manual de utilização (Usuário) .....	63
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	75
<b>SV</b> Bruksanvisning (Användare) .....	87
<b>NO</b> Bruksanvisning (Bruker) .....	99
<b>PL</b> Instrukcja użytkowania (Użytkownik) .....	110
<b>CS</b> Návod k použití (Uživatel) .....	122
<b>EL</b> Οδηγίες χρήσης (Χρήστης) .....	133
<b>RU</b> Руководство по применению (Пользователь) .....	146

## INFORMATION

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### INFORMATION

Date of last update: 2021-04-01

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.


The product "12K100N=\* DynamicArm and 12K110N=\* DynamicArm Plus" is referred to simply as the product/elbow joint below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

### 2.2 Function

The product is a myoelectrically controlled elbow joint driven by an electric motor. It helps the wearer perform everyday tasks in combination with other prosthetic components.

A Li-Ion battery integrated into the product provides the energy required. The "ElbowSoft/Elbow-Soft TMR" software is used for the patient-specific adjustment of the product.

If the elbow joint is off or the battery is depleted, the elbow joint can still be locked and unlocked in any position, even when loaded, by means of the pull cable.

The product can lift a load of up to 5 kg at a forearm lever length of 305 mm. Once this load is exceeded, the joint will be locked. Flexion and extension are only possible again after the load is reduced.

### 2.2.1 Definition of terms

#### AFB (Automatic Forearm Balance)

AFB is a flexion assist in the form of a mechanical gear mechanism in the forearm of the elbow joint. The Automatic Forearm Balance provides for harmonious movement control, free swing and reduction of energy consumption. It stores the energy that is released when the arm is extended and uses this to assist with flexion. The degree of flexion assistance can be adjusted to the individual weight of the prosthetic forearm and to different clothing by means of a dial.

The free swing of the forearm corresponds to the natural movement behaviour of the upper limb while walking. After extension of the forearm, the vario-gear automatically decouples the forearm from the drive unit so it can swing freely. The AFB slightly dampens the free swing, thereby making a close approximation of the forearm's physiological movement behaviour possible.

The product does not consume electrical energy during the free swing. A short muscle signal to flex the product ends the free swing.

#### Humeral rotation feature

The humeral rotation feature is an upper arm rotation joint with bilateral stops. It is used for rotation of the forearm. The resistance of the humeral rotation feature can be changed using the adjusting screw.

#### Lock

Releasing and locking the elbow joint is automated according to operating situations, the applied load and the control signal generated by the patient by means of the electronically controlled locking drive.

The product has a load capacity of up to 230 N in the locked state with a forearm lever length of 305 mm. The lock will slip under higher loads.

#### Mechanical release

With the product switched off or when the battery is drained, the elbow can be manually unlocked and locked again with a slight pull on the unlock cable. In this way, the forearm can be brought into the desired position. Mechanical release can also take place under load.

#### On/off switch

Press the on/off switch to operate it (arrow).

A signal sounds and the product briefly vibrates once when it is turned on. Turning the elbow joint off during longer passive breaks is recommended. This increases the period of use before charging the battery.

## 3 Intended use

### 3.1 Indications for use

The product is to be used **exclusively** for exoprosthetic fittings of the upper limbs and serves as an anatomical replacement for the elbow joint and forearm.

### 3.2 Conditions of use

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment. Permissible ambient conditions are described in the technical data (see page 24).

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The fitting of a user with the product may only be carried out by O&P professionals who have been authorised through the corresponding Ottobock training.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

<p> <b>WARNING</b></p> <p><b>The heading describes the source and/or the type of hazard</b></p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none"> <li>&gt; E.g.: Consequence 1 in the event of failure to observe the hazard</li> <li>&gt; E.g.: Consequence 2 in the event of failure to observe the hazard</li> <li>▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.</li> </ul>
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### 4.3 General safety instructions

<p> <b>WARNING</b></p> <p><b>Non-observance of safety notices</b></p> <p>Personal injury/damage to the product due to using the product in certain situations.</p> <ul style="list-style-type: none"> <li>▶ Observe the safety notices and the stated precautions in this accompanying document.</li> </ul>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p> <b>WARNING</b></p> <p><b>Use of the product while operating a vehicle or machinery</b></p> <ul style="list-style-type: none"> <li>&gt; Risk of accident due to unexpected behaviour of the product.</li> <li>&gt; Risk of injury due to faulty control or malfunction of the product.</li> <li>▶ Applicable national legal regulations must be complied with when operating vehicles of any kind or when operating machinery with the product.</li> <li>▶ Have your driving ability examined and verified by an authorised test centre.</li> <li>▶ Switch off the product before driving.</li> </ul>
--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**⚠ WARNING**

**Charging the prosthesis without taking it off**

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

**⚠ WARNING**

**Operating the prosthesis near active implanted systems**

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**⚠ CAUTION**

**Overloading due to unusual activities**

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Independent user manipulation of system components**

Injury due to faulty control or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

**⚠ CAUTION**

**Signs of wear and tear on the product components**

Injury due to faulty control or malfunction of the product

- ▶ Observing the specified service intervals is mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty (see the section "Maintenance").

**⚠ CAUTION**

**Mechanical stress on the prosthesis**

Injury due to faulty control or malfunction of the prosthesis.

- ▶ Do not subject the prosthesis to mechanical vibrations or impacts.
- ▶ Check the prosthesis for visible damage before each use.

**⚠ CAUTION**

**Penetration of dirt and humidity into the product**

Risk of injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.
- ▶ Do not expose the product, and especially the elbow joint, to dripping or splashing water.
- ▶ Wear the product, and especially the elbow joint, under suitably resistant clothing in the rain.

**⚠ CAUTION**

**Changing prosthetic components when switched on**

Injury due to faulty control or malfunction of the prosthetic components.

- ▶ Before changing prosthetic components (e.g. terminal device) switch the product off.

**⚠ CAUTION**

**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

**⚠ CAUTION**

**Incorrect donning and doffing of the prosthesis**

Risk of injury due to sudden flexion of the prosthesis.

- ▶ Switch off the prosthesis before donning or doffing.
- ▶ The prosthesis must always be flexed for donning and doffing.

**⚠ CAUTION**

**Manual unlocking of elbow lock under load**

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

**⚠ CAUTION**

**Overheating of the product**

Injury due to loss of product functionality.

- ▶ Do not expose the product to direct sunlight or heat sources.
- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range. See the section "Technical Data".

**⚠ CAUTION**

**Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.99 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points,...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

**⚠ CAUTION**

**Overheating of drive unit due to continuous activity**

Injury due to touching the overheated drive unit.

- ▶ Avoid touching the housing of the drive unit in case of overheating.

**⚠ CAUTION**

**Unintentional unlocking of the terminal device**

Risk of injury due to releasing the terminal device from the forearm (e.g. while carrying objects).

- ▶ When using a quick-disconnect wrist unit, position the terminal device so that slight twisting does not release the terminal device from the forearm.

**⚠ CAUTION**

**Risk of pinching in the joint flexion area**

Injuries due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the joint.

**⚠ CAUTION**

**Charging the product with damaged power supply unit/charger/charger cable**

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

**NOTICE**

**Coating, gluing or painting the prosthesis**

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

**NOTICE**

**Use of incorrect power supply unit/battery charger**

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see scope of delivery).

## 5 Scope of delivery

- 1 pc. 12K110N=\* DynamicArm Plus
- or
- 1 pc. 12K100N=\* DynamicArm
  - 1 pc. 757L24 charger
  - 1 pc. Instructions for use (user)
  - 1 pc. 646D644 service pass

## 6 Charging the battery

The product is powered by an integrated high-quality Li-Ion battery with sufficient capacity to operate for one full day of usual everyday activities. Turning the product off during longer periods of passive use (e.g. air or rail travel, visit to a theatre or cinema, etc.) will make the battery last longer. We recommend charging the product once a day when used by the patient on a daily basis. Electronic battery management provides the patient with information about the battery charge level (battery management).

The charging unit consists of the charging plug and the power cord. The battery charger has an input voltage range of 100–240 V and may be operated in a mains frequency range of 50–60 Hz.

### INFORMATION

Please note the corresponding instructions for use of the battery charger.

### 6.1 Charging process

#### ⚠ WARNING

#### Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

1. Extend and switch off the product.
2. Lock the product with the pull cable.
3. Take off the elbow joint.
4. Connect the power cord to the battery charger.
5. Insert the charging plug into the charging receptacle (see figure).  
Do not force it in!
6. Plug the power cord into a wall socket. If the LED shows an orange light, the battery is charging.  
If the LED shows a green light but the battery is drained, the charging plug was not inserted correctly.
7. If the LED shows a green light after charging, this indicates that the battery is fully charged.
8. Disconnect the power cord of the battery charger from the wall socket.

### 6.2 Charging times

Charging time	Capacity	Operating time
4 hours	100%	approx. 18 hours

Charging time	Capacity	Operating time
1.5 hours	80%	approx. 14 hours
20 minutes	40%	approx. 4 hours

### 6.3 Display of the current charge level during the charging process

The battery charger has a battery capacity indicator LED:

LED shows orange light	Battery is charging.
LED shows green light	Charging process is finished and the battery is fully charged.

### 6.4 Battery management

Electronic battery management helps the patient maintain the function of individual prosthetic components as long as possible as the battery charge level drops.

Level 1	The product and the connected prosthetic components move with maximum speed and power.
Level 2	The flexion force becomes weaker, the product "fatigues". The functions of the connected prosthetic components remain intact.
Level 3	The patient is informed of the low battery charge level by a vibration signal and a sound sequence. The product stops functioning. The elbow lock can be released and locked manually. The functions of the connected prosthetic components remain intact.
Level 4	A sound sequence informs the patient that the battery is empty. Now all other prosthetic components also cease to function.

### 6.5 Charging plug

The charging plug is suitable for bilateral amputees. It is designed so it can be inserted into the charging receptacle with a prosthetic hand or by mouth.

The following points must be observed when charging the battery:

- We recommend charging the product once a day when used by the patient on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.
- Note the permissible temperature range for charging the battery (see page 24).

#### INFORMATION

The product can remain connected to the charger after charging is complete and the green LED indicator lights up. Damage to the battery caused by overcharging is not possible. Partial charging also does not decrease the service life of the battery (no memory effect).

## 7 Cleaning and Care

#### NOTICE

##### Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

- 1) Clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock Derma Clean) when needed.  
Make sure that no liquids get into the product and product components.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

## 8 Maintenance

Regular maintenance (service inspections) every 24 months is mandatory in the interest of your own safety, in order to maintain operating reliability, to protect the extended warranty and to maintain basic safety.

The grace period is no more than one month before or three months after maintenance is due.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components always have to be submitted to the O&P professional for maintenance and repairs:

The prosthesis, battery charger and power supply unit.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 9.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address:  
<http://www.ottobock.com/conformity>

### 9.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and

2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution: Exposure to Radio Frequency Radiation.**

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP  
3820 West Great Lakes Drive  
Salt Lake City, Utah 84120-7205 USA  
Phone + 1-801-956-2400  
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road  
L7L 5N5 Burlington, Ontario  
Canada  
Phone + 1-800-665-3327

---

**Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.**

## 10 Technical data

<b>Ambient conditions</b>	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 93% relative humidity, non-condensing
Charging temperature	+5 °C/+41 °F to +40 °C/+104 °F
Operation	+5 °C/+41 °F to +45 °C/+113 °F Max. 93% relative humidity, non-condensing

<b>General information</b>	
Reference number	12K100N=*/12K110N=*
Weight (dependent on forearm length)	approx. 1,000 g
Max. lifting force	50 N
Flexion angle	Approx. 15°–145°
Service life	5 years

<b>Battery of the product</b>	
Battery type	Li-Ion
Output voltage	approx. 3.7 V
Charging voltage	approx. 4.2 V
Capacity	1,880 mAh
Dimensions of battery cells	33.8 x 48.8 x 10.5 mm
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500
Weight	38.5 g (battery without options)
Charging time until battery is fully charged	4.0 h

<b>Power supply</b>	
Reference number	757L24
Storage (with and without packaging)	-25 °C/-13 °F to +70 °C/+158 °F 10% to 95% relative humidity
Transport (with and without packaging)	-25 °C/-13 °F to +70 °C/+158 °F 10% to 95% relative humidity
Operation	-25 °C/-13 °F to +40 °C/+104 °F Max. 95% relative humidity, non-condensing
Input voltage	90 V~ to 264 V~
Mains frequency	47 Hz to 63 Hz

## 11 Appendices

### 11.1 Symbols Used

Declaration of conformity according to the applicable European directives

Manufacturer

Compliance with the requirements according to "FCC Part 15" (USA)

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Medical device

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Serial number (YYYY WW NNN)

YYYY – year of manufacture

WW – week of manufacture

NNN – sequential number

Non-ionising radiation

Article number

## 11.2 Operating states/error signals

### 11.2.1 Warnings/error signals

Beep signal	Vibration signal	Error	Required action
1 x long	1 x	Critical error (e.g. a sensor is not operational)	Turn the product off and on again or, with the product turned on, connect the battery charger. Contact Ottobock Service if the error has not been resolved.
1 x long	5 x	Severe error (e.g. temperature of the lifting motor too high)	
2 x	1 x	Malfunction (e.g. incompatible component connected)	

## 11.2.2 Signals for operating states

Sound sequence		Event	Elbow joint	Electric wrist rotator	Hand/Greifer
	Low, low, high, low	The product has been switched on.	Function is present	Function is present	Function is present
	Very high, high, low, very low	Decreasing battery charge level.	No function	Function is present	Function is present
	6 x low	Battery charger was connected during operation.	No function	No function	No function
	Low, high	Prosthetic component is being controlled.	Function of the prosthetic component in question is not possible.		

### INFORMATION

It is possible that only individual or several components can no longer be controlled.

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- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

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The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website <http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road

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