



## Axon-Bus Prosthetic System

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	24
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	46
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	69
<b>ES</b> Instrucciones de uso (Usuario) .....	91
<b>PT</b> Manual de utilização (Usuário) .....	112
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	134
<b>SV</b> Bruksanvisning (Användare) .....	156
<b>DA</b> Brugsanvisning (Bruger) .....	176
<b>NO</b> Bruksanvisning (Bruker) .....	197
<b>FI</b> Käyttöohje (Käyttäjä) .....	218
<b>CS</b> Návod k použití (Uživatel) .....	239
<b>JA</b> 取扱説明書(ユーザー用) .....	260

**INFORMATION**

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### INFORMATION

Date of last update: 2020-09-15

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.


The Axon-Bus prosthesis solution is referred to as the product in the following.

## 2 Product description

### 2.1 Design

The Axon-Bus prosthetic system comprises the following components:

1. Elbow joint  
(12K500 AxonArm Hybrid/12K501 Axon-Arm Ergo)
2. Charging receptacle
3. Rotation unit  
(9S501 AxonRotation adapter/9S503 Axon-Rotation)
4. Gripping component  
(8E500=\* Michelangelo Hand/8E600=\* AxonHook) with connected wrist unit

### 2.2 Function

#### 2.2.1 Elbow joint

##### 12K500 AxonArm Hybrid

The 12K500 AxonArm Hybrid is a mechanically locked elbow joint.

Locking and unlocking is carried out mechanically, e.g. by means of body harnesses.

##### 12K501 AxonArm Ergo

The 12K501 AxonArm Ergo is a myoelectrically locked elbow joint.

Locking and unlocking is controlled myoelectrically by means of electrode signals or with a switch.

#### 2.2.2 Battery charger

The 757L500 AxonCharge Integral battery charger is used to charge the battery. LEDs on the charger provide information about the following:

- The battery status of the charger
- Operating states of the Axon-Bus prosthetic system
- Battery charge level

#### 2.2.3 Battery

Power is supplied to the Axon-Bus prosthesis solution by the battery integrated in the socket (757B501 AxonEnergy Integral).

The charging receptacle on the socket has the following functions:

- Contact for charging the battery
- LEDs to indicate the operating state and battery charge level
- Button to turn the Axon-Bus prosthesis solution on/off
- Beeper for audible feedback on operating states

#### 2.2.4 Rotation unit

##### 9S501=\* AxonRotation Adapter (mechanical rotation unit)

The gripping component can be manually rotated continuously through 360°. There are 24 ratchet positions for this purpose, one at every 15°.


##### 9S503=\* AxonRotation (active rotation unit)

The gripping component can be actively rotated 160° from the neutral position to pronation and supination.

## 2.2.5 Gripping component

### 8E500 Michelangelo Hand

Complex gripping kinematics combined with an anatomical appearance and low weight support everyday activities. The following gripping options are possible:


#### Neutral Position

Natural, physiological appearance in the rest position.


#### Lateral Power Grip

The thumb moves laterally to the index finger; this laterally fixates medium-sized objects when the thumb is in the half-open position.


#### Lateral Pinch

The thumb moves laterally to the index finger; this laterally fixates flat objects when the thumb is in the closed position.

#### Open Palm

In the Open Palm position the thumb is positioned far to the outside; this achieves a flat hand position with a completely opened thumb position.


#### Opposition Power Grip

The opening width makes it possible to hold cylindrical objects with large diameters and with an open thumb position.

#### Tripod Pinch

Together with the index finger and the middle finger, the thumb forms a three-point support; small objects can then be securely fixated with the thumb in closed position.

### **Finger Abduction/Adduction**

By spreading the fingers, a flat, thin object (< 3 mm / < 0.12 inch) can be held between the fingertips when closing the hand.

### **8E600 Axon Hook**

The gripping component consists of two hooks ("hook tips"), one hook being fixed and the second hook movable. Objects can be grasped by opening and closing the hook tips.

### **2.3 Combination possibilities**

This product is to be combined exclusively with components in the Axon-Bus prosthetic system. Components of the Ottobock MyoBock system or components from other manufacturers cannot be used with this product.

## **3 Intended use**

### **3.1 Indications for use**

The Axon-Bus prosthetic system is to be used **exclusively** for exoprosthetic fittings of the upper limbs.

### **3.2 Conditions of use**

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see page 39).

### **3.3 Indications**

- Amputation level transradial and transhumeral
- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm
- The user must be able to understand usage instructions and safety notices and put them into practice.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

### **3.4 Contraindications**

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### **3.5 Qualification**

The fitting of a user with the product may only be carried out by O&P professionals who have been authorised through the corresponding Ottobock training.

## **4 Safety**

### **4.1 Explanation of warning symbols**

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
---------------------------------------------------------------------------------------------------	-----------------------------------------------------------------

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

## 4.2 Structure of the safety instructions

### The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

## 4.3 General safety instructions

### Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

### Operating a vehicle with the prosthesis

Accidents due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

### Operating machines with the prosthesis

Injury due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of industrial machines or motor-driven equipment.

### Operating the prosthetic system near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the prosthetic system.

- ▶ When operating the prosthetic system in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

### Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.

- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

**⚠ WARNING**

**Charging the prosthesis without taking it off**

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

**⚠ CAUTION**

**Unauthorised manipulation of Axon-Bus components**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Manipulations to the Axon-Bus prosthetic system other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The Axon-Bus prosthetic system and any damaged Axon-Bus components may only be opened and repaired by authorised, qualified Ottobock personnel.

**⚠ CAUTION**

**Storage of the gripping component in the closed state**

Injury due to faulty control or malfunction of the gripping component because of damage to the sensors and mechanism.

- ▶ Store the gripping component in the neutral position or open state only.

**⚠ CAUTION**

**Signs of wear and tear on the product components**

Injury due to faulty control or malfunction of the product

- ▶ In the interest of preventing injuries and to maintain product quality, we recommend performing regular service.
- ▶ Please contact your O&P professional for more information.

#### 4.4 Information on Proximity to Certain Areas

**⚠ CAUTION**

**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the prosthetic system.

**⚠ CAUTION**

**Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

#### CAUTION

##### **Wearing in extreme ambient temperatures**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see page 39).

## 4.5 Information on Use

#### CAUTION

##### **Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)**

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

#### CAUTION

##### **Overloading due to unusual activities**

Injury due to unexpected behaviour of the Axon-Bus prosthetic system because of malfunction.

- ▶ The Axon-Bus prosthetic system was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the Axon-Bus prosthetic system and the Axon-Bus components not only increases their service life, but is primarily in the interest of your own safety.
- ▶ If the Axon-Bus prosthetic system and the Axon-Bus components have been subjected to extreme loads (e.g. through a fall or similar occurrence) the prosthesis system must be inspected by a prosthetist immediately for possible damage. If necessary, the prosthetist will forward the prosthesis system to an authorised Ottobock Service Centre.

#### CAUTION

##### **Mechanical stress on the Axon-Bus prosthetic system**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Do not subject the Axon-Bus components to mechanical vibrations or impacts.
- ▶ Check the Axon-Bus components for any visible damage before each use.

#### CAUTION

##### **Changing Axon-Bus gripping components when switched on**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Power down the Axon-Bus prosthetic system by pressing the button in the charging receptacle before changing any Axon-Bus components (e.g. Axon-Bus gripping component).

#### CAUTION

##### **Penetration of dirt and moisture in the Axon-Bus components**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the Axon-Bus prosthetic system or Axon-Bus components (e.g. the Axon-Bus gripping component).

**⚠ CAUTION**

**Unintentional unlocking of the Axon-Bus gripping component**

Injury due to releasing the Axon-Bus gripping component from the forearm (e.g. while carrying objects).

- ▶ The two release buttons should only be activated in order to change the Axon-Bus gripping component, deliberately and with due consideration of the respective situation.

**⚠ CAUTION**

**Risk of pinching between the hook tips**

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the hook tips when using or rotating the product.

**⚠ CAUTION**

**Gripping objects with different gripping force when using the "AxonHook" gripping component**

Injury due to unexpected behaviour of the gripping component.

- ▶ Note that the gripping force can vary depending on the consistency (soft/hard) of the object being gripped.

**⚠ CAUTION**

**Use of the Axon-Bus prosthetic system with dirty or damaged contacts of the rotation unit and gripping component**

Injury due to unexpected behaviour of the Axon-Bus prosthetic system.

- ▶ Make sure that the contacts of the rotation unit and gripping component are always clean and free of grease.
- ▶ Clean the electrical contacts regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

**NOTICE**

**Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

- ▶ Use only the following products for cleaning/disinfecting the inner socket:

**Cleaning:** 453H10=1 Ottobock DermaClean

**Disinfection:** Colourless, conventional medical disinfectants

## 4.6 Information on the Power Supply/Battery Charging

**⚠ CAUTION**

**Charging the Axon-Bus prosthetic system with soiled or damaged contacts**

Injury due to unexpected behaviour of the Axon-Bus prosthetic system caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

**NOTICE****Charging the Axon-Bus prosthetic system with an incorrect power supply/battery charger**

Damage to the Axon-Bus prosthetic system due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for the Axon-Bus prosthetic system by Ottobock (see instructions for use and catalogues).

**NOTICE****Contact of the charging plug with magnetic data carriers**

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

**INFORMATION**

The Axon-Bus prosthetic system cannot be used as long as the charging plug is connected to the charging receptacle. The Axon-Bus prosthetic system is deactivated for the duration of charging.

**4.7 Notes on Using an Elbow Joint (e.g.12K501 AxonArm Ergo, 12K500 AxonArm Hybrid)****⚠ CAUTION****Risk of pinching in elbow joint flexion region**

Injury due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the elbow joint.

**⚠ CAUTION****Manual unlocking of elbow lock under load**

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

**⚠ CAUTION****Penetration of dirt and moisture in the prosthetic arm**

Injury due to unexpected behaviour of the Axon-Bus prosthetic system or malfunction.

- ▶ Do not let foreign particles or liquids get into the prosthetic arm.
- ▶ Do not expose the prosthetic arm, especially the elbow joint, to dripping or splashing water.
- ▶ Wear the prosthetic arm and especially the elbow joint under suitably resistant clothing in the rain.

**NOTICE****Coating, gluing or painting the prosthesis**

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

## 5 Scope of Delivery and Accessories

### 5.1 Scope of delivery

The product is delivered preinstalled. The scope of delivery includes:

- 1 pc. Axon-Bus prosthetic system
- 1 pc. 757L500 AxonCharge Integral
- 1 pc. Instructions for use (user)

### 5.2 Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- 4X500 AxonCharge Mobile
- 757S500=\* AxonCharge adapter

## 6 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.
- Use the 757L500 battery charger to charge the battery.

### 6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round plug of the power supply to the receptacle on the battery charger so that the plug locks into place (see fig. 2).

**INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.**

- 3) Plug the power supply unit into the wall socket.

### 6.2 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle of the product.  
→ The correct connection of the battery charger to the product is indicated by feedback (see page 43).
- 2) The charging process starts.  
→ The Axon-Bus prosthetic system is switched off automatically.
- 3) Disconnect the product after the charging process is complete.

### 6.3 Display of the current charge level

#### Power supply and battery charger operating states

During the charging process, the current battery charge level (Status signals) is indicated by the LEDs on the battery charger lighting up in sequence (see fig. 3, items 2-5).

#### Displaying the charge level on the product

The battery charge level can be queried at any time.

- 1) With the prosthetic system switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

## 6.4 Safety Shutoff

The purpose of the Axon-Bus prosthetic system safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

## 7 Use

### 7.1 Switching the Axon-Bus prosthesis system on and off

- ▶ 1) Push and hold the charging receptacle button until a confirmation sound is heard (at least 1 second).  
→ The prosthesis and prosthetic components are switched on. Repeat this process to switch the prosthesis and prosthetic components off.

#### INFORMATION

Switching off the Axon-Bus prosthetic system during extended passive pauses (e.g. air or rail travel, visit to a theatre or cinema, etc.) will make the battery last longer before requiring a recharge. Only the entire Axon-Bus prosthetic system with all connected components can be switched off. Individual Axon-Bus components cannot be switched off separately.

### 7.2 Opening the Axon-Bus gripping component in an emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open.  
→ A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

### 7.3 Changing the Axon-Bus gripping component

#### ⚠ CAUTION

#### Changing Axon-Bus gripping components when switched on

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Power down the Axon-Bus prosthetic system by pressing the button in the charging receptacle before changing any Axon-Bus components (e.g. Axon-Bus gripping component).

#### Attaching the Axon-Bus gripping component to the socket

- 1) Slide the terminal device on to the socket (rotation anchored in the socket) until it locks in place audibly.
- 2) Verify proper attachment by pulling on the terminal device.

### Separating the Axon-Bus gripping component from the socket

- 1) To switch off the Axon-Bus prosthetic system, press the charging receptacle button and hold for more than one second.
- 2) Press both unlock buttons on the wrist unit (see fig. 4).

**INFORMATION: For safety reasons the Axon-Bus gripping component cannot be removed if only one unlock button is pressed.**

- 3) Pull off the Axon-Bus gripping component.

### 7.4 Proper handling of the Axon-Bus gripping component

#### INFORMATION

The handling procedures demonstrated in this section are independent of the gripping component used. The handling procedures are to be used for all gripping components.

#### Hand flexion, hand extension procedure

This function enables the flexion or extension of the Axon-Bus gripping component.

##### Flexible mode

- 1) Push the unlock lever (see fig. 5, see fig. 6) of the Axon-Bus gripping component (to the stop) until it engages.

→ The wrist unit can be moved without engaging at the ratchet positions.

- 2) When the unlock lever is pushed again, the wrist unit engages at the next available position.

##### Rigid mode

- 1) Push the unlock lever (see fig. 5, see fig. 6) of the Axon-Bus gripping component lightly (not to the stop) and move the wrist unit to the desired position.

- 2) When the unlock lever is released, the wrist unit engages at the next available position.

#### Pronation, supination procedure

#### CAUTION

##### Risk of pinching between the hook tips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the hook tips when using or rotating the product.

#### INFORMATION

The following steps apply only to the product "9S501 AxonRotation adapter" (mechanical rotation). When using the product "9S503 AxonRotation", pronation and supination are controlled by the associated muscle signals.

The wrist unit supports pronation and supination.

- ▶ Grasp the gearbox and turn the wrist unit to the desired position.

#### Storage of the gripping component

#### CAUTION

##### Storage of the gripping component in the closed state

Injury due to faulty control or malfunction of the gripping component because of damage to the sensors and mechanism.

- ▶ Store the gripping component in the neutral position or open state only.

### Supporting yourself with the gripping component

**⚠ CAUTION**

**Mechanical stress on the Axon-Bus prosthetic system**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Do not subject the Axon-Bus components to mechanical vibrations or impacts.
- ▶ Check the Axon-Bus components for any visible damage before each use.

The flexion/extension of the wrist unit must be used when supporting yourself on the product.

## 8 Cleaning

### NOTICE

#### Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

- 1) Clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean) when needed.  
Make sure that no liquids get into the product and product components.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

## 9 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

The grace period is no more than one month before or three months after maintenance is due.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components always have to be submitted to the O&P professional for maintenance and repairs:

The prosthesis, battery charger and power supply unit.

## 10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

### 10.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution: Exposure to Radio Frequency Radiation.**

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP  
3820 West Great Lakes Drive  
Salt Lake City, Utah 84120-7205 USA  
Phone + 1-801-956-2400  
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website <http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road

**Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.**

**11 Technical data**

**Michelangelo Hand gripping component**

<b>Ambient conditions</b>	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing

<b>General information</b>	
Reference number	8E500
Opening width	120 mm/4.72 inch
Weight of Michelangelo hand alone, without the AxonRotation Adapter and without a prosthetic glove	approx. 510 g/18 oz
Flexion of the wrist	75°/four locking positions
Extension of the wrist	45°/three locking positions
Expected lifetime given compliance with the recommended maintenance intervals	5 years

The following gripping forces and load limits apply only when the battery of the Axon-Bus prosthetic system is fully charged and at room temperature.

<b>Maximum gripping forces</b>	
Gripping force "Opposition Mode"	Approx. 70 N
Gripping force "Lateral Mode"	Approx. 60 N
Gripping force "Neutral Mode"	Approx. 15 N

<b>Load limits</b>	
Maximum vertical load on the palm with locked wrist joint (e.g. when holding a sphere)	10 kg / 22.1 lbs
Maximum load on the actively driven fingers (index finger, middle finger) in the fully open hand position (e.g. holding a plate)	10 kg / 22.1 lbs
Maximum load on the actively driven fingers (index finger, middle finger) in the closed hand position (e.g. carrying bags)	20 kg / 44.1 lbs
Maximum vertical load on the knuckles (e.g. bracing oneself on the fist)	60 kg / 132 lbs

<b>Load limits</b>	
Weight of objects (typical diameter 19 mm/0.75 inch) before they slide out of the hand ("Power Grip" grip type)	18 kg / 39.6 lbs

### **AxonHook gripping component**

<b>Ambient conditions</b>	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing

<b>General information</b>	
Reference number	8E600=*
Min. opening width	130 mm/5.12 inch
Min. closing time	0.75 sec
Max. closing time	10 sec
Weight	440 g/15.52 oz
Expected lifetime given compliance with the recommended maintenance intervals	500,000 gripping cycles/5 years

<b>Maximum gripping forces</b>	
Gripping force at the hook tips	approx. 110 N

<b>Load limits</b>	
Maximum lateral load on the hooks with locked wrist unit (e.g. holding an object/horizontal forearm)	16 kg/35.2 lb
Maximum load on the hooks in the closed hook position (e.g. carrying bags/vertical forearm)	50 kg/110.2 lb
Maximum vertical weight on the hooks (e.g. bracing)	60 kg/132.3 lb

### **AxonRotation wrist unit**

<b>General information</b>	
Reference number	9S503
Expected lifetime given compliance with the recommended maintenance intervals	5 years
Active pronation from the neutral position	160°
Active supination from the neutral position	160°
Maximum torque	1.5 Nm
Idle speed	25 rpm
Angular rotation	150°/sec.
Weight	approx. 140 g/4.94 oz

<b>General information</b>	
Reference number	9S501
Expected lifetime given compliance with the recommended maintenance intervals	5 years
Pronation	360° (24 locking positions)
Supination	360° (24 locking positions)
Weight	approx. 90 g/3.17 oz

## Elbow joint

<b>General information</b>	<b>AxonArm Hybrid</b>	<b>AxonArm Ergo</b>
Reference number	12K500=*	12K501=*
Expected lifetime given compliance with the recommended maintenance intervals	5 years	5 years
Weight	approx. 840 g/29.6 oz	approx. 840 g/29.6 oz
Flexion angle	15°–145°	15°–145°
Maximum vertical load with locked elbow joint and a forearm length of 305 mm/12 inch	23 kg/50.7 lb	

## Axon-Bus prosthetic system

<b>Prosthesis battery</b>	
Battery type	Li-Ion
Battery service life	2 years
Charging time until battery is fully charged	4 hours
Behaviour of the product during the charging process	The product is non-functional
Operating time of the product with fully charged battery	1 day with average use

## 12 Appendices

### 12.1 Symbols Used

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Serial number (YYYY WW NNN)  
 YYYY – year of manufacture  
 WW – week of manufacture  
 NNN – sequential number

**MD**

Medical device

Legal manufacturer

**IP33**

Protection against penetration of solid foreign objects with a diameter &gt; 2.5 mm, protection against water spray

**LOT**

Lot number (PPPP YYYY WW)

PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

**REF**

Article number

**12.2 Operating states/error signals**

The product indicates operating states audibly and visually.


**12.2.1 Signals for operating states****Battery charger connected/disconnected**

Beep signal	Event
1 x short	Battery charger connected
2 x short	Charging begins
1 x long	Charging ends




**Switching the Axon-Bus prosthetic system on/off**

Charging receptacle	Beeper	Event
	2 x short	Switch on the Axon-Bus prosthetic system
	1 x long	Switch off the Axon-Bus prosthetic system


**12.2.2 Warnings/error signals****Error during use**

Charging receptacle	Beeper	Explanation/correction
	1x long	<b>Error in the Axon-Bus prosthetic system</b> <ul style="list-style-type: none"> <li>Switch the Axon-Bus prosthetic system off/on</li> <li>Open/close the Axon-Bus prosthetic system</li> <li>Consult prosthetist if error persists</li> </ul>
	-	<b>Axon-Bus prosthetic system outside operating temperature range</b> <ul style="list-style-type: none"> <li>Allow the Axon-Bus prosthetic system to warm up or cool down (Technical data)</li> </ul>

Charging receptacle	Beeper	Explanation/correction
-	3 x short	<b>Axon-Bus prosthetic system switches itself off</b> <ul style="list-style-type: none"> <li>Charge battery</li> </ul>
-	Pulsating	<b>Opening the Axon-Bus gripping component in an emergency</b> <ul style="list-style-type: none"> <li>System is switched off, system can be used without restriction after being switched on again</li> </ul>



Battery charger	Explanation/correction
  	<b>An error has occurred in the Axon-Bus prosthetic system</b> <ul style="list-style-type: none"> <li>Consult prosthetist</li> </ul>
  	<b>Axon-Bus prosthetic system requires service</b> <ul style="list-style-type: none"> <li>Consult prosthetist</li> </ul>




### Error while charging the product

Battery charger	Explanation/correction
  	<b>Axon-Bus prosthetic system is not charging</b> <ul style="list-style-type: none"> <li>Disconnect and then reconnect the charging plug.</li> <li>Clean charger contacts.</li> </ul>
      	<b>Charger is overheated or overcooled</b> <ul style="list-style-type: none"> <li>Disconnect charging plug and allow charger to cool down/warm up.</li> </ul>
  	<b>Battery charger or power supply unit is defective</b> <ul style="list-style-type: none"> <li>Consult prosthetist</li> </ul>

### 12.2.3 Status signals

#### Battery charger is connected

Battery charger	Explanation
  	Power supply and battery charger operational
  	Battery is charging
  	Battery is charging, battery 25% charged

Battery charger	Explanation
  	Battery is charging, battery 50% charged
  	Battery is charging, battery 75% charged
  	Battery is fully charged, battery 100% charged

### Battery charge level

Charging receptacle	Event
	Battery is fully charged Battery charge level is 100%
	Battery charge level is over 50%
	Battery charge level is less than 50%

### 12.3 Glossary

The term "Axon" stands for **A**daptive **e**xchange **o**f **n**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

### 12.4 Directives and manufacturer's declaration

#### 12.4.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (see page 29).

#### Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

## Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air,
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	$\pm 2$ kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage drops	IEC 61000-4-11	0% $U_T$ ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% $U_T$ ; 1 period and 70% $U_T$ ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% $U_T$ ; 250/300 periods

## Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM $\pm 5$ kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, GSM 800/900, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 to 2,570	Bluetooth WLAN 802.11 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

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- Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/ TV technician for help.
- Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
 5470 Harvester Road  
 L7L 5N5 Burlington, Ontario  
 Canada  
 Phone + 1-800-665-3327

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**The 8E500 Michelangelo is covered by the following patents:**

Australia: AU 2006 332 253; AU 2006 332 292; AU 2006 332 315; AU 2006 332 317; AU 2006 332 318; AU 2006 332 316;

Canada: CA 2 631 970; CA 2 631 982; CA 2 632 241; CA 2 632 551

China: CN 101 340 864; CN 101 340 865; CN 101 340 866; CN 101 340 867; CN 101 346 107; CN 101 346 109; CN 101 681 704

Japan: JP 4 843 055; JP 5 242 409; JP 5 123 206; JP 5 074 414; JP 5 155 183; JP 5 389 667; JP 5 074 415

Mexico: MX 285 670; MX 286 932; MX 291 871; MX 291 872; MX 292 784

Russia: RU 2 387 412; RU 2 414 871; RU 2 416 379; RU 2 423 952; RU 2 427 348; RU 2 429 803; RU 2 469 429

South Korea: KR 101 131 692; KR 101 169 834; KR 101 178 679; KR 101 178 646

USA: US 7 867 287; US 8 016 969; US 8 257 446; US 8 188 835; US 8 579 991

Patents pending in Australia, Brazil, Canada, China, EPA, Germany, India, Japan, Mexico, South Korea, Russia, Taiwan and USA.

**The 8E600 AxonHook is covered by the following patents:**

Russia: RU 2 472 470

European Patent EP 2175808 in DE, FR, GB, IT, NL, SE, TR

Patents pending in Canada, Germany, Taiwan and USA.

**The 8E500 Michelangelo or parts thereof are covered by the following registered designs and design patents:**

Australia: 317214; 317789

Canada: ©Ottobock No. 122162; 122163

China: ZL 200730154423.X; ZL 200730154429.7; ZL 201130050582.1; ZL 201130050654.2

European Design: No.000786421; No.000786694; No.001824004

Germany: 40701345.8; 40701357.1

Japan: Registered Design No. 1 365 277; 1 365 278

Russia: 69 071; 70 542

Taiwan: R.O.C. Design Patent D 128 170; D 128 171

USA: Patent US D 595,854; D 597,672; D 694,189

Design Patents pending in India.

Michelangelo is a tradename of Otto Bock Healthcare GmbH.

Michelangelo is a registered trademark in many countries of the world, beside others registered at the US Patent and Trade-mark Office, Reg.-No. 4008171.

Otto Bock Healthcare Products GmbH  
Brehmstraße 16 · 1110 Wien · Austria  
T +43-1 523 37 86 · F +43-1 523 22 64  
info.austria@ottobock.com · www.ottobock.com