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

INFORMATION

Date of last update: 2020-10-06

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The 8E550=R/L Michelangelo Hand Transcarpal is referred to as the product / Axon-Bus terminal device / Michelangelo Hand below.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Construction and Function

The Axon-Bus prosthetic system comprises the following components:

2.1.1 Battery charger

The 757L500 AxonCharge Integral battery charger is used to charge the battery. LEDs on the charger provide information about the following:

- The battery status of the charger
- Operating states of the Axon-Bus prosthetic system
- Battery charge level

2.1.2 Battery

Power is supplied to the Axon-Bus prosthetic system by the battery integrated in the socket (757B501 AxonEnergy Integral).


The charging receptacle on the socket provides the following functions:

- Contact for charging the battery
- LEDs to indicate the operating state and battery charge level
- Button to switch the Axon-Bus prosthetic system on and off
- Beeper for audible feedback on operating states

2.1.3 Gripping component


8E550 Michelangelo Hand Transcarpal

Complex gripping kinematics combined with an anatomical appearance and low weight support everyday activities. The following gripping options are possible:


Neutral Position

Natural, physiological appearance in the rest position.


Lateral Power Grip

The thumb moves laterally to the index finger; this laterally fixates medium-sized objects when the thumb is in the half-open position.

Lateral Pinch


The thumb moves laterally to the index finger; this laterally fixates flat objects when the thumb is in the closed position.

Open palm

In the open palm position, the thumb is positioned all the way out; this allows for a flat hand position with a completely opened thumb position. The degree of mobility for internal and external rotation of the forearm depends on the length of the residual limb. A compensatory movement from the upper arm may be required in order to hold a flat object in a horizontal position.

INFORMATION: The starting position of the rotation setting can be adjusted by the O&P professional within a range of $\pm 15^\circ$.


Opposition Power Grip

The opening width makes it possible to hold cylindrical objects with large diameters and with an open thumb position.

Tripod Pinch

Together with the index finger and the middle finger, the thumb forms a three-point support; small objects can then be securely fixated with the thumb in closed position.

Finger Abduction/Adduction

By spreading the fingers, a flat, thin object (< 3 mm / < 0.12 inch) can be held between the fingertips when closing the hand.

3 Intended use

3.1 Indications for use

The 8E550 Michelangelo Hand Transcarpal is to be used **exclusively** for upper limb exoprosthetic fittings.

3.2 Conditions of use

The Axon-Bus prosthetic system was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the socket connection (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the Axon-Bus prosthetic system should not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The prosthetic system is intended **exclusively** for use on **one** user. Use of the prosthetic system by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see page 31).

3.3 Indications

The 8E550 Michelangelo Hand Transcarpal can be used for unilateral or bilateral amputees from a transradial to transcarpal amputation level or in the case of dysmelia of the forearm.

- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 **WARNING**

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 **WARNING**

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

 **WARNING**

Operating a vehicle with the prosthesis

Accidents due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

 **WARNING**

Operating machines with the prosthesis

Injury due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of industrial machines or motor-driven equipment.

 **WARNING**

Operating the prosthetic system near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the prosthetic system.

- ▶ When operating the prosthetic system in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ WARNING

Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

⚠ CAUTION

Unauthorised manipulation of Axon-Bus components

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Manipulations to the Axon-Bus prosthetic system other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The Axon-Bus prosthetic system and any damaged Axon-Bus components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Signs of wear and tear on the product components

Injury due to faulty control or malfunction of the product

- ▶ In the interest of preventing injuries and to maintain product quality, we recommend performing regular service.
- ▶ Please contact your O&P professional for more information.

4.4 Information on Proximity to Certain Areas

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the prosthetic system.

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

⚠ CAUTION

Wearing in extreme ambient temperatures

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see page 31).

4.5 Information on Use

⚠ WARNING

Operating the prosthesis near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ CAUTION

Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

⚠ CAUTION

Risk of pinching between the fingertips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the fingertips when using the product.

⚠ CAUTION

Overloading due to unusual activities

Injury due to unexpected behaviour of the Axon-Bus prosthetic system as a result of a malfunction.

- ▶ The Axon-Bus prosthetic system was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the socket connection (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the Axon-Bus prosthetic system and the Axon-Bus components not only increases their service life, but is primarily in the interest of your own safety.
- ▶ If the Axon-Bus prosthetic system and the Axon-Bus components have been subjected to extreme loads (e.g. through a fall or similar occurrence) the prosthesis system must be inspected by an O&P professional immediately for possible damage. If necessary, the prosthetist will forward the prosthetic system to an authorised Ottobock Service Centre.

⚠ CAUTION

Mechanical stress on the Axon-Bus prosthetic system

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Do not subject the Axon-Bus components to mechanical vibrations or impacts.
- ▶ Check the Axon-Bus components for any visible damage before each use.

⚠ CAUTION

Penetration of dirt and moisture in the Axon-Bus components

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the Axon-Bus prosthetic system or Axon-Bus components (e.g. the Axon-Bus gripping component).

NOTICE

Improper care of the Axon-Bus prosthetic system

Damage to the Axon-Bus prosthetic system due to the use of incorrect cleaning agents.

- ▶ Clean the Axon-Bus prosthetic system with a damp cloth and mild soap only (e.g. 453H10=1 Ottobock DermaClean).

- ▶ Use only the following products for cleaning/disinfecting the inner socket:

Cleaning: 453H10=1 Ottobock DermaClean

Disinfection: 453H16 Ottobock DermaDesinfect.

4.6 Information on the Power Supply/Battery Charging

⚠ WARNING

Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

⚠ CAUTION

Charging the Axon-Bus prosthetic system with soiled or damaged contacts

Injury due to unexpected behaviour of the Axon-Bus prosthetic system caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

NOTICE

Charging the Axon-Bus prosthetic system with an incorrect power supply/battery charger

Damage to the Axon-Bus prosthetic system due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for the Axon-Bus prosthetic system by Ottobock (see instructions for use and catalogues).

NOTICE

Contact of the charging plug with magnetic data carriers

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

INFORMATION

The Axon-Bus prosthetic system cannot be used as long as the charging plug is connected to the charging receptacle. The Axon-Bus prosthetic system is deactivated for the duration of charging.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

The product is delivered preinstalled. The scope of delivery includes:

- 1 pc. Axon-Bus transcarpal prosthetic system
- 1 pc. 757L500 AxonCharge Integral
- 1 pc. Instructions for use (user)
- 1 pc. prosthesis passport
- 1 pc. cosmetic case for battery charger and power supply

5.2 Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- 4X500 AxonCharge Mobile
- 757S500=* AxonCharge adapter
- 8S501=* AxonSkin Natural for men (different skin colours)
- 8S502=* AxonSkin Natural for women (different skin colours)
- 8S500=* AxonSkin Visual (sheer, translucent)
- 8S511=* AxonSkin Silicone for men (different skin colours)
- 8S512=* AxonSkin Silicone for women (different skin colours)

6 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.
- Use the 757L500 battery charger to charge the battery.

6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round plug of the power supply to the receptacle on the battery charger so that the plug locks into place (see fig. 2).

INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.

- 3) Plug the power supply unit into the wall socket.

6.2 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle of the product.
→ The correct connection of the battery charger to the product is indicated by feedback (see page 34).
- 2) The charging process starts.
→ The Axon-Bus prosthetic system is switched off automatically.
- 3) Disconnect the product after the charging process is complete.

6.3 Display of the current charge level

Power supply and battery charger operating states

During the charging process, the current battery charge level (Status signals) is indicated by the LEDs on the battery charger lighting up in sequence (see fig. 3, items 2-5).

Querying battery charge level on prosthetic system

The battery charge level can be queried at any time.

- 1) With the prosthetic system switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

6.4 Safety Shutoff

The purpose of the Axon-Bus prosthetic system safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

7 Use

7.1 Switching the Axon-Bus prosthesis system on and off

The product is switched on and off by pressing the button in the charging receptacle on the outside of the socket.

INFORMATION

Switching off the Axon-Bus prosthetic system during extended passive pauses (e.g. air or rail travel, visit to a theatre or cinema, etc.) will make the battery last longer before requiring a recharge. Only the entire Axon-Bus prosthetic system with all connected components can be switched off. Individual Axon-Bus components cannot be switched off separately.

7.2 Opening the Axon-Bus gripping component in an emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open.
→ A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

INFORMATION

If emergency opening of the Michelangelo Hand cannot be performed, the entire prosthesis must be removed from the residual limb if necessary.

7.3 Proper handling of the Axon-Bus gripping component

Storage of the gripping component

⚠ CAUTION

Storage of the gripping component in the closed state

Injury due to faulty control or malfunction of the gripping component because of damage to the sensors and mechanism.

- ▶ Store the gripping component in the neutral position or open state only.

⚠ CAUTION

Mechanical stress on the Axon-Bus prosthetic system

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Do not subject the Axon-Bus components to mechanical vibrations or impacts.
- ▶ Check the Axon-Bus components for any visible damage before each use.

8 Cleaning

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

► Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1-N DermaClean) when needed.
Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

9 Maintenance

- Inspect the entire prosthesis for wear during normal consultations.

9.1 Service

Since all moving mechanical components are subject to wear and tear, regular service is required within the guarantee period. The entire product is inspected by Ottobock Service (Ottobock Myo-Service). Wear parts are replaced as required.

Failure to adhere to the service schedule invalidates the warranty.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

10.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and

2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA
Phone + 1-801-956-2400
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.
5470 Harvester Road
L7L 5N5 Burlington, Ontario
Canada
Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

11 Technical data

| Environmental conditions | |
|---|---|
| Storage and transport in original packaging | -20 °C/-4 °F to +40 °C/+104 °F |
| Storage and transport without packaging | -20 °C/-4 °F to +40 °C/+104 °F max. 80% relative humidity, non-condensing |
| Operation | -10 °C/+14 °F to +60 °C/+140 °F max. 80% relative humidity, non-condensing |
| Charging the battery | +5 °C/+41 °F to +40 °C/+104 °F |

| General information | |
|--|------------------------|
| Reference number | 8E550=L-M or 8E550=R-M |
| Opening width | 120 mm/4.72 inch |
| Weight of Michelangelo Hand without prosthetic glove | approx. 460 g / 16 oz |
| Product service life | 5 years |

The following gripping forces and load limits apply only when the battery of the Axon-Bus prosthetic system is fully charged and at room temperature.

| Maximum gripping forces | |
|----------------------------------|--------------|
| Gripping force "Opposition Mode" | Approx. 70 N |
| Gripping force "Lateral Mode" | Approx. 60 N |
| Gripping force "Neutral Mode" | Approx. 15 N |

| Load limits | |
|--|----------------|
| Maximum vertical load on the palm (e.g. when holding a sphere) | 10 kg/22.1 lbs |
| Maximum load on the actively driven fingers (index finger, middle finger) in the fully open hand position (e.g. holding a plate) | 10 kg/22.1 lbs |
| Maximum load on the actively driven fingers (index finger, middle finger) in the closed hand position (e.g. carrying bags) | 20 kg/44.1 lbs |
| Maximum vertical load on the knuckles (e.g. bracing oneself on the fist) | 60 kg/132 lbs |
| Weight of objects (typical diameter 19 mm/0.75 inch) before they slide out of the hand ("Power Grip" grip type) | 18 kg/39.6 lbs |

| Data transfer | |
|-----------------------------|------------------------|
| Wireless technology | Bluetooth 2.0 |
| Range | approx. 10 m / 32.8 ft |
| Frequency range | 2402 MHz to 2480 MHz |
| Modulation | GFSK |
| Data rate (over the air) | 704 kbps |
| Maximum output power (EIRP) | -1.78 dBm (0.66 mW) |

| Prosthesis battery | |
|---------------------------|--------|
| Battery type | Li-Ion |

| Prosthesis battery | |
|--|-------------------------------|
| Battery service life | 2 years |
| Charging time until battery is fully charged | 4 hours |
| Behaviour of the product during the charging process | The product is non-functional |
| Operating time of the product with fully charged battery | 1 day with average use |

12 Appendices

12.1 Symbols Used

Manufacturer

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)

YYYY – year of manufacture

WW – week of manufacture

NNN – sequential number

Lot number (PPPP YYYY WW)

PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

Medical device

12.2 Operating states/error signals

The product indicates operating states audibly and visually.

12.2.1 Signals for operating states

Battery charger connected/disconnected


| Beep signal | Event |
|--------------------|---------------------------|
| 1 x short | Battery charger connected |
| 2 x short | Charging begins |
| 1 x long | Charging ends |

Switching the Axon-Bus prosthetic system on/off

| Charging receptacle | Beeper | Event |
|--|-----------|---|
|  | 2 x short | Switch on the Axon-Bus prosthetic system |
| | 1 x long | Switch off the Axon-Bus prosthetic system |


12.2.2 Warnings/error signals


Error during use

| Charging receptacle | Beeper | Explanation/correction |
|--|-----------|---|
|  | 1x long | Error in the Axon-Bus prosthetic system <ul style="list-style-type: none"> Switch the Axon-Bus prosthetic system off/on Open/close the Axon-Bus prosthetic system Consult prosthetist if error persists |
|  | - | Axon-Bus prosthetic system outside operating temperature range <ul style="list-style-type: none"> Allow the Axon-Bus prosthetic system to warm up or cool down (Technical data) |
| - | 3 x short | Axon-Bus prosthetic system switches itself off <ul style="list-style-type: none"> Charge battery |
| - | Pulsating | Opening the Axon-Bus gripping component in an emergency <ul style="list-style-type: none"> System is switched off, system can be used without restriction after being switched on again |

| Battery charger | Explanation/correction |
|--|--|
|  | An error has occurred in the Axon-Bus prosthetic system <ul style="list-style-type: none"> Consult prosthetist |
|  | Axon-Bus prosthetic system requires service <ul style="list-style-type: none"> Consult prosthetist |







Error while charging the product

| Battery charger | Explanation/correction |
|---|---|
|  | Axon-Bus prosthetic system is not charging <ul style="list-style-type: none"> Disconnect and then reconnect the charging plug. Clean charger contacts. |
|   | Charger is overheated or overcooled <ul style="list-style-type: none"> Disconnect charging plug and allow charger to cool down/warm up. |

| Battery charger | Explanation/correction |
|--|--|
|  | |
|  | Battery charger or power supply unit is defective <ul style="list-style-type: none"> Consult prosthetist |

12.2.3 Status signals

Battery charger is connected

| Battery charger | Explanation |
|--|--|
|  | Power supply and battery charger operational |
|  | Battery is charging |
|  | Battery is charging, battery 25% charged |
|  | Battery is charging, battery 50% charged |
|  | Battery is charging, battery 75% charged |
|  | Battery is fully charged, battery 100% charged |

Battery charge level

| Charging receptacle | Event |
|--|--|
|  | Battery is fully charged Battery charge level is 100% |
|  | Battery charge level is over 50% |
|  | Battery charge level is less than 50% |

12.3 Glossary

The term "Axon" stands for **Adaptive exchange of neuroplacement data**. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

12.4 Directives and manufacturer's declaration

12.4.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (see page 23).

Electromagnetic emissions

| Interference measurements | Compliance | Electromagnetic environment directive |
|---|---|--|
| HF emissions according to CISPR 11 | Group 1/class B | The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely. |
| Harmonics according to IEC 61000-3-2 | Not applicable – power below 75 W | – |
| Voltage fluctuations/flicker according to IEC 61000-3-3 | Product meets the requirements of the standard. | – |

Electromagnetic interference immunity

| Phenomenon | EMC basic standard or test procedure | Interference immunity test level |
|---|--------------------------------------|--|
| Electrostatic discharge | IEC 61000-4-2 | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air, |
| High-frequency electromagnetic fields | IEC 61000-4-3 | 10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz |
| Magnetic fields with rated power frequencies | IEC 61000-4-8 | 30 A/m 50 Hz or 60 Hz |
| Electrical fast transients/bursts | IEC 61000-4-4 | ± 2 kV 100 kHz repetition rate |
| Surges Line against line | IEC 61000-4-5 | ± 0.5 kV, ± 1 kV |
| Conducted interference induced by high-frequency fields | IEC 61000-4-6 | 3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz |
| Voltage drops | IEC 61000-4-11 | 0% U _T ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees |
| | | 0% U _T ; 1 period and 70% U _T ; 25/30 periods Single phase: at 0 degrees |
| Voltage interruptions | IEC 61000-4-11 | 0% U _T ; 250/300 periods |

Interference resistance against wireless communication devices

| Test frequency [MHz] | Frequency band [MHz] | Radio service | Modulation | Maximum power [W] | Distance [m] | Interference immunity test level [V/m] |
|----------------------|----------------------|--|--|-------------------|--------------|--|
| 385 | 380 to 390 | TETRA 400 | Pulse modulation 18 Hz | 1.8 | 0.3 | 27 |
| 450 | 430 to 470 | GMRS 460, FRS 460 | FM ± 5 kHz deviation 1 kHz sine | 1.8 | 0.3 | 28 |
| 710 | 704 to 787 | LTE band 13, 17 | Pulse modulation 217 Hz | 0.2 | 0.3 | 9 |
| 745 | | | | | | |
| 780 | | | | | | |
| 810 | 800 to 960 | GSM 800/90- 0, TETRA 800, iDEN 820, CDMA 850, GSM 800/90- 0, LTE band 5 | Pulse modulation 18 Hz | 2 | 0.3 | 28 |
| 870 | | | | | | |
| 930 | | | | | | |
| 1,720 | 1,700 to 1,990 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS | Pulse modulation 217 Hz | 2 | 0.3 | 28 |
| 1,845 | | | | | | |
| 1,970 | | | | | | |
| 2,450 | 2,400 to 2,570 | Bluetooth WLAN 802.1- 1 b/g/n, RFID 2450 LTE band 7 | Pulse modulation 217 Hz | 2 | 0.3 | 28 |
| 5,240 | 5,100 to 5,800 | WLAN 802.1- 1 a/n | Pulse modulation 217 Hz | 0.2 | 0.3 | 9 |
| 5,500 | | | | | | |
| 5,785 | | | | | | |

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The 8E550 Michelangelo is covered by the following patents:

Australia: AU 2006 332 253; AU 2006 332 292; AU 2006 332 315; AU 2006 332 317; AU 2006 332 318; AU 2006 332 316;

Canada: CA 2 631 970; CA 2 631 982; CA 2 632 241; CA 2 632 551

China: CN 101 340 864; CN 101 340 865; CN 101 340 866; CN 101 340 867; CN 101 346 107; CN 101 346 109; CN 101 681 704

Japan: JP 4 843 055; JP 5 242 409; JP 5 123 206; JP 5 074 414; JP 5 155 183; JP 5 389 667; JP 5 074 415

Mexico: MX 285 670; MX 286 932; MX 291 871; MX 291 872; MX 292 784

Russia: RU 2 387 412; RU 2 414 871; RU 2 416 379; RU 2 423 952; RU 2 427 348; RU 2 429 803; RU 2 469 429

South Korea: KR 101 131 692; KR 101 169 834; KR 101 178 679; KR 101 178 646

USA: US 7 867 287; US 8 016 969; US 8 257 446; US 8 188 835; US 8 579 991

Patents pending in Australia, Brazil, Canada, China, EPA, Germany, India, Japan, Mexico, South Korea, Russia, Taiwan and USA.

The 8E550 Michelangelo or parts thereof are covered by the following registered designs and design patents:

Australia: 317214; 317789

Canada: ©Ottobock No. 122162; 122163

China: ZL 200730154423.X; ZL 200730154429.7; ZL 201130050582.1; ZL 201130050654.2

European Design: No.000786421; No.000786694; No.001824004

Germany: 40701345.8; 40701357.1

Japan: Registered Design No. 1 365 277; 1 365 278

Russia: 69 071; 70 542

Taiwan: R.O.C. Design Patent D 128 170; D 128 171

USA: Patent US D 595,854; D 597,672; D 694,189

Design Patents pending in India.

Michelangelo is a tradename of Otto Bock Healthcare GmbH.

Michelangelo is a registered trademark in many countries of the world, beside others registered at the US Patent and Trade-mark Office, Reg.-No. 4008171.

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com