
# ottobock.

## 17B205=\*

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for use .....	5
<b>FR</b>	Instructions d'utilisation .....	7
<b>IT</b>	Istruzioni per l'uso .....	9
<b>ES</b>	Instrucciones de uso .....	11
<b>PT</b>	Manual de utilização .....	13
<b>NL</b>	Gebruiksaanwijzing .....	15
<b>SV</b>	Bruksanvisning .....	17
<b>DA</b>	Brugsanvisning .....	18
<b>PL</b>	Instrukcja użytkowania .....	20
<b>HU</b>	Használati utasítás .....	22
<b>CS</b>	Návod k použití .....	24
<b>HR</b>	Upute za uporabu .....	25
<b>TR</b>	Kullanma talimatı .....	27
<b>EL</b>	Οδηγίες χρήσης .....	29
<b>RU</b>	Руководство по применению .....	31
<b>JA</b>	取扱説明書 .....	33
<b>ZH</b>	使用说明书 .....	35

**INFORMATION**

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**INFORMATION**

Last update: 2016-02-09

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

## 1 Indications for use

With the E-MAG Control, the 17B205=\* Medial Support generally has to be used starting from a patient weight of 55 kg. With the E-MAG Active, the medial support does not have to be used up to a body weight of 85 kg (subject to the restrictions on use in the 647G639 Installation and Service Manual). Using the medial support is mandatory for a patient weight over 85 kg.

## 2 Delivery condition

The Medial Support is delivered with a dummy for assembly into the orthosis.

## 3 Assembly

The joint is integrated into the to-be-fabricated orthosis with a 605P8=\* Aluminum System Bar.

The dummy is used for assembly into the orthosis.

For trial fittings, the Medial Support can be fixed to the bars with screws.

Before the final delivery of the orthosis, the Medial Support must be glued with the 636W28 Adhesive.

## 4 Service

The Medial Support should be serviced every 6 months along with the EMAG System. The maintenance plan included with the Assembly and Service Instructions of the E-MAG Joint System can be used as a guideline.

The Medial Support should be serviced every 6 months along with the EMAG System. The maintenance plan included with the Assembly and Service Instructions of the E-MAG Joint System can be used as a guideline.

For maintenance works, use the 17BS205 Maintenance Set.

The bushing with collar (2) and 4 x bearing washer (1) in different colours are enclosed to the maintenance set. When replacing the bearing washers, please make sure to use always washers of the same colour.

Washer thickness	Colour
0.4	blue
0.5	white
0.6	black
0.7	red

Tighten the joint screw to a torque of 35 Nm.

## 5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH