

## 17PA1=\*

<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use .....	10
<b>FR</b>	Instructions d'utilisation .....	15
<b>IT</b>	Istruzioni per l'uso .....	20
<b>ES</b>	Instrucciones de uso .....	25
<b>PT</b>	Manual de utilização .....	31
<b>NL</b>	Gebruiksaanwijzing .....	37
<b>SV</b>	Bruksanvisning .....	43
<b>DA</b>	Brugsanvisning .....	47
<b>NO</b>	Bruksanvisning .....	53
<b>PL</b>	Instrukcja użytkowania .....	58
<b>HU</b>	Használati utasítás .....	64
<b>CS</b>	Návod k použití .....	70
<b>HR</b>	Upute za uporabu .....	75
<b>TR</b>	Kullanma talimatı .....	80
<b>RU</b>	Руководство по применению .....	86


1

2

3

## INFORMATION

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### INFORMATION

Last update: 2017-03-29

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide you with important information on the processing of the 17PA1=\* system ankle joints made of plastic material.

## 2 Product description

### 2.1 Components/design

## 3 Intended use

### 3.1 Indications for use

The system ankle joints are intended to be used **exclusively** in pairs for the fitting of lower limb orthoses for patients.

Reference number	Max. body weight
17PA1=14	45 kg
17PA1=20	100 kg

### 3.2 Indications

Partial or total paralysis of the leg muscles.

Orthopaedic conditions

Indications must be determined by the physician.

### 3.3 Further Usage Restrictions

The product is not suitable for use in sports.

## 4 Safety

### 4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

## 4.2 Safety instructions

### CAUTION

#### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).

### CAUTION

#### **Incorrect alignment or assembly**

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

### CAUTION

#### **Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

### CAUTION

#### **Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

### CAUTION

#### **Unallowable use of lubricants**

Risk of injury and loss of functionality as well as damage to the product

- ▶ The product is designed to be free of lubricants. Avoid contact with lubricants (e.g. grease, oil).

- ▶ After contact with lubricants, promptly clean the product or parts of the product with a degreasing cleaning agent (e.g. acetone or isopropyl alcohol).

#### NOTICE

#### **Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)**

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

#### NOTICE

#### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Preparation for use

#### INFORMATION

- ▶ Check the exact parallel alignment of the system joints before delivering the orthosis.

#### INFORMATION

In order to facilitate the parallel alignment of the system joints, each pair of system joints is supplied with an adjustment aid. The adjustment aid consists of 2 screws, 1 parallel truing tube and, depending on the size of the system joints, 2 adapter pieces.

The system ankle joints can be equipped with a choice of stop screws or compression springs in order to limit movement, e. g. for plantar flexion. Due to varying alignment parameters, such as the generated forefoot leverage or the flexion position of the orthosis, the system ankle joints are subjected to increased strain.

## Assembling the system joints

### INFORMATION

Stops and compression springs can be used for the individual adjustment of the system joint. The system joint can therefore be adapted to the patient's requirements.

- 1) Slide the bearing bushings (5) onto the system ankle joint (6).
- 2) Secure the foot stirrup (7) with the joint bolt (4), the joint nut with internal thread (8) and the cap screw (3).
- 3) **INFORMATION: Different set screws are used depending on the size of the system joint.** Compression spring: Install the ball (12), cylinder pin (13) and compression spring (14) with the set screw (1).
- 4) Stop: Insert the cylinder pin (2) into the system ankle joint and screw in the set screw (1).

### Installing the system bar in the system ankle joint

- 1) Insert the system bar (18) into the system bar receiver (11) and mark the bore on the system bar (18).
- 2) Remove the system bar (18) from the system bar receiver (11) and drill a  $\varnothing$  4 mm bore in the system bar (18).
- 3) Secure the system bar (18) and the system ankle joint with the set screw (9).

### Gluing the system bar

For stability reasons, the system bar (18) must be glued into the system bar receiver (11) of the system ankle joints using 636W28 Ottobock special adhesive and hardener.

- 1) Use a degreasing cleaner (such as 634A1) to degrease the bonding surfaces of the system bar (18) and the system bar receiver (11).
- 2) Mix the special adhesive and hardener thoroughly and apply it to the components being glued.
- 3) Insert the system bar (18) and secure it with the set screw (9). The final bonding strength of the adhesive is reached after **16 h**.
- 4) Prior to delivering the orthosis:  
Secure the joint bolt (4) and the joint nut (8) with 636K13 Loctite 241.  
Verify the proper functionality of the product.

## **Removing the system bar**

- 1) Unscrew the set screw (9).
- 2) Tap out and remove the system bar (18).

## **6 Cleaning**

After contact with water containing salt, chlorine or soap, or when they get dirty, the system joints must be cleaned promptly:

- 1) Rinse the system joints with clean fresh water.
- 2) Dry the system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

The system joints must be cleaned promptly after coming into contact with lubricants (e.g. grease, oil):

- 1) Moisten a cloth with a degreasing cleaning agent (e.g. acetone or isopropyl alcohol).
- 2) Wipe down the system joints or components with the cloth.

## **7 Maintenance**

The manufacturer recommends a semi-annual inspection to verify functionality and check for wear of the system joints.

Spare parts are listed in the section "Components/design".

## **8 Disposal**

Dispose of the product in accordance with national regulations.

## **9 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **9.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### **9.2 CE Conformity**

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH