## 4F34

<b>DE</b>	Gebrauchsanweisung .....	2
<b>EN</b>	Instructions for Use.....	3
<b>FR</b>	Instructions d'utilisation .....	5
<b>IT</b>	Istruzioni per l'uso .....	6
<b>ES</b>	Instrucciones de uso .....	8
<b>PT</b>	Manual de utilização.....	9
<b>SV</b>	Bruksanvisning.....	11
<b>NL</b>	Gebruiksaanwijzing.....	13
<b>DA</b>	Brugsanvisning.....	14
<b>PL</b>	Instrukcja użytkowania.....	16
<b>HU</b>	Használati utasítás .....	17
<b>CS</b>	Návod k použití.....	19
<b>TR</b>	Kullanma talimatı.....	20
<b>EL</b>	Οδηγίες χρήσης .....	22

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English

Date of the last update: 2013-07-25


- Please read this document carefully.
- Follow the safety instructions.

The 4F34 Lock Control Unit is to be used **exclusively** in lower limb components.

## 2 Description and function

A safe, easily operated manual knee lock control is very important for the geriatric patient.

The 4F34 Lock Control Unit may be used with cable controlled knee lock mechanisms that utilize a cable travel of no more than 11 mm maximum.

**Lock Control Unit without cover,  
left side, lock lever anterior, pressure activated**

In the above illustration the knee lock is released when the lock lever (2) is pushed distally. Re-routing the cable (5) through the other entrance hole (6) will cause release of the lock when the lever is pulled proximally.

To change the position of the lock lever (2) from anterior to posterior, first pull off the cover (not shown). Remove the oval head screw (3) and pull the lock lever out of the housing (1). Rotate the lock lever 180° and reassemble in the reverse order.

The cable is secured by two set screws in the bushing (4) of the lock lever. A hexagon wrench for use in securing and adjusting cable length is enclosed for your convenience.

### **3 Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

### **4 CE conformity**

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH