

749F16

DE Gebrauchsanleitung (Benutzer / Techniker)	3
EN Instructions for use (user / technician)	6
FR Mode d'emploi (Utilisateur / Technicien)	10
IT Istruzioni d'uso (utilizzatore / tecnici)	14
ES Instrucciones de uso (usuario / técnico)	18

INFORMATION

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INFORMATION

Date of the last update: 2017-02-21

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

1.1 Construction and Function

749F16 Silicone Sanding Cones are designed as tools for polishing and rounding off the edges of medical devices (e. g. prosthetic sockets) made

with thermoplastics, for use on a socket router machine or on sanding or milling machines.

The tool is suitable for processing the following materials:

Material		Processing	
Chemical composition	Otto Bock ThermoLyn® (article number)	Machine speed (rpm)	Contact pressure
PE-LD	Trolen (616T3)	1500	Slight
PE-HD 200	PE 200 (616T19, 616T58, 616T60, 616T61, 616T95, 616T495)	1500	Slight
PE-HD 500	RCH 500 (616T22, 616T43, 616T44)	1500	Slight
PE-HD 1000	RCH 1000 (616T16)	1500	Slight
PP-C	PP-C (616T20, 616T229, 616T221, 616T222)	1500	Slight
PP-H	PP-H (616T20, 616T420)	1500	Slight
Polystyrene	Rigid (616T52, 616T252)	3000	Normal
Copolyester	Clear (616T83)	3000	Normal
	PETG Clear (616T183)	3000	Normal
Ionomer	Flexible (616T39, 5Z3)	3000	Normal
EVA	Soft (616T53, 616T69, 616T253, 616T269)	3000	Normal
	Supra Soft (616T59)	3000	Normal
EVA/LDPE	EVA/LDPE Sylvershield® (616T200)	3000	Normal

2 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figure 1).

The following components are included in the scope of delivery in the specified quantities and / or are available for reordering as single components (■), single components with minimum order quantities (▲), single-component packs (●) or accessories:

749F16=5/8 Silicone Sanding Cone Set – thread 5/8"					
Fig.		Item no.	Quantity	Designation	Reference number
–	■	–	1	Instructions for Use	647G778
1	■	1	1	Silicone sanding cone, conical	749F16=5/8-1
1	■	2	1	Silicone sanding cone, conical and pointed	749F16=5/8-2

749F16=5/8 Silicone Sanding Cone Set – thread 5/8"					
Fig.		Item no.	Quantity	Designation	Reference number
1	■	3	1	Silicone sanding cone, round	749F16=5/8-3
1	■	4	1	Silicone sanding cone, round with grooves	749F16=5/8-4

749F16=M16 Silicone Sanding Cone Set – thread M16					
Fig.		Item no.	Quantity	Designation	Reference number
–	■	–	1	Instructions for Use	647G778
1	■	1	1	Silicone sanding cone, conical	749F16=M16-1
1	■	2	1	Silicone sanding cone, conical and pointed	749F16=M16-2
1	■	3	1	Silicone sanding cone, round	749F16=M16-3
1	■	4	1	Silicone sanding cone, round with grooves	749F16=M16-4

749F16=1/2 Silicone Sanding Cone Set – thread 1/2"-13					
Fig.		Item no.	Quantity	Designation	Reference number
–	■	–	1	Instructions for Use	647G778
1	■	1	1	Silicone sanding cone, conical	749F16=1/2-1
1	■	2	1	Silicone sanding cone, conical and pointed	749F16=1/2-2
1	■	3	1	Silicone sanding cone, round	749F16=1/2-3
1	■	4	1	Silicone sanding cone, round with grooves	749F16=1/2-4

3 Explanation of warning symbols

 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.

4 Use

4.1 Choosing a Tool

CAUTION

Using improper tools

Injuries, damage to the tool and/or the material to be processed

- ▶ Use only a tool that is suitable for the material to be processed (see section "Construction and Function").
- ▶ Replace any severely damaged tools (see Fig. 2, item 3) with new tools (see Fig. 2, item 1).

4.2 Mounting/Dismounting

CAUTION

Improper mounting/dismounting of the tool

Injuries

- ▶ Only mount the tool on a machine that has the corresponding thread on the tool holding fixture.
- ▶ Mount/dismount the tool only when the machine is switched off.

4.3 Handling

CAUTION

Careless working

Injuries and burns, damage to the tool and/or the material to be processed

- ▶ Observe the Operating Instructions of the machine.
 - ▶ Avoid any contact with the tool as it heats up during processing.
- 1) Use the values for speed and contact pressure given in the table in the section "Construction and Function".
 - 2) Switch on the machine and set the required speed (see Operating Instructions of the machine).
 - 3) Firmly hold the material to be processed in your hand and guide it against the tool with the required contact pressure.

4.4 Cleaning and Care

Inspection prior to each use

INFORMATION

Process-related air bubbles may appear in the silicone at the edge towards the tool holding fixture, however these do not impair the tool's functionality.

- 1) Examine the silicone surface of the tool for any damage.
- 2) **Normal signs of wear** (rough surface, see Fig. 2, item 2): Clean the tool if necessary.
- 3) **Severe signs of wear** (cracks or holes in the surface, see Fig. 2, item 3): Dispose of the tool and replace it with a new one.

Cleaning as necessary

- 1) Clean the tool with clear fresh water.
- 2) Dry the tool using a soft cloth to swab away the water.

5 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

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Otto Bock HealthCare GmbH