
ottobock.

17B98=L16, =R16, =L20, =R20

| | | |
|----|----------------------------------|----|
| DE | Gebrauchsanweisung | 2 |
| EN | Instructions for use | 4 |
| FR | Instructions d'utilisation | 6 |
| IT | Istruzioni per l'uso | 8 |
| ES | Instrucciones de uso | 10 |
| PT | Manual de utilização | 12 |
| NL | Gebruiksaanwijzing | 14 |
| SV | Bruksanvisning | 16 |
| PL | Instrukcja użytkowania | 17 |
| CS | Návod k použití | 19 |
| EL | Οδηγίες χρήσης | 21 |

INFORMATION

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INFORMATION

Last update: 2015-01-06

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1 Indications for use

The 17B98 System Ankle Joints with Lamination Adapters **are to be used exclusively** for lower limb orthotic fittings using specified lamination techniques.

2 Construction / Components

The 17B98 System Ankle Joints with Lamination Adapters can be used together with the 17B97 System Knee Joints with Lamination Adapters (please refer to the 647H160 Package Insert).

The foot stirrups have a modified joint head and the ankle joints include hardened bearing nuts. As with the System lamination joint bars, the joint's widened surface provides improved anchoring to the fiber composite material. The joint may be attached to reinforcement material by stitching through the bore holes.

The two 17Y128=16x80 System Lamination Joints (also available for 20 mm joint width and 100 mm length) provide the proximal connection of dorsiflexion assist ankle joints using the provided sidebars.

2.1 Components (shown on cover page)

- (1) 17B63=L16 System Dorsiflexion Assist Ankle Joints, 16 mm side bar width or
17B63=R16 System Dorsiflexion Assist Ankle Joints, 16 mm side bar width or
17B63=L20 System Dorsiflexion Assist Ankle Joints, 20 mm side bar width or
17B63=R20 System Dorsiflexion Assist Ankle Joints, 20 mm side bar width
- (2) 17B99=16 System Foot Stirrups (Lamination Pieces) or
17B99=20 System Foot Stirrups (Lamination Pieces)
- (3) 17Y128=16x80 System Lamination Joints, 3x80 mm, 2 pieces or
17Y128=20x80 System Lamination Joints, 3x80 mm, 2 pieces

504H1=7 Speed Rivets, 2 pieces (not illustrated)
636W28 Ottobock Special Adhesive (not illustrated)

3 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

3.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

3.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH