

Quality for life

743L30=230/110 LaserLine

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Date of the last update: 2014-05-08

- Please read this document carefully.
- Follow the safety instructions.

1 Introduction

The LaserLine must be properly handled in order to ensure successful operation and to represent no risk to individuals.

Warranty claims become void, if the apparatus has not been used as indicated, – and in the case of damages, liability will have to be assumed by the user.

2 Usage

The 743L30 LaserLine is an apparatus for checking the alignment of prostheses and orthoses.

3 Description of LaserLine

LaserLine consists of a triangle shaped basic unit with four built-in accumulators and a set-up accommodating two diode lasers which permit accurate alignment of laser light projection within the room (figs. 1 and 2).

It is operated by rechargeable accumulators and both lasers can be simultaneously activated for about 4 hours without being connected to the mains. If only one laser is in operation, use of the apparatus increases to 8 hours. Along with the LaserLine a charger is supplied, and the accumulators are charged by this via the charging socket (1).

The lasers should not be activated during the charging procedure because the diodes might be damaged due to current fluctuations.

The accumulators are charged when supplied.

The diode lasers may be turned on separately by two flip switches (2 and 3). In figure 1 the laser light outlets are marked by A and B.

- switch “I” (2) (sign for “vertical”) for the lower vertical laser, outlet A
- switch “>” (3) (sign for “angle”) for the upper, rotatable laser, outlet B.

The housing with the laser diodes can be rotated in two directions in relation to the basis. The horizontal alignment of the housing is provided by means of the adjustment screw (4) and the lock button (5) as well as the longitudinal water level (6) and the cross water level (7).

By means of the rotation button (8) the upper beam can be turned to the left by 180° and to the right by 45° about a horizontal axis.

The lower beam can be directed within the horizontal plane by means of antagonistic spring control (9) for about $\pm 9^\circ$.

In order to apply LaserLine properly, Ottobock recommends to use the apparatus only together with the stand. The stand was especially developed for this purpose (article no. 743X30) and can be ordered as accessory.

4 Transport and Installation

LaserLine is secured in a transport box when being delivered. Installation is to be provided as follows:

Remove LaserLine from the transport box: When LaserLine is to be used together with the Ottobock Stand (art. no. 743X30), attach it to the holding device of the stand by tightening the ratchet screw located at the underside of the apparatus. Make sure, that the guide pins at the underside are placed in the long-hole as well. These pins ensure that when the apparatus requires to be moved to the right or left, the direction of projection is maintained. Move the apparatus to a medial position of the long-hole.

Place the LaserLine in front of the object to be measured ensuring that the direction of projection is vertically arranged in relation to the object plane. Level the apparatus using the two operating elements – lock button (5) and adjustment screw (4), so that the air bubbles are centered in the water levels each.



Before turning the lasers on, please observe the following safety instructions.

5 Safety Instructions

Operation of laser apparatuses requires special precautionary measures.

Laser radiation can be dangerous to individuals, especially to their eyes. The following instructions have been based on the regulations for Laser Class 2. LaserLine has been classified according to Laser Class 2. In the case of apparatuses of class 2, the eyes are protected by the optico-facial winking reflex when for a short time they accidentally come to look into the laser radiation.

This laser apparatus of class 2 may be used without further protective measures, if it is ensured that the eyes look neither intentionally into the laser radiation for a time exceeding 0,25s nor do they look repeatedly into the laser radiation and/or directly reflecting laser radiation, respectively. If the optico-facial winking reflex is suppressed or delayed (for example due to medical treatment or illness), the laser apparatus should not be used.

The strong fanning of the laser beam as well as the reduction of light intensity with increasing working distance ensure, that LaserLine, when properly used, can be considered in general to be completely harmless.

According to the Accident Prevention Regulations VBG 93 "Laser Radiation" the maximally admissible radiation has been determined. From the following minimum distances radiation to the eyes is harmless for respective radiation durations:

- in the case of radiation time of 100 sec. from 15 cm
- in the case of radiation time of 1000 sec. from 45 cm

The above indications apply to the radiation with one laser beam.

Both laser beams must form a congruent line only for adjustment, because power intensity doubles.

5.1 Safety Recommendations

Make sure, that the laser beams do not meet the eyes of the patients or other persons being present.

- Avoid unintentional switching-on of the apparatus.
- The laser light comes out of the openings A and B (fig. 1).

- Make sure that the beams – if possible – will not reach the level of the patient's eyes. You may achieve this for example by coming closer to the patient with the apparatus levelled.
- If possible, project to the patient from dorsal or lateral aspects.
- Keep a minimum distance of about 1 m – if possible. See also above: minimum distances at corresponding radiation times.
- Make sure, that the beams do not meet reflecting objects (mirrors, bright surfaces).
- Do not place optical equipment (lenses, concave mirrors) into the radiation channel. Such equipment may change the beam (focussing) in a way, that it will become more dangerous, and thus the apparatus will have to be assigned to a higher laser class.
- Inform the patient and other attending persons about the danger of looking into the beam.

6 Prerequisites for Correct Application

Please make sure that the projection of light is “vertically” arranged on the object. This can be best ensured in the following way:

- The LaserLine must be horizontally aligned by means of a longitudinal water level (6) and a cross water level (7), i. e. both water bubbles should be in the center of the corresponding water level.
- The surface of the object to which the light hits, should be vertically positioned in the room, and the LaserLine should be “vertically” placed in front of it, i.e. the projection plane and the main direction of projection must be vertically positioned in relation to each other.

If these two requirements are met, the following may be stated:

- The lower laser projects a plumb line to the plane of the object.
- If the angle scale of the upper laser is in the “zero” position, this laser also projects a plumb line to the object. The two lasers then project parallel lines onto the object. The distance between the parallel lines can be changed by the antagonistic spring control (9).
- If the upper beam is turned in relation to the lower one, the angle between the beams is displayed at the angle scale on the rotation button (8).

7 Range of Application / Proper Use

In the following the different possibilities of application are described in more detail:

7.1 Initiating Steps

For most applications first the following steps are to be executed:

- Make sure that the patient takes a normal posture, or, in the same way, place the object (prosthesis, orthosis, plaster model, etc.) perpendicularly within the room.
- LaserLine should be accurately situated laterally or dorsally (depending on the kind of application) in front of the patient or the object – at a distance of approx. 1 m.
- Turn the stand until the direction of laser projection will vertically meet the patient or the object.
- Adjust LaserLine horizontally and turn the upper laser to the “zero” position using the rotation button (8).
- After making sure that the patient will not look into the beams, switch on both lasers.
- Using the antagonistic spring control (9), move the lower laser until both beams are congruently arranged.

The following adjustments of the apparatus are for different types of application.

7.2 Applications – Plumb Lines

Examples:

- Projection of the plumb line during negative plaster taking for lower extremity prostheses, orthoses and spinal orthoses (figs. 3 and 4).
- Projection of the plumb line when trying on the plaster socket and prosthetic socket or during balance (fig. 8).
- Projection of the plumb line during alignment and trial fitting of lower extremity prostheses (fig. 6).

Method:

Please, observe the initiating steps mentioned before.

After loosening the ratchet screw at the underside of the apparatus, move LaserLine to the side until the vertical laser beam will run through the requested point of the object.

7.3 Applications – Shifting

Examples:

- Demonstration of the medio-lateral shifting of the foot in relation to the knee center – from the dorsal aspect.
- Checking the deviation of the trunk (lateral deviation from C7 in relation to the anal fold) – from the dorsal aspect (fig. 7).

Method:

For these applications the plumb line – as described above – is to be projected from the dorsal aspect through the knee joint and C7, respectively.

7.4 Applications – Control of Angles

Examples:

- Flexion, abduction and adduction angles of body sections in relation to the plumb line – from the frontal or lateral aspect. (fig. 8)
- Trunk positions, for example for spinal orthoses.

Method:

Beginning from two vertical, congruent laser lines – an intersection point is provided by turning the upper laser beam coming from the “zero” position. This intersection point is the vertex of the angle to be checked.

Take this point to the requested place by vertical and lateral moving of LaserLine on the stand, for example to the hip joint or knee joint rotation point.

Turn now the upper laser to the requested direction. The angle in relation to the plumb line can now be checked by means of the angle scale of the rotation button (8).

Further examples:

- Control of the knee angle in the sagittal plane
- Control of abduction/adduction angles, in the frontal plane, of upper or lower legs relative to knee axis

Method:

In the above cases not the angle in relation to the plumb line is to be checked but that between two directions. To achieve this, the scale ring at the rotation button (8) must first be disengaged by loosening the scale ring lock screw.

Proceed as described above. After having rotated the upper laser to the direction requested first, rotate the scale in a way, that 0° is according to this direction.

Then proceed rotating the laser up to the second direction requested. Check the angle on the scale.

Note:

For correct control the scale must be placed at zero with the laser line being vertically arranged. (In case of a levelled apparatus, the laser line of the upper laser is a plumb line, only if it is parallelly placed to the line of the lower laser).

7.5 Applications with Horizontal Lines

Examples:

- Water level function for control of the pelvic position (fig. 3).
- Water level function for knee joint axes (fig. 5).
- Water level function for control of the shoulder position (fig. 7)

Method:

By turning the upper laser to the horizontal position (90°) a level line is obtained. Now move LaserLine along the stand until the horizontal line has reached the requested height.

LaserLine provides alignment analysis for orthopedic appliances. Heights, displacements and angles serve the illustration and documentation of abnormal positions and therapy results.

8 Maintenance and Repair

The apparatus requires no special maintenance work. If the LaserLine functions are not according to your expectations, it may be possible to correct these by one of the following measures. Otherwise, please contact your Ottobock Service Workshop.

8.1 Unfocused Projection of Lines

In case the apparatus does not project sharp lines, clean the lenses carefully using a cotton bud soaked with pure alcohol. If the defect cannot be eliminated in this way, the apparatus should be returned to the supplier.

8.2 Lasers do not Light

Charge the apparatus using the enclosed charger. Remove the charger plug from the charging socket (1). The lasers cannot be operated as long as the charger is connected.

8.3 Replacement of Defective Accumulators

Defective accumulators can be replaced after 4 hexagon screws of the bottom cover have been loosened. Make sure that during exchange the cables are not damaged, that polarity of the accumulators are as indicated on the accumulator pack, and that the plug connection at the accumulator pack is correctly reinstalled. Use only 625B4 Ottobock Accumulator Pack. **Do not use non-rechargeable batteries.**

9 Technical Data

Length×width×height	mm	100×110×115
Power requirements	V/Hz	230/50-60 or 110/60
Weight net/gross	kg	1.4/2.0
Color		blue – anodized

– Red light of light emitting diodes (635 nm), projected laser lines are visible during day light as well, laser class 2, output power 1 mW

- Opening angle of the laser beams: 100°
- Upper laser can be turned by 225° (+45°, -180°), with angle scale for control of the rotation angle, - 2° division.
- Lower laser with vertically projected line, can be horizontally turned by approx. ± 9°.
- Power supply: 4.8 V Accumulator Pack
- Connection thread W 1/4"
- Charger voltage 230 V/50-60 Hz or 110 V/60 Hz

10 Standard Equipment

The standard equipment consists of:

- LaserLine (1 pc.)
- Charger for 230 V/50-60 Hz or 110 V/60Hz (1 pc.)
- Transport box (1 pc.)

11 Accessories and Service Parts

625B4	Accumulator Pack
743X30	Stand

12 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

13 Declaration of Conformity

Ottobock as manufacturer with sole responsibility declares that the 743L30 LaserLine conforms to the requirements of Guideline 2004/108/EG + 2006/95/EG.

The design and construction of the apparatus were carried out in accordance with the following Standards: EN 60825-1, EN 50081-1, EN 50082-2; EN 55011; EN 60335.

The conformity of the LaserLine is guaranteed.

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Otto Bock HealthCare GmbH