
14A1

DE	Gebrauchsanweisung	3
EN	Instructions for use	5
FR	Instructions d'utilisation	7
IT	Istruzioni per l'uso	9
ES	Instrucciones de uso	11
PT	Manual de utilização.....	14
SV	Bruksanvisning	16
PL	Instrukcja użytkowania	18
HU	Használati utasítás	20
CS	Návod k použití	22
TR	Kullanma talimatı.....	24

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English

Last update: 2015-03-19

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1 Application

The 14A1 Lock Set is to be used **solely** with Ottobock Silicone ArmLiners as an accessory component for upper limb prosthetic fittings.

2 Components

- (1) 14A110 Lock with Release Button
- (2) 14A107 Pin (23 mm)
- (3) 14A108 Lamination Ring including Lamination Protection

3 Description

The 14A110 Lock with Release Button attaches the 14Y1/14Y5 Ottobock Silicone ArmLiner to the prosthesis.

Pressing a release button on the lock disengages the 14A107 Pin from the lock.

Important!

During socket fabrication, please make sure that the residual limb is healthy and free of any swelling, redness, or similar physical impairments. The socket must be fabricated with the standard pin included in the delivery.

If the length of the 14A107 Pin is insufficient for engaging in the lock due to short-term fluctuations of residual limb volume, it can be replaced. A longer pin may be used only in exceptional cases for a short time. Use of a longer pin must be examined and evaluated by the responsible prosthetist in advance on a case-by-case basis. In addition, the prosthetist must make sure that the prosthesis fits securely and can be used without restrictions with a longer pin. Never use a longer pin to compensate for a poorly fitted socket. In such cases a new socket must be fabricated immediately.

The following lengths are available (not included in the delivery):

14A107=1 (28 mm)

14A107=2 (33 mm)

For myoelectric prostheses use only a standard pin. A change in pin length could alter the position of the electrodes with the result that the optimal control of the myoelectric prosthetic components can no longer be guaranteed.

4 Fabrication of Inner and Outer Socket

The 14A111 Dummy Set is required for socket fabrication.

Fabrication instructions are described in Technical Information Bulletin 3.3.1.

To finish the prosthesis, apply one drop of 635K13 Loctite or equal to the threads of the release button and screw it on completely by hand.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

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This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH