


617R11=PE/ 617R11=PP

DE Gebrauchsanweisung	3
EN Instructions for Use	5
FR Instructions d'utilisation	7
IT Istruzioni per l'uso	9
ES Instrucciones de uso	11
SV Bruksanvisning	13

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Date of the last update: 2014-09-16

- Please read this document carefully.
- Follow the safety instructions.

1 Description

The **617R11=PE** and **617R11=PP** products are thermoplastic prepregs 3 × 20 × 1200 mm in size. In combination with the appropriate thermoplastic sheets, these rectangular cross-sections can be incorporated into a custom orthosis in a single fabrication process. **They provide partial reinforcement** of the orthosis for optimal fitting. Furthermore, the non-reinforced areas can be shaped as needed, while also avoiding any distortion of parts.

The thermoplastic prepregs can be used, for example, for manufacture of:

- dorsiflex-assist ankle joint braces
- knee braces
- lower limb braces (even as test braces)
- body braces
- wrist braces


Fig. 1: Sectional view of thermoplastic pre-preg

The prepregs consist of thermoplastic (PE or PP) with unidirectional 'endless' glass fiber threads (Fig. 1). The glass fiber rovings are blended into the cover layer of the section and **unevenly** distributed along either side. The heat-weldable side contains a higher proportion of fibers than the opposite side to counteract distortion of the plastic sheet material. This side can be easily recognized by its 2-mm radius and "Ottobock" logo.

2 Processing

First, pull 2 stockinette layers over the model.

Then mark the positions of the preprep sections and measure their length (Fig. 2).

Cut off the preprep sections to the length and round their edges.

Remove the ends of the glass fiber rovings from the cut surface.

Clean and degrease the sheets and preprep sections with isopropyl alcohol before heating.

Now heat up the sheet cut-offs. The optimal processing temperature for PE (617R11=PE) is 190°C (375 F). The recommended plastic sheet material is 616T95 PE 200 (Fig. 3).

The processing temperature for PP (616R10=PP) is about 250°C (480 F). The recommended plastic sheet material is 616T20.

These temperatures, which exceed the normal processing conditions, allow for a shorter heating time, and are necessary to create the proper weld between the plastic materials. The plastic materials have been stabilized for such processing temperatures, and must be processed immediately after the processing temperature has been reached or they may become brittle due to their lengthy exposure to 250 °C (480 F) during the heating process. Because the prepreg sections heat more quickly than the sheets, they can be placed into the oven later than the sheets.

When the prepreg sections and sheet cut-offs are completely heated, position the self-adhesive profiles on the model (Fig. 4).

Make sure that the tops and bottoms are correctly oriented! (with "Ottobock" writing on the top)

Finally put the sheet cut-off on the model with the prepreg sections and turn on the vacuum (Fig. 5).

Once the vacuum is applied, the prepreg sections and the cut-off sheets will automatically be welded together.

Quick application of the vacuum will ensure a quality weld.

After the material cools down, continue the orthosis fabrication as usual (Fig. 6 and 7).

3 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

4 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH