
## C-Brace joint unit 17K01=\*

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	15

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## INFORMATION

Date of last update: 2022-03-24

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

### Note the instructions for use

This document contains the sections "Intended use" and "Safety", and does not replace the instructions for use. The instructions for use are available in digital form (PDF format) (see page 2 of this document). To view PDF files, you need an Internet browser or other program capable of opening this file format.

You may request a printed copy of the instructions for use under the link on page 2. This will be delivered to the specified address within 7 days after the order is received.

The product "C-Brace joint unit 17KO1=\*" is referred to as the product/component/orthosis/joint unit below.

Basic UDI-DI: 40644110000000017KO1LU

## 2 Intended use

### 2.1 Indications for use

The product is intended **solely** for orthotic fittings of the lower limbs.

### 2.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sports activities that involve jumping, sudden movements or fast sequences of steps (e.g. basketball, badminton, sporty riding).

Permissible ambient conditions are described in the technical data (Technical data).

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

### 2.3 Indications

- Unilateral or bilateral lower limb paresis or flaccid paralysis, e.g. due to post-polio syndrome, traumatic paresis, incomplete paraplegia.
- Physical prerequisites such as muscle status, joint mobility and possible axis deviations are crucial and proper control of the orthosis must be guaranteed.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations
- The existing muscle strength of the hip extensors and flexors must permit the controlled swing-through of the limb (compensation using the hip is possible).

### 2.4 Contraindications

#### 2.4.1 Absolute Contraindications

- Flexion contracture in the knee and/or hip joint in excess of 10°

- Knee varus/valgus malposition in excess of 10°
- Severe spasticity
- Body weight over 125 kg / 275 lbs

### 2.4.2 Relative Contraindications

- Moderate spasticity

## 2.5 Qualification

### 2.5.1 Qualification of the O&P professional

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

### 2.5.2 Qualifications of the therapist or nursing staff

The therapists or nursing staff must be trained in handling the product. Training must be carried out by the authorised orthopaedic technician.

## 3 Safety

### 3.1 Explanation of warning symbols

**⚠ WARNING** Warning regarding possible serious risks of accident or injury.

**⚠ CAUTION** Warning regarding possible risks of accident or injury.

### 3.2 Structure of the safety instructions

#### **⚠ WARNING**

#### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 3.3 General safety instructions

#### **⚠ WARNING**

#### **Non-observance of safety notices**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

#### **⚠ WARNING**

#### **Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

**⚠ CAUTION**

**Failure to observe warning/error signals**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (Status display on the control panel) and corresponding change in damping settings must be observed.

**⚠ CAUTION**

**Penetration of dirt and moisture into components of the orthosis**

- > Falling due to unexpected behaviour of the orthosis as a result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that neither solid particles, foreign bodies nor liquids penetrate the orthosis components.
- ▶ The orthosis components are protected against splashed water from all directions.
- ▶ However, the orthosis components are not protected against submersion, jets of water and steam.
- ▶ If water has penetrated the components of the orthosis, remove the straps and pads if possible and allow the components to dry.
- ▶ In case of a malfunction after drying, the orthosis must be inspected by an authorised Ottobock Service Centre.
- ▶ Always close the charging receptacle with the protective cap after disconnecting the charging cable.

**⚠ CAUTION**

**Independent manipulations on the joint unit and orthosis components**

Falling due to the breakage of load-bearing components or malfunction of the orthosis.

- ▶ Manipulations on the joint unit and orthosis components other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The joint unit and orthosis components may only be opened and repaired by authorised, qualified Ottobock personnel.

**⚠ CAUTION**

**Use of the product when battery charge level is too low**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the product if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

**⚠ CAUTION**

**Mechanical stress during transport**

- > Falling due to unexpected product behaviour as a result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Only use the transport packaging for transportation.

**⚠ CAUTION**

**Signs of wear and tear on the product components**

Falling due to damage or malfunction of the product.

- ▶ Regular service inspections (maintenance) are mandatory in the interest of patient safety and in order to maintain operating reliability and protect the warranty.

**⚠ CAUTION**

**Use of unapproved accessories**

> Falling due to product malfunction as a result of reduced interference resistance.

- > Interference of other electronic devices due to increased emissions.
- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the sections "Scope of delivery" (Scope of delivery) and "Accessories" (Accessories).

**3.4 Information on the Power Supply/Battery Charging**

**⚠ CAUTION**

**Charging the product with damaged power supply unit/charger/charger cable**

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

**⚠ CAUTION**

**Charging while wearing the product**

> Falling due to walking and getting caught on the connected battery charger.

- > Falling due to unexpected product behaviour because of changed damping behaviour.
- ▶ Inform the patient that the product may only be charged while sitting down.

**3.5 Information on Alignment/Adjustment**

**⚠ CAUTION**

**Improper assembly of the screw connections**

Falling due to breakage or loosening of the screw connections.

- ▶ Clean the threads before every installation.
- ▶ Apply the specified tightening torque values for installation (see the section "Technical data").
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

**⚠ CAUTION**

**Incorrectly secured screws**

Falling due to breakage of load-bearing components caused by screw connections coming loose.

- ▶ Before installing the joint unit and joint mechanism, the screws must be secured before they are tightened to the specified tightening torque (see the section "Technical data" Technical data).

**⚠ CAUTION**

**Operator errors during the adjustment process with the adjustment app**

Falling due to unexpected orthosis behaviour.

- ▶ Do not charge the orthosis battery during the adjustment process since the orthosis is not functional while the battery is being charged.
- ▶ The orthosis must not remain unattended during the adjustment process while connected to the adjustment software and being worn by the patient.
- ▶ Observe the maximum range of the Bluetooth connection and note that obstacles may limit this range.
- ▶ During the data transfer (tablet to orthosis), the user should sit still or stand with support.
- ▶ If only temporary changes to the settings are to be made while connected to the adjustment app, these changes must be reversed before the adjustment app is closed.  
It must also be ensured that the patient does not leave the range of the Bluetooth connection with settings that have been changed temporarily.
- ▶ Inform the patient immediately if the data connection is accidentally interrupted during the adjustment process.
- ▶ The connection to the orthosis must always be disconnected after adjustments have been completed.
- ▶ Participation in an Ottobock product training course is mandatory prior to the initial use of the product. Additional product training courses may be required to qualify for software updates.
- ▶ Correctly entering the patient data and especially the body weight is an important criterion for the quality of the fitting. If the values are too high, the orthosis may not switch to the swing phase. If the values are too low, the orthosis may trigger the swing phase at the wrong time.
- ▶ When entering the weight and dimensions, note the settings for the units. To change the units, see the section "Changing the units for weight and length" (Changing the units for weight and length).
- ▶ If the patient uses walking aids (e.g. crutches or walking canes) during the adjustment process, readjustment is required as soon as the patient stops using these walking aids.
- ▶ Do not pass on your personal access data.

**⚠ CAUTION**

**Use of unsuitable orthosis components**

Falling due to unexpected behaviour of the product or breakage of load-bearing components.

- ▶ Use the product only in combination with components listed in the section "Combination possibilities" (Combination possibilities).

**⚠ CAUTION**

**Disassembling the joint unit on a standing patient**

- > Risk of falls due to lack of support from the product.
- > Damage to the medial support.
- ▶ Disassemble the joint unit only on a sitting patient or with the orthosis taken off.

**⚠ CAUTION**

**Reusing the screws to assemble the joint unit on the frame after maintenance**

Falling due to breakage of screws that are reused to assemble the joint unit.

- ▶ The mounting screws of the joint unit have to be replaced each time after maintenance. This is required at least every 24 months or after reaching one million steps.
- ▶ To assemble the loaner unit and subsequently the joint unit after maintenance, use the screw set that is included with the loaner unit or can be ordered under article number 17KO1A=SET-1.

### 3.6 Notes on applying the product

#### **⚠ CAUTION**

##### **Foreign objects between the leg and orthosis shells**

Pressure points on the leg due to foreign objects at the contact points between the leg and orthosis shells.

- ▶ Smooth out wrinkles in the padding material and clothing.
- ▶ Check the leg for pressure points.

#### **⚠ CAUTION**

##### **Improper fit of the orthosis**

Falling/skin irritation due to insufficient support/hold by the orthosis shells.

- ▶ Take the orthosis off immediately and put it back on.
- ▶ Observe the instructions for applying and removing.

#### **⚠ CAUTION**

##### **Pinching skin in the area of the closures**

Injuries and swelling due to disruption of blood circulation to the skin in the area of the closures.

- ▶ Do not over-tighten the closures when applying the product.

#### **⚠ CAUTION**

##### **Volume fluctuations of the leg or problems with the fit of the orthosis**

Injuries, friction and local pressure caused by improper fit of the orthosis shells (too firm/too loose).

- ▶ Sensitivity problems and skin damage require special attention with regard to the fit. Instruct the patient to inspect the skin areas daily.
- ▶ The patient should consult a doctor or O&P professional if there are any signs of skin damage, even if they are only minor.
- ▶ If the improper fit is the result of weight gain or loss, new thigh shells and lower leg shells have to be fabricated based on a new plaster cast.
- ▶ Check the leg for pressure points.

### 3.7 Information on Proximity to Certain Areas

#### **⚠ CAUTION**

##### **Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

#### **⚠ CAUTION**

##### **Operating the product in very close proximity to other electronic devices**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.

- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

#### **⚠ CAUTION**

##### **Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Ensure that the patient is not in the vicinity of sources of strong magnetic and electrical interference during trial fitting (such as theft prevention systems, metal detectors...). If this cannot be avoided, ensure at least that the patient has a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.

#### **⚠ CAUTION**

##### **Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)**

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of strong magnetic fields.
- ▶ Make sure that the patient takes off the product before entering the room or area and stores the product outside this room or area.
- ▶ Damage to the product caused by exposure to strong magnetic fields cannot be repaired.

#### **⚠ CAUTION**

##### **Remaining in areas outside the allowable temperature range**

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Ensure that the patient is not in areas outside the permissible temperature range (Technical data) during trial fitting.

### **3.8 Information on Use**

#### **⚠ CAUTION**

##### **Improper mode switching**

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that the patient stands securely during all switching processes.
- ▶ Inform the patient that the changed damping characteristics have to be verified after switching and feedback from the acoustic signal emitter must be observed.
- ▶ Switching back to basic mode is mandatory once the activities in MyMode have been completed.
- ▶ If required, take the weight off the product and correct the switching.

#### **⚠ CAUTION**

##### **Risk of pinching in the joint flexion area**

- > Injury due to pinching of body parts.
- > Damage to garments due to pinching in the mechanism of the joint unit or medial support.
- ▶ Ensure that no body parts or garments are in this area when flexing the product.

**⚠ CAUTION**

**Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)**

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- ▶ Inform the patient that they must pay attention when pulsating vibration signals start. They indicate the risk of overheating.
- ▶ Inform the patient that the activity level has to be reduced immediately when these pulsating vibration signals begin, so the hydraulic unit can cool down.
- ▶ Inform the patient that full activity may be resumed after these vibration signals stop.
- ▶ If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. In this case, the product has to be inspected by an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Overloading due to unusual activities**

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sports activities that involve jumping, sudden movements or fast sequences of steps (e.g. basketball, badminton, sporty riding).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Using the product while switched off**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Prior to using the product, briefly press the **ⓘ** button on the control panel to check whether it is switched on. The button symbol lights up green **ⓘ** when the product is switched on.

**3.9 Notes on the movement patterns**

**⚠ CAUTION**

**Improper use of the stance function**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that the patient is standing safely when using the stance function and checks the lock of the knee joint before placing their full weight on the orthosis.
- ▶ Instruct the patient in the correct use of the stance function. Information on the stance function Stance function.

**⚠ CAUTION**

**Walking up stairs**

- > Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that the patient uses the handrail when walking up stairs and places most of the sole of the foot on the stair surface.
- ▶ Always step up with the less affected leg, followed by the leg with the product.
- ▶ If there is a connection to a foot component, support by bouncing with the product extended while walking up stairs is not allowable and therefore has to be avoided.
- ▶ Particular caution is required when carrying children up the stairs.

**⚠ CAUTION**

**Walking down stairs with moveable orthosis foot component**

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Inform the patient that the handrail always has to be used when walking down stairs, and that most of the sole of the foot has to be set onto the stair surface.
- ▶ The warnings and error signals have to be observed (Status display on the control panel).
- ▶ Notify the patient that resistance in the flexion and extension direction can change in case of warnings and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

**⚠ CAUTION**

**Walking down stairs with rigid orthosis foot component**

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Inform the patient that the handrail always has to be used when walking down stairs, and that the patient has to roll over the edge of the step with the middle of the shoe.
- ▶ The warnings and error signals have to be observed (Status display on the control panel).
- ▶ Notify the patient that resistance in the flexion and extension direction can change in case of warnings and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

**⚠ CAUTION**

**Walking down stairs and ramps**

Falling due to unexpected, increased stance phase damping during the transition from walking on the level to stairs or ramps, for example when the orthosis is in a MyMode.

- ▶ Inform the patient of the changed product behaviour, especially in the MyModes.
- ▶ Inform the patient that changed stance phase damping has to be verified before walking on stairs or ramps.

**⚠ CAUTION**

**Lack of support from the product when walking down stairs**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Inform the patient that it is essential to check whether corresponding flexion of the joint is possible before walking down stairs. If this is not possible, either the Cockpit app or turning the product off/on must be used to switch back to basic mode.

**⚠ CAUTION**

**Improper use of the MyMode "Training mode"**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that the patient is standing safely when using this MyMode and checks the lock of the knee joint before placing their full weight on the orthosis.
- ▶ Inform the patient that the knee joint is locked in the flexion direction in this MyMode. For information about this mode, MyMode "Training mode".
- ▶ Switching back to basic mode is mandatory once the activities in this MyMode have been completed.

**⚠ CAUTION**

**Improper use of the MyMode "Freeze position"**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that the patient is standing safely when using this MyMode and checks the lock of the knee joint before placing their full weight on the orthosis.
- ▶ Inform the patient that the knee joint is locked in the both the flexion and extension direction in this MyMode. For information about this mode, MyMode "Freeze position".
- ▶ Switching back to basic mode is mandatory once the activities in this MyMode have been completed.

**⚠ CAUTION**

**Lack of support from the product while sitting down**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Inform the patient that it is essential to check whether corresponding flexion of the joint is possible before sitting down. If this is not possible, either the Cockpit app or turning the product off/on must be used to switch back to basic mode.

**⚠ CAUTION**

**Quickly pushing the hip forward with the orthosis extended**

> Falling due to unexpected activation of the swing phase.

- ▶ Note that the joint may flex unexpectedly when the hip is pushed forward quickly while the orthosis is extended.
- ▶ With the patient, practice movement patterns where this behaviour can occur. If the patient participates in sports where this movement pattern can occur, configure a corresponding "user-defined MyMode" using the adjustment app.

**⚠ CAUTION**

**Changing the gait pattern influences activation of the swing phase**

> Falling due to unexpected activation of the swing phase.

- ▶ Inform the patient that changing the gait pattern influences activation of the swing phase. Adjustment by the O&P professional has to be repeated in this case.

### 3.10 Notes on the safety modes

**⚠ CAUTION**

**Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Using the product when it is defective is prohibited.

- ▶ The product must be inspected by an authorised Ottobock Service Centre.

#### **⚠ CAUTION**

##### **Safety mode cannot be deactivated**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ If safety mode cannot be deactivated by recharging the battery, a permanent error has occurred.
- ▶ Using the product when it is defective is prohibited.
- ▶ The product must be inspected by an authorised Ottobock Service Centre.

#### **⚠ CAUTION**

##### **Safety signal occurs (ongoing vibration)**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (Operating states/error signals) have to be observed.
- ▶ After the safety signal has been emitted, further use of the product is prohibited.
- ▶ The product must be inspected by an authorised Ottobock Service Centre.

#### **⚠ CAUTION**

##### **Using the product in safety mode**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (Operating states/error signals) have to be observed.
- ▶ Particular caution is necessary when using a bicycle without a freewheel (with a fixed gear).

### **3.11 Information on the use of a mobile device with the cockpit app**

#### **⚠ CAUTION**

##### **Improper use of the mobile device**

Falling due to altered damping behaviour as a result of unexpected switching into a MyMode.

- ▶ Use the instructions for use (user) to instruct the patient on the proper handling of the mobile device with the Cockpit app.

#### **⚠ CAUTION**

##### **Independently applied changes or modifications made to the mobile device**

Falling due to altered damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Do not make any independent changes to the hardware of the mobile device on which the app is installed.
- ▶ Do not make any independent changes to the software/firmware of the mobile device that are not included in the update function of the software/firmware.

#### **⚠ CAUTION**

##### **Improper mode switching with the mobile device**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that the patient stands securely during all switching processes.
- ▶ Inform the patient that the changed damping characteristics have to be verified after switching, and feedback from the acoustic signal emitter and the mobile device display must be observed.
- ▶ Switching back to basic mode is mandatory once the activities in the MyMode have been completed.

A series of horizontal lines spaced evenly down the page, providing a template for writing or drawing.

Legal Manufacturer of C-Brace joint unit 17KO1=\* and C-Brace Setup 560X17-  
Otto Bock Healthcare Products GmbH  
Brehmstraße 16 · 1110 Wien · Austria  
T +43-1 523 37 86 · F +43-1 523 22 64  
info.austria@ottobock.com · www.ottobock.com