

Academy Tech Tips

Pairing the Remote Control (4X250) with a Compact Knee

When you first use the 4X250 remote for the Compact knee, it will need to be paired with the knee. A paper clip works well to accomplish this task.



Given that only one joint should be paired with a remote control, it must be ensured that no other joint is within 3 meters during the pairing.

Looking at the face of the remote, you will notice there is a small hole at the top of the remote. Holding the remote at least 12 to 18 inches away and no more than 3 feet away from the knee, insert the paper clip into the hole and you will feel a button. Apply slight pressure and hold the button down until you hear 5 beeps indicating that the remote is now paired to the knee.

If you have any further questions, please contact Ottobock at 800 328 4058 and ask to speak with a member of our Professional and Clinical Services department.