

## 3R80, 3R80=ST

<b>DE</b> Gebrauchsanweisung (Anwender) .....	3
<b>EN</b> Instructions for use (user) .....	7
<b>FR</b> Instructions d'utilisation (utilisateur) .....	11
<b>IT</b> Istruzioni per l'uso (utente) .....	15
<b>ES</b> Instrucciones de uso (usuario) .....	19
<b>PT</b> Instruções de utilização (utilizador) .....	23
<b>NL</b> Gebruiksaanwijzing (gebruiker) .....	27
<b>SV</b> Bruksanvisning (användare) .....	32
<b>PL</b> Instrukcja używania (użytkownik) .....	35
<b>HU</b> Használati útmutató (felhasználó) .....	40
<b>CS</b> Návod k použití (uživatel) .....	44
<b>RO</b> Instrucțiuni de utilizare (utilizator) .....	48
<b>HR</b> Upute za uporabu (korisnik) .....	52
<b>TR</b> Kullanım kılavuzu (kullanıcı) .....	56
<b>RU</b> Руководство по применению (для пользователей) .....	59
<b>UK</b> Інструкція щодо використання (користувачі) .....	64
<b>JA</b> 取扱説明書(使用者) .....	68
<b>ZH</b> 使用说明书(用户) .....	72


## 1 Technical data

**1** max. 150 kg  
max. 330 lbs



**3** 3R80  
1196 g



**4** 3R80=ST  
1191 g

## 2 Using the lock

### 1 Activating the lock

### 2 Deactivating the lock

## INFORMATION

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### INFORMATION

Date of last update: 2025-02-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The "3R80/3R80=ST" prosthetic knee joint is referred to as the product in the following. This document is intended for users of the product.

## 2 Safety-related information

### CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the maximum lifetime of the product (see page 9).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Observe the operating conditions of the product (see page 8) and do not expose the product to prohibited environmental conditions (see page 8).
- ▶ Do not use the product if it is damaged or in a questionable condition. Clean the product if it gets dirty. If necessary, have it inspected, repaired or replaced by the manufacturer or a specialist workshop.
- ▶ To avoid the risk of pinching, do not reach into the joint mechanism.
- ▶ Do not touch hot prosthetic components.
- ▶ If there are changes in functionality, reduce all activities so the prosthetic components can cool down.
- ▶ Have the prosthetic components checked by qualified personnel in the event of functional changes.
- ▶ In the event of leaks, stop using the product and have it repaired.
- ▶ The product may only be used by one person; reuse on another person is prohibited.
- ▶ Only use the approved cleaning agents and disinfectants.
- ▶ Observe the instructions for cleaning and care (see page 9).
- ▶ Follow the maintenance instructions (see page 10).

### Signs of changes in or loss of functionality during use

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

## 3 Product description

The product has the following key features:

- Monocentric prosthetic knee joint with rotation hydraulics
- Product components for stance phase stability:
  - Adjustable stance phase flexion resistance (hydraulic damping)

- Adjustable trigger threshold for switching between swing phase resistance to stance phase flexion resistance
- Lock (activated and deactivated by the user)
- Product components for swing phase control:
  - Adjustable swing phase flexion resistance (hydraulic damping)
  - Adjustable swing phase extension resistance (hydraulic damping)
  - Hydraulic unit with extension assist function (the extension assist function can be enhanced by installing an additional spring)

## 4 Intended use

### 4.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 4.2 Operating conditions

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

Allowable field of application	
Recommended mobility grade: <b>3 + 4</b> (see fig. 1 item ①)	Everyday prosthesis
Allowable body weight: <b>≤ 150 kg</b> (see fig. 1 item ②)	

### 4.3 Environmental conditions

Storage and transport
<ul style="list-style-type: none"> <li>• Storage temperature: -20 °C to +60 °C</li> <li>• Relative humidity: 20 % to 90 %</li> <li>• No mechanical vibrations or impacts</li> </ul>

#### INFORMATION

Clean the product after contact with humidity/chemicals/solids in order to avoid increased wear and tear and damage (see page 9).

Allowable environmental conditions
<ul style="list-style-type: none"> <li>• Temperature range: -10 °C to +45 °C</li> <li>• Relative humidity: 20 % to 90 %</li> </ul>
<ul style="list-style-type: none"> <li>• Contact with soap solution (e.g. when taking a shower)</li> <li>• Contact with perspiration</li> <li>• Submersion in fresh water/chlorinated water (e.g. swimming pool): maximum immersion depth 2 m</li> </ul>
<ul style="list-style-type: none"> <li>• Submersion in salt water with a low salt content (14 days/year, 1 hour/day)</li> <li>• Permissible salt content: max. 3.5 % (e.g. ocean)</li> <li>• Maximum immersion depth swimming pool: 1 m</li> <li>• Maximum immersion depth sea: 0.5 m</li> </ul>
<ul style="list-style-type: none"> <li>• Contact with dust, foam cover particles, drifting sand (e.g. walking on the beach)</li> </ul>
<ul style="list-style-type: none"> <li>• Contact with salty, condensing air (e.g. in a 40 °C steam sauna)</li> </ul>
<ul style="list-style-type: none"> <li>• UV resistant</li> </ul>

### Unallowable environmental conditions

- Contact with liquid-binding particles (e.g. talcum)
  - Contact with large amounts of sand and dust (e.g. digging, kneeling in sand, construction site work)
  - Contact with acids
  - Contact with urine
- 
- Submersion in salt water with excessive salt content
  - Salt content: > 3.5 % (e.g. brine bath, Dead Sea)
- 
- Cleaning agents and disinfectants with solvents, chlorine and phosphate
- 
- High water pressure (e.g. diving, jumping into water)

## 4.4 Lifetime

This prosthetic component has been load tested by the manufacturer according to ISO 10328. The maximum lifetime is 5 years.

## 5 Use

**Use lock** (see fig. 2)

### CAUTION

#### Locking function

Risk of injury due to unintentional activation or deactivation of the locking function

- ▶ Check that the joint is functioning correctly each time you activate or deactivate the lock.
- ▶ Make sure that cosmetic covers and clothing cannot activate the lock.

To activate or deactivate the lock, put weight on the prosthetic knee joint, push against the extension stop and activate the lock lever.

- 1 Activate: Push the lock lever down (e.g. for relaxed standing)
- 2 Deactivate: Pull the lock lever up (e.g. for walking)

**Use in wet areas**

### CAUTION

#### Use in wet areas

Risk of injury due to unlocked joint or walking too fast on wet surfaces

- ▶ Lock the joint before using it in wet areas.
- ▶ Walk carefully on wet or slippery surfaces.

- 1) Observe the environmental conditions prior to use in wet areas (see page 8).
- 2) Activate the lock.
- 3) Walk carefully on wet and slippery surfaces.
- 4) Deactivate the lock after leaving wet areas.
- 5) Dry the joint.

## 6 Cleaning

### Permissible cleaning agents and disinfectants

#### Requirements:

Free of solvents, chlorine and phosphates

> **Required materials:** 2 microfibre cloths, isopropyl alcohol 634A58, tap water

- 1) **In case of light soiling (e.g. small deposits):** Moisten the microfibre cloth with water.
- 2) **In case of heavy soiling (e.g. tough deposits):** Moisten the microfibre cloth with isopropyl alcohol.

- 3) **NOTICE! Verify the material compatibility when using isopropyl alcohol. Do not remove lubricant from bearing points.**
- 4) Clean the product with the moistened cloth.
- 5) Dry the product with a dry microfibre cloth.
- 6) Allow to air dry in order to remove residual moisture.

#### **Cleaning after contact with salt water**

> **Required materials:** Microfibre cloth, tap water

- 1) **NOTICE! Remove salt crystals from the prosthetic knee joint to prevent corrosion.**
- 2) **Cleaning the prosthetic knee joint:** Immerse in tap water.
- 3) **Rinse out dirt and salt crystals:** Move the joint in the water (flex and extend several times).
- 4) **Remove residual water:** Remove the joint from the water and allow it to drain.
- 5) **Drying:** Dry with a microfibre cloth and allow to air dry.
- 6) **Perform a functional test:** Test the joint for proper function (flexing, extending, locking).
- 7) **Checking by qualified personnel:** After maximum use in salt water (14 days/year, 1 hour/day), have it checked by qualified personnel.

#### **Cleaning with Disinfectants**

> **Required materials:** 2 microfibre cloths, colourless disinfectant

- 1) **NOTICE! Verify material compatibility when using disinfectants. Do not remove lubricant from bearing points.**
- 2) Apply disinfectant to the microfibre cloth.
- 3) Disinfect the product with the moistened cloth.
- 4) Dry the product with a dry microfibre cloth.
- 5) Allow to air dry in order to remove residual moisture.

## **7 Maintenance**

- ▶ **NOTICE! Do not lubricate and grease the prosthetic joint.**
- ▶ **NOTICE! Repair work must be performed exclusively by manufacturer service.**
- ▶ Check (inspect) the prosthetic components after 30 days of initial use.
- ▶ Have annual safety checks carried out.
- ▶ **NOTICE! The protector may partially open or become damaged when kneeling. If contact is unavoidable, close the protector again after kneeling. If the protector no longer closes, please have it checked by qualified personnel.**

## **8 Disposal**

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

## **9 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **9.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### **9.2 CE conformity**

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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