Oslo Bock®

QUALITY FOR LIFE

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Dorso Arexa – Fitting made easy

With valuable practical tips

Fitting of the Dorso Arexa made easy

The Dorso Arexa covers a fitting spectrum for waist circumferences between 60 and 150 cm with one single size. It can therefore be adjusted to the needs and requirements of each individual patient.

Here are a few practical tips and tricks in order to make it easy for you to perform the fitting.

1

2

3

Height adjustment

Patient-specific height adjustment is performed via a continuously adjustable sternal and symphysis support.

- Please loosen the screws in the base plate, adjust the sternal and symphysis supports continuously to the required height and re-tighten the screws (1).
- Always observe the maximum extension length (see marking) (2).

Fitting the sternal and symphysis supports

The sternal and symphysis supports can be fitted through uniform bending of both tubes by hand.

- Leave the sternal and symphysis supports in the base plate.
- Put the base plate on a table edge and bend the tubes (3).

1

2

3

4

5

Fitting the circumference

The circumference is adjusted by removing or replacing the extension elements as well as by shortening the waistband. The following steps are necessary in order to fit the Dorso Arexa to your patient.

1 Preparation

- Attach the two-part waistband to the back support with screws (countersunk head screws); the waistband is adjustable at three height levels (1).
- Loosen the ratchet closures on the waistband (2).
- Now slide the back support underneath the patient's back (3).
- Push the ratchet closures onto the toothed belt and leave one tooth exposed on each side (4).
- Place the base plate on the patient's stomach (5).

1

2

3

2 Shortening the waistband

A patient's waist circumference can be different between standing and lying down. This especially applies to patients with larger waist circumferences. When shortening the waistband take care to ensure that you don't cut it too short if you adjust the orthosis while the patient is lying down.

Tip:

For patients with larger waist circumferences we can therefore recommend to leave two additional waistband segments on each side of the waistband.


- Measuring the patient's waist circumference.

Tip: There are two possible measurement approaches:

1. Measure the waist circumference from the end of one ratchet closure to the end of the other ratchet closure (1).
- Shorten the waistband at both ends to one half of the measured circumference (3).
2. Measure the patient's full circumference and subtract 18 cm from the measured value (Distance between the ratchet closures) (2).
- Shorten the waistband at both ends to one half of the identified circumference (3).

1

2

3

4

5

6

3 Attaching the ratchet closures

- Make sure that the ratchet closures are flush with the waistband end and snap the ratchet closures as well as the toothed belt cover in place (1).
- After making final adjustment to the seat and the fit, secure the closures with the screws (2).

Tip:

Hand-tighten the screws. If tightened too much, they will be difficult to move later.

4 Shortening and attaching the waistband pads

- Shorten the pads to the required length (3).
- Attach the waistband pad to the waistband using clips (4).

5 Removing the extension elements

- Remove the same number of unnecessary extension elements from both sides (5). In doing so, please note that the strap guide element is always the final element (6).

1

2

3

4

5

Applying the Dorso Arexa

Tip:

For very slim patients, replacing the dark grey (hard) extension elements with the light grey (soft) ones is recommended. To do so, loosen the dark grey elements which are attached directly to the base plate and insert the light grey elements (1). This results in a higher degree of flexibility (2), so that the Dorso Arexa rests firmly and optimally against the body even for slim patients.

- Shorten the base plate pad to the required length (3).

Tip:

To insert the waistband in the strap guide element please proceed as follows:

Set the waistband at an angle on the lower small raised part of the strap guide element (4). Then click the waistband onto the upper small raised part (5).

Otto Bock®

QUALITY FOR LIFE

Otto Bock HealthCare GmbH

Max-Nieder-Strabe 15 · 37115 Duderstadt/Germany · Phone +49 5527 848-1424 · Fax +49 5527-72330 · export@otobock.de · www.otobock.com