ottobock.

Quickstart

E-MAG Active 17B203

Quality for life

E-MAG Active

Gangzyklus

Gait cycle

Loopcyclus

Cycle de la marche

Ciclo de marcha

Ciclo de marcha

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Gångcykel

Cykl chodu

Cyklus chůze

Járásciklus

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Gangcyklus

Gangesyklus

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EN E-MAG Active Patient Information

Your orthotist fitted you with a complete leg orthosis.

The E-MAG Active knee joint system is a special feature of your orthosis. This joint promotes a natural gait pattern and offers you a high degree of safety in the stance phase. The use of the orthosis depends on whether you were previously fitted with a knee joint that included a so-called Swiss lock or a free-motion knee joint with posterior positioning, or if you have not had a prior fitting at all. In any case, ask your orthotist or therapist to instruct you regarding the correct use of the orthosis.

The purpose of this Quick Start Guide is to provide a brief explanation of key points you need to consider when using the E-MAG Active as well as its features. An electronic control unit incorporated into the E-MAG Active measures your leg position on the fitted side as you walk. When you are in the stance phase – which is the phase from heel impact to toe-off – the joint remains locked. Only after the toes leave the ground does the joint open automatically, allowing you to swing the leg forward freely. Shortly before the heel strikes the ground again, the joint is once again locked when the extension position is reached. **1** The E-MAG Active therefore results in a more physiological gait pattern and helps you save energy while walking. In order to use the orthosis safely, it is essential to recognise the limits of the joint technology. It depends on proper control for trouble-free functionality.

What do you need consider while walking?

The most important moment during the gait cycle is the point of toe-off, which is when the joint is unlocked so the leg can swing forward freely. The joint is designed so that the knee is securely locked during the stance phase. However, this design also requires the load to be taken off the lock (extension) before initiating the swing-through phase. If the load is not taken off the joint, it cannot unlock and remains engaged. In this case the leg cannot swing forward! **2** To facilitate the use of the orthosis your orthotist may have selected an ankle joint, or angle between the lower leg and foot, that enables you to reach the so-called joint extension position.

How to practice using the orthosis

Practice using your orthosis at your orthotist's or therapist's facility. In most cases parallel bars are available, allowing you to walk freely while supporting yourself with both hands at any time. Treadmills are also highly suitable, as the repetitive nature of this type of training is very effective. Once you can control the rhythm required for movement yourself, you can continue training with the aid of a crutch. However, please keep in mind that this rhythm itself requires practice and cannot be executed perfectly from the first step. You will require several training sessions before you can easily release the lock.

3 Operating the electronic lock / unlock function

A Permanent locking of the joint.

You can use Button A to enhance your safety on uneven surfaces by deactivating the joint. It remains locked and no longer responds to your pattern of movement. In order to do so, push the top button as shown in illustration 3 on page 5.

B One-time unlocking of the joint.

Button B allows you to electronically unlock the joint so you can sit down. However, depending on your level of expertise, you can also initiate automatic unlocking. You should first practice this procedure with your orthotist or therapist. Push the bottom button as shown in illustration 3 on page 5. Remember that the load also needs to be taken off the joint prior to sitting down (as described in the section “What do you need to take into account while walking?”).

ATTENTION: For safety reasons, both buttons require a “double click” for activation, just like using a computer mouse.

4 Control of manual unlocking

- 1 To permanently unlock the joint, the lock can be released manually. This allows the joint to be used for special activities such as cycling.
- 2 The lock must be reactivated for active support during the stance phase.

Power supply - how to install and replace the battery

The battery should last at least one day, provided you charge it regularly using the supplied charger. Charging the battery after each use of the orthosis is recommended. Please do not hesitate to ask your orthotist if you have any questions. The electronic control unit emits a warning tone when the battery is low! The system switches to safe mode once the battery is drained. The joint then remains locked at all times (see „Operating the mechanical release“).

Maintenance and care

Please bring the orthosis to your orthotist for maintenance at least once every 6 months. Your orthotist has a maintenance schedule with a checklist for the inspection of all relevant components.

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Otto Bock HealthCare GmbH