

Otto Bock®

QUALITY FOR LIFE



## 4R119 Lamination Anchor

Adapter for optimizing prosthetic alignments

Correct alignment is essential to allow optimal functionality of a prosthesis and has a significant influence on the quality of the amputee's fitting!

Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt/Germany · Phone +49 5527 848-0 · Fax +49 5527 848-1524 · export@ottobock.com · www.ottobock.com


# 4R119 Lamination Anchor

## The 4R119 Lamination Anchor supports you on the road to an optimized prosthetic alignment as follows:

- Suitable for use in the transfemoral area
- Dependent on existing clearance the 4R119 establishes a direct connection to the Otto Bock knee joint with adjustment pyramid
- The Lamination Anchor 4R119 is provided with a dorsally positioned angled arm requiring less bending work. This allows for easier installation and provides an optimized prosthetic alignment through consideration of the following aspects:
  - A properly positioned socket connection creates the basis for correct prosthetic alignment
  - The stump position has to be considered for positioning
  - Frequently the stump position is characterized by a flexion angle of 3–5° (observe the individual situation, e.g. with hip flexion contraction)
- When correctly positioned the 4R119 considers the flexion angle, provides for proper positioning of the knee joint and allows for a proper prosthetic alignment in relation to the alignment reference line



4R119 Lamination Anchor with rotatable pyramid receiver and angled arm




### Advantages at a glance

Prosthetist	Amputee
<ul style="list-style-type: none"> <li>• Timesaving               <ul style="list-style-type: none"> <li>- Bending of dorsal anchor arm is less often/rarely required</li> <li>- Less time-consuming fitting sessions (e.g. faster adjustment of knee joint)</li> </ul> </li> <li>• Optimized functionality of prosthesis               <ul style="list-style-type: none"> <li>- Customer satisfaction</li> <li>- Reduced risk of occurrence of liability cases</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>• Relief in daily activities               <ul style="list-style-type: none"> <li>- Correct alignment allows to better take advantage of the functionality of installed components</li> </ul> </li> <li>• Amputee benefits from more safety offered by a prosthesis with optimized functionality</li> </ul>

### Technical data

<b>Article number</b>	4R119
<b>Weight</b>	165 g
<b>System height</b>	44 mm
<b>Material</b>	INOX*
<b>Max. body weight</b>	150 kg

\*Stainless steel



max. 150 kg