# Anwenderinformation

DE

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# User Information

EN

## **Dyna Ankle – Dynamic Ankle Orthosis designed by Prof. Dr. Grifka**

The dynamic ankle orthosis – Dyna Ankle – was designed for the functional treatment of sprains and ruptures of the deltoid ligament. The placement of the orthosis and its dynamic tether promotes healing by protecting the injured ligament from excess stress.

Please be sure to follow the instructions of the physician who has prescribed the orthosis. Avoid walking or moving your ankle without the protection of wearing the orthosis. Supportive shoes, such as sports shoes, are recommended for use with this orthosis.

### **To apply the orthosis, proceed as follows:**

Place the foot in the plastic shell, sliding your heel to the back (fig. 1).

Pull the lower hook and pile strap from the inner side of your foot across the top of the foot, through the outer plastic guide and tighten the strap so that the heel will snugly fit to the shell (fig. 2 and 3).

Pull the upper hook and pile strap through the plastic loop and fasten it (fig. 4).

Fasten the hook and pile strap on the forefoot as well (fig. 5).

Slightly stretch the elastic functional strap and fasten it (fig. 6). The shape of the orthosis and the tension of the strap ensure correct positioning of the joint.

When putting on the shoe, open the laces completely to ensure easy donning of the orthosis.

For new injuries, the orthosis should be worn day and night – generally for 6 weeks. If swelling increases, or pain develops, contact your physician immediately.

After the doctor says the injury has healed, you may use the Dyna Ankle orthosis during vigorous activities for added protection.

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