
## 8385 Genu Therma Tendon 8386 Genu Therma Chondro

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	9
<b>FR</b> Instructions d'utilisation .....	13
<b>IT</b> Istruzioni per l'uso .....	18
<b>ES</b> Instrucciones de uso .....	23
<b>NL</b> Gebruiksaanwijzing .....	28
<b>SV</b> Bruksanvisning .....	33
<b>DA</b> Brugsanvisning .....	38
<b>NO</b> Bruksanvisning .....	42
<b>FI</b> Käyttöohje .....	47
<b>PL</b> Instrukcja użytkowania .....	51
<b>CS</b> Návod k použití .....	56
<b>TR</b> Kullanma talimatı .....	61
<b>ET</b> Kasutusjuhhis .....	65
<b>RU</b> Руководство по применению .....	70
<b>ZH</b> 使用说明书 .....	75

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## INFORMATION

Last update: 2014-07-29

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 8385 Genu Therma Tendon and 8386 Genu Therma Chondro knee supports.

## 2 Application

### 2.1 Intended Use

The support is intended **exclusively** for fittings of the lower limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

### 2.2 Indications

- Patellofemoral pain syndrome (e.g. chondromalacia patellae, chondro-pathia patellae, retropatellar osteoarthritis)
- Patellar tip syndrome (jumper's knee)
- Osgood-Schlatter disease

Indications must be determined by the physician.

### 2.3 Contraindications

#### 2.3.1 Absolute Contraindications

Not known.

#### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

### 2.4 Effects

The support relieves pressure on the patellar tendon insertion as well as the retropatellar joint surface and can relieve pain.

## 3 Safety

### 3.1 Explanation of Warning Symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General Safety Instructions

 **CAUTION**

**Using the product without following the instructions for use**

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

 **CAUTION**

**Use of a product with neoprene**

Allergic reactions due to neoprene or heat allergies

- ▶ Do not wear products with neoprene for more than 3 – 4 hours without interruption (unless explicitly instructed otherwise by the doctor).
- ▶ Do **not** use the product in case of a known neoprene or heat allergy.
- ▶ **Inform the patient.**

 **CAUTION**

**Contact with heat, embers or fire**

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

 **CAUTION**

**Reuse on other patients and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

#### NOTICE

##### **Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Only use this product carefully and for its intended purpose.
- ▶ Do not make any improper changes to the product.
- ▶ **Inform the patient.**

#### NOTICE

##### **Contact with salves, lotions or other products that contain oils or acids**

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

## 4 Handling

#### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Size Selection

- 1) Measure the circumference of the lower leg 15 cm below the centre of the patella.
- 2) Determine the size of the support (see size chart).

### 4.2 Fitting and Application

#### CAUTION

##### **Incorrect or excessively tight application**

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the support is applied properly and fits correctly.
- ▶ Do not use the product in case of a known heat allergy.

- ▶ When in doubt, do not continue to use the product if skin irritation occurs.
- ▶ **Inform the patient.**

#### NOTICE

#### Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

### 8385 Genu Therma Tendon

> Unfasten the hook-and-loop closure of the support.

> The patient is sitting on the front edge of a chair.

1) Wrap the support around the knee from behind (see Fig. 1).

**INFORMATION: The support must be applied so that the small cushion on the inner side is located in the upper part of the support.**

2) Position the small cushion precisely under the kneecap.

3) Close the hook-and-loop closure (see Fig. 2).

### 8386 Genu Therma Chondro

> Open all hook-and-loop straps on the support.

> The patient is sitting on the front edge of a chair.

1) Hold the support by the top edge on both sides and pull it over the knee until the patella cutout is centred on the patella (see Fig. 3).

2) Position the small cushion precisely under the kneecap.

3) Fasten the hook-and-loop straps first above and then below the hollow of the knee (see Fig. 4).

## 4.3 Cleaning

#### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

1) Fasten all hook-and-loop closures.

2) Recommendation: Use a laundry bag or net.

- 3) Wash the support in warm water at **40 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 5 Disposal

Dispose of the product in accordance with national regulations.

## 6 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH