


## 3R67

<b>DE</b> Gebrauchsanweisung .....	12
<b>EN</b> Instructions for use .....	19
<b>FR</b> Instructions d'utilisation .....	27
<b>IT</b> Istruzioni per l'uso .....	34
<b>ES</b> Instrucciones de uso .....	42
<b>PT</b> Manual de utilização .....	50
<b>NL</b> Gebruiksaanwijzing .....	57
<b>SV</b> Bruksanvisning .....	65
<b>DA</b> Brugsanvisning .....	72
<b>NO</b> Bruksanvisning .....	80
<b>FI</b> Käyttöohje .....	87
<b>PL</b> Instrukcja użytkowania .....	94
<b>HU</b> Használati utasítás .....	102
<b>CS</b> Návod k použití .....	109
<b>HR</b> Upute za uporabu .....	116
<b>TR</b> Kullanma talimatı .....	124
<b>RU</b> Руководство по применению .....	131
<b>ZH</b> 使用说明书 .....	139


**1** 5 Lieferumfang  
**1** 5 Scope of delivery

**3R67**


4Z59=5.5X6



710H10=2X3


506G3=M6X12


501S42=M6X18

647G1098

**2** **6 Technische Daten**  
**1** **6 Technical data**


510 g

65 g

151 mm


79 mm


① Systemhöhe (SH)  
 System height (SH)

97 mm

32 – 65 mm

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### **INFORMATION**

Last update: 2017-04-19

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

## INFORMATION

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings.

### 1.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 1.2 Meanings of pictograms in the illustration

	Numbering for the illustrations		Reference to the number of the corresponding section
	Numbering for the parts of an illustration		

## 2 Product description

The product (3R67) has the following key features:

- Polycentric prosthetic knee joint
- Product components for stance phase stability
  - Polycentric kinematics
- Product components for swing phase control:
  - Separately adjustable, hydraulic damping for extension and flexion

## 3 Intended use

### 3.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

### 3.2 Area of application

#### **CAUTION**

##### **Excessive strain on the product**

Fall due to breakage of load-bearing components

- ▶ Only use the product according to its allowable field of application.

##### **Allowable amputation types**

Knee disarticulation, Transfemoral amputation

##### **Allowable field of application**

Children's prosthesis (Everyday prosthesis)

##### **Positioning in the prosthesis**

Only according to the sections and illustrations in the instructions for use

##### **Allowable body weight**

See item number ④ in Figure ② "Technical data" on page 3

##### **Allowable body height**

See item number ⑤ in Figure ② "Technical data" on page 3

##### **Allowable foot size**

See item number ⑥ in Figure ② "Technical data" on page 3

### 3.3 Combination possibilities

#### **⚠ CAUTION**

##### **Failure to observe the manufacturer requirements regarding combination possibilities**

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- ▶ Contact the manufacturer with questions.

#### **INFORMATION**

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

### 3.4 Environmental conditions

#### **⚠ CAUTION**

##### **Use under unallowable environmental conditions**

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

#### **Allowable environmental conditions**

##### **Use**

Temperature range: -10 °C – +45 °C

##### **Transportation and storage**

Temperature range: -20 °C – +60 °C

Relative humidity: ≤100%

Air pressure: 750 mbar<sup>(1)</sup> – 1013 mbar<sup>(2)</sup>

<sup>(1)</sup>: elevation approx. 2500 m = min. aircraft cabin pressure

<sup>(2)</sup>: sea level

##### **Miscellaneous**

Penetration of particles:

Particles of cosmetics

#### **Unallowable environmental conditions**

Not waterproof

Not corrosion-resistant

##### **Penetration of particles**

Dust, sand, highly hygroscopic particles (e.g. talcum)

##### **Transportation and storage**

Mechanical vibrations or impacts

### 3.5 Reuse and service life

#### **⚠ CAUTION**

##### **Reuse on another patient**

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

#### **⚠ CAUTION**

##### **Exceeding the service life**

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life defined in this section is not exceeded.

This prosthetic component has been load tested by the manufacturer. Maximum service life is 2 years.

### 4 General safety instructions

#### **⚠ CAUTION**

##### **Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

#### **⚠ CAUTION**

##### **Overheating of the hydraulics due to excessive strain**

Burns, injuries caused by falling due to changes in functionality and damage to prosthetic components

- ▶ Do not touch overheated prosthetic components.
- ▶ Reduce all activities if there are changes in functionality, in order to allow the overheated prosthetic components to cool down.
- ▶ If overheating or changes in functionality occur, have the prosthesis checked by authorised qualified personnel.

#### **⚠ CAUTION**

##### **Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

##### **Signs of changes in or loss of functionality during use**

Noticeable changes in functionality may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

### 5 Scope of delivery

The scope of delivery is shown in Figure **1** on page 2.

In the illustration, a reference number shown in bold stands for a package in which the product components shown in the frame can be reordered.

Only product components with reference numbers in the illustration can be reordered separately.

## 6 Technical data

The technical data are shown in Figure **2** on page 3.

①: System height (SH); ②: Proximal system height (PSH); ③: Distal system height (DSH)

## 7 Preparation for use

### 7.1 Information on fabrication of a prosthesis

#### **⚠ CAUTION**

##### **Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

#### **⚠ CAUTION**

##### **Initial use of the prosthesis by the patient**

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

### 7.2 Bench alignment of the prosthesis

- ▶ **CAUTION!** To enable the patient to stand safely, align the prosthesis using the alignment images for the sagittal plane (Figure **3** on page 4) and for the frontal plane (Figure **4** on page 5) as well as the instructions for use of all prosthetic components that were used.

#### **Installing the tube adapter ( **4** on page 5 – ①)**

#### **⚠ CAUTION**

##### **Incorrect mounting of the tube**

Risk of injury due to breakage of load-bearing components

- ▶ Clean the contact surfaces of the tube and the tube adapter using a degreasing cleaning agent.
  - ▶ Position the tube in the tube adapter so the minimum insertion depth is met.
- 
- ▶ Use the smallest value for the allowable insertion depth as the minimum insertion depth.
    - ②: Allowable insertion depth:  $\geq 30$  mm
    - ③: Unallowable insertion depth:  $< 30$  mm

#### **⚠ CAUTION**

##### **Improper assembly of the screw connections**

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

- ▶ **CAUTION!** Use thread lock at the latest when finishing the prosthesis according to Figure [4] on page 5.

### 7.3 Optimising the static alignment

- ▶ To optimise the static alignment, position the patient on the measuring device.
- ▶ In the sagittal plane, check the course of the load line relative to the marked alignment reference point and correct this if needed (see Figure [5] on page 6).
- ▶ In the frontal plane, check the course of the load line relative to the marked points (prosthetic foot, prosthetic knee joint, spina iliaca anterior superior) and correct if needed (see Figure [6] on page 7).

### 7.4 Optimising during dynamic trial fitting

#### **⚠ CAUTION**

##### **Adjusting the Settings**

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

##### **Checking the factory settings ([7] on page 8 – ①)**

- ▶ Prior to the dynamic trial fitting, restore the factory settings (see Figure [7] on page 8).
  - ②: Flexion damping – ③: Middle setting
  - ④: Extension damping – ⑤: Minimum setting
- ▶ Optimise the prosthesis settings in the sequence according to this section.

##### **Setting flexion damping ([8] on page 9 – ①)**

- ▶ For setting flexion damping, see Figure [8] on page 9.
  - ②: Increase flexion damping
  - ③: Decrease flexion damping
- ▶ **CAUTION!** The lower leg of the prosthesis must not swing through too far in the dorsal direction at the beginning of the swing phase.  
Set flexion damping so that the lower leg of the prosthesis swings through in the dorsal direction according to the patient's requirements, but reaches full extension in a timely manner for the next heel strike.

##### **Setting extension damping ([9] on page 10 – ①)**

- ▶ For setting extension damping, see Figure [9] on page 10.
  - ②: Increase extension damping
  - ③: Decrease extension damping
- ▶ **CAUTION!** The prosthetic knee joint has to reach full extension before the heel strike, even at a slow walking speed. Observe the information under the heading "Fine-tuning the settings".
- ▶ Adjust extension damping so that the prosthetic knee joint does not swing too hard against the extension stop and so that the prosthetic foot swings through sufficiently according to the needs of the patient.

##### **Fine-tuning the settings**

- ▶ **CAUTION!** If the settings cannot be fully adjusted to the comfort needs of the patient, the settings must be established primarily based on safety aspects. The extension position has to be reached with every step and at any walking speed.
- ▶ Adapt the settings of the prosthetic knee joint to the patient by means of fine-tuning and exercises.
- ▶ Check the adjustment of the prosthesis settings during normal consultations and every 3 months during the safety inspection.  
Advise the patient and guardian to have the prosthesis inspected if changes in function occur.

## 7.5 Attaching the cosmetic cover

### CAUTION

#### **Use of highly hygroscopic particles (grease-absorbent substances such as talcum)**

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not allow the product to come into contact with highly hygroscopic particles.
- ▶ Note the following points for attaching the cosmetic cover:
  - Adapt the cosmetic cover to the patient and mount it on the prosthesis.
  - Check the settings from the dynamic trial fitting.
- ▶ To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.

## 7.6 Finishing the prosthesis

- ▶ **CAUTION!** To avoid a risk of falling, replace set screws that are too short with longer ones.
- ▶ Replace set screws that are too long and can cause product damage (for instance through collision or friction) with shorter ones.
- ▶ Before the patient leaves the fitting area, finish the prosthesis by tightening all screw connections. In doing so, note the instructions for use for all prosthetic components regarding installation torque values and thread lock.

## 8 Use

### 8.1 Information for use

#### NOTICE

#### **Mechanical overload**

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### 8.2 Cleaning

#### CAUTION

#### **Use of unsuitable cleaning agents or disinfectants**

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
- ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.

#### **Cleaning light soiling**

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

#### **Cleaning heavier soiling**

- > **Required materials:** cleaning cloth, 634A58 isopropyl alcohol, soft cloth
- ▶ **NOTICE!** Verify material compatibility. Do not remove lubricant from bearing points.
- ▶ Clean the product with a cleaning cloth and isopropyl alcohol.
- ▶ Dry the product with a cloth.

## Cleaning with Disinfectants

- > **Required materials:** colourless, alcohol-free disinfectant (verify material compatibility!), soft cloth
- ▶ Disinfect the product with the disinfectant.
- ▶ Dry the product with a cloth.
- ▶ Allow to air dry in order to remove residual moisture.

## 9 Maintenance

### CAUTION


#### **Failure to follow the maintenance instructions**

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ **NOTICE!** Do not lubricate and grease the prosthetic joint.
- ▶ **NOTICE!** Repair work must be performed exclusively by the manufacturer's technical service.
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Perform a safety inspection every 3 months.
- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.
- ▶ In case of changes in prosthesis function or relevant patient changes (such as the activity level, body weight, body height), perform unscheduled inspections and implement measures according to the situation.

#### **Inspecting and replacing the stop on page 11 –**

- ▶ Check the stop shown in Figure  on page 11 for wear and integrity every time the prosthesis is inspected.
- ▶ Replace a damaged or missing stop in the illustrated sequence.

## 10 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

## 11 Legal information

### 11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 11.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH