
## 50K324-US Agilium Reactive

<b>EN</b> Instructions for use .....	4
<b>FR</b> Instructions d'utilisation .....	9

6

**Material**

Al, PP, PA, PU, PES, EL, PE Schaum/PE foam, Stahl/steel

# 1 Foreword

## INFORMATION

Date of last update: 2020-05-15

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K324 Agilium Reactive knee joint brace.

## 2 Product description

Fig. 1/ Item:	Designation	Fig. 1/ Item:	Designation
1	Brace	2	Allen wrench
3	Screwdriver	4	Joint stop for flexion limitation
5	Joint stop for extension limitation	6	Strap guide loop
7	Condyle pad		

## 3 Intended use

### 3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 3.2 Indications

- Unicompartmental knee osteoarthritis
- Post-operative fitting after meniscus reconstruction
- Ligament injuries requiring unilateral relief

The indication must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

None known.

#### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent, swollen scars; reddening and hyperthermia of the treated limb; pronounced varicose veins, especially with impaired return flow; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area where the device will be applied; sensory and circulatory disorders in the legs, e. g. associated with diabetic neuropathy.

- Bicompartamental osteoarthritis of the knee

#### 3.4 Mechanism of Action

The brace relieves the medial or lateral compartment of the knee joint using the three-point principle. In the process, corrective force is applied to the contralateral condyle.

## 4 Safety

### 4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

 **CAUTION**

#### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

 **CAUTION**

#### **Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**⚠ CAUTION**

**Breakage of the brace frame or joint due to improper shaping**

Injuries caused by sharp edges, loss of function

- ▶ Avoid adapting the brace repeatedly.
- ▶ Do not make any improper changes to the brace.
- ▶ **Inform your patient.**

**⚠ CAUTION**

**Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

**NOTICE**

**Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

**NOTICE**

**Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

**NOTICE**

**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 5 Handling

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

### 5.1 Application

- > The patient should be seated.
- > The knee joint should be slightly bent and relieved of pressure.
- 1) Slide the brace over the leg from the front. When doing this, the brace joint must be on the lateral side of the affected compartment.
- 2) Position the brace on the leg in such a way that the brace joint is at the height of the middle of the patella (see fig. 2).
- 3) Fasten the top strap for the lower leg. When doing so, be sure that the tibia pad is under the fastening strap.
- 4) Close the rest of the straps in the following order: the bottom lower leg strap, the bottom upper leg strap, the top upper leg strap. Do this by pulling the straps through the strap guide loops and fastening them.
- 5) Ensure that the brace is in a firm and comfortable position.

### 5.2 Adaptation

#### Adjusting the corrective force

- > **Prerequisite:** The patient should be seated.
- > **Prerequisite:** The knee joint should be relieved of pressure and slightly bent.
- 1) Set the corrective force on the socket screw above the brace joint (see fig. 3). This changes the varus/valgus angle of the brace (see fig. 6).
- 2) Check the corrective force while standing and walking.
- 3) If necessary, adapt the corrective force until the desired adjustment has been made.

### **Optional: shorten the bottom lower leg and upper leg straps**

- 1) Remove the Y-hook-and-loop.
- 2) Shorten the straps to the required length.
- 3) Stick the Y-hook-and-loop onto the strap.

### **Optional: shorten the top lower leg strap**

- 1) Loosen the hook-and-loop closure on the strap guide loop.
- 2) Remove the hook-and-loop circle.
- 3) Shorten the strap to the desired length.
- 4) Fasten the hook-and-loop circle to the top lower leg strap.

### **Optional: limiting the knee flexion and extension**

The extension limitation at delivery is **0°**. The flexion is not limited.

- 1) Remove the condyle pad.
- 2) To limit the extension, remove the anterior screw. To limit the flexion, remove the posterior screw.
- 3) Choose the joint stops appropriate for the desired flexion and extension.
- 4) Insert the joint stops into the brace knee joint. Insert the joint stop for extension limitation at the front and the joint stop for flexion limitation at the back. Ensure that the heads of the joint stops are pointing outwards (see fig. 4).
- 5) Secure the joint stop with the screw (see fig. 5).
- 6) Re-attach the condyle pad.

## **5.3 Removal**

- > **Prerequisite:** The patient should be seated.
  - > **Prerequisite:** The knee joint should be relieved of pressure and slightly bent.
- 1) Open and unfasten all straps.
  - 2) Pull the brace forwards off the leg.

## **5.4 Cleaning**

### **NOTICE**

#### **Use of improper cleaning agents**

Damage to the product due to use of improper cleaning agents

- Only clean the product with the approved cleaning agents.

- 1) Remove the straps and pads from the brace.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.

- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 5) If the brace comes into contact with salt water or dirt, rinse with clear water and allow to air dry.

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

### 7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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