17B203=* E-MAG Active

DE Gebrauchsanweisung (Benutzer)	3
EN Instructions for use (user)	8
FR Instructions d'utilisation (Utilisateur)	14
IT Istruzioni per l'uso (Utilizzatore)	20
ES Instrucciones de uso (Usuario)	25
PT Manual de utilização (Usuário)	31
NL Gebruiksaanwijzing (Gebruiker)	37
SV Bruksanvisning (Användare)	43
DA Brugsanvisning (Bruger)	48
NO Bruksanvisning (Bruker)	54
FI Käyttöohje (Käyttäjä)	59
CS Návod k použití (Uživatel)	65
SK Návod na používanie (Užívateľ)	70
RU Руководство по применению (Пользователь)	76
ZH 使用说明书 (用户)	82

1 Vorwort

Deutsch

INFORMATION

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INFORMATION

Last update: 2016-01-29

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Have yourself instructed by qualified personnel in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

Your orthopaedic technician has fabricated an individual leg orthosis for you. The integrated E-MAG Active is a knee joint system with stance phase control, which offers you a free swing phase and locks the knee joint for a secure stance phase prior to heel strike. The knee joint determines the position of your leg with the help of the electronics and enables the swing phase accordingly. Read the "Quick Start", 646D517.

The E-MAG Active is a sealed system, which means that the joint is less susceptible to contamination and also prevents clothing from becoming trapped. Not only is this gentle on your clothing but, above all, offers you enhanced safety.

Two clearly distinguishable switches on a button control unit allow you to electronically unlock the joint (once) or engage the lock (permanently). In addition, the mechanical (temporary) unlock feature supports a larger area of application, for example when cycling. These instructions for use describe the functions and operation of the E-MAG Active as well as how to charge it and care for it.

Please do not hesitate to ask your orthopaedic technician if you have any technical questions.

2 Safety

2.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

2.2 General safety instructions

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

CAUTION

Explosion hazard

Do not throw the battery into fire. Never expose it to heat. Avoid direct sunlight.

CAUTION

Risk of injury in case of continued use after malfunctions

Never use the E-MAG Active knee joint if there is a malfunction. Please contact your specialist workshop immediately.

CAUTION

Risk of injury due to electromagnetic interference

Malfunctions may occur near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). This can cause the patient to fall.

Avoid proximity to strong sources of magnetic and electric interference (e.g. transformer stations, transmitters).

CAUTION

Risk of falling due to improper operation of the service switches

The small service switches A and T must not be activated or manipulated. They are exclusively for use by qualified personnel in your specialist workshop.

CAUTION

Risk of falling while using the control unit

To avoid the risk of falling, you may only activate the locking and sitting function while standing and not while walking.

CAUTION

Risk of falling due to use of grease (or other lubricants)

This could contaminate the joint and jam the lock. The bearing washers and joint locking mechanism should never be lubricated.

NOTICE

Risk of damage due to water/moisture

Water will destroy the E-MAG Active. Do not wear your orthosis while taking a shower, bathing, or swimming, and do not stay in environments with high humidity with your orthosis.

NOTICE

Risk of damage due to incorrect environmental conditions

Solid particles and fluid in the joint, lock or electronics of the E-MAG Active may lead to malfunctioning. Do not expose the components to excessive smoke or dust in the environment, mechanical vibrations or shocks.

⚠ CAUTION

Exposure of the product to unsuitable environmental conditions

Patient injury, damage, brittleness or destruction due to improper handling

- ▶ Do not expose the product to condensing ambient humidity or liquids.
- ▶ Do not expose the product to abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

3 Initial operation

INFORMATION

The E-MAG Active can only be operated with the intended 317B20 battery. It was developed especially for the E-MAG Active. Both the receptacle and lock unit are designed to fit with the form of this battery.

INFORMATION

Before putting the E-MAG Active into operation, read the instructions for the 317B20 battery and the 317L20 battery charger. In particular, please be sure to take note of the safety information provided.

This section will describe how the E-MAG Active can be activated/deactivated and how it can be reactivated from the energy-saving mode. The section also explains how the electronic system starts and which signals the system emits during activation.

3.1 Charging the battery

INFORMATION

After charging, the E-MAG Active can be used for around 5,000 steps (corresponds to about 5 km) with one battery charge. Nevertheless, charging the battery daily (ideally at night) is recommended.

INFORMATION

In daily use, the charging time depends on the remaining capacity and maintenance of the battery. The charger can charge two batteries at the same time.

INFORMATION

The battery should always be removed when the orthosis is not in use.

INFORMATION

The LEDs on the charger must not be illuminated prior to inserting the battery. If they are illuminated, please submit the charger to an authorised Ottobock Service Centre.

First plug the battery charger into a power outlet and then insert the battery into the charger (see fig. 1). One full charging cycle takes 8 to 10 hours. Place the charged battery into its receptacle; it should snap easily into the holder. The contacts have to face up and in. In order to do so, insert the battery into the receptacle at an angle and press it against the orthosis. One can hear the battery lock into place when it is properly inserted.

3.2 Activating/deactivating the E-MAG Active

CAUTION

Risk of falling

Missing signals during activation indicate a failure of the E-MAG Active. Remove the battery immediately and then reinsert. Please contact your specialist workshop immediately if the error occurs again.

INFORMATION

Failure to unlock the knee joint for more than one hour automatically causes the control electronics to switch to stand-by mode. To reactivate the system, simply push the top button on the control electronics.

Activation is performed by inserting the battery. Every time the E-MAG Active is activated, the electronic control unit runs a short self-test that emits the following signals:

Signal	Visual feedback (LEDs)	Audible feedback	Tone / length
System ready > Knee joint locked	Short flashing of all LEDs (white)		Rising
System not ready > Knee joint flexed (open)	Continuous flashing of both LEDs (white and green)		Tone-pause-tone series

When the E-MAG Active is flexed, the tone-pause-tone series is emitted as a warning signal until the E-MAG Active is extended/locked. Deactivation is performed by removing the battery.

4 Operation

This section explains how to operate the knee joint control unit by button and tells you about the warnings used by the feedback system of the E-MAG Active to inform you of malfunctions.


4.1 Control unit

CAUTION

Risk of falling due to unlocking / locking

Both permanent and one-time unlocking as well as permanent locking of the system can result in an elevated risk of falling. The patient should not use these functions whilst walking.

The electronic control unit measures the position of the orthosis while walking. It unlocks the joint before toe-off and locks it prior to heel strike at the end of the swing phase, so that the patient can step down safely. Your orthopaedic technician adapted the electronic control unit of the E-MAG Active to your needs during the first trial fitting. Two manual switching possibilities are available to you to give you full control over the E-MAG Active (see fig. 2):


4.1.1 Electronic locking (permanent)

If you feel unsafe, you can lock the joint at any time using the top button for locking the knee joint. This may be the case when walking on a slope or if you want to stand up for an extended period of time without focusing on the electronic control unit. In order to do so, briefly press the upper button twice in quick succession (like double-clicking a computer mouse). This function is confirmed by an audible signal. It is deactivated with a double mouse click on the upper button until you hear an audible tone.

4.1.2 Electronic unlocking (once)

With the lower button, the patient can unlock the joint once (e.g. in order to sit down). In order to do so, briefly press the lower button twice in quick succession (like double-clicking a computer mouse). This function is not confirmed by an audible signal.

4.1.3 Mechanical unlocking (temporary)

For certain activities (e.g. cycling), it may be beneficial to temporarily unlock the E-MAG Active joint. Mechanical unlocking of the system is indicated by the symbol of a cyclist (switch to "🚲").

In order to restore the functionality of the E-MAG Active joint, cancelling the unlocking feature is required. In order to secure the system again, the switch is operated again in the opposite direction. The normal state of the joint is indicated by the symbol of a walking person (switch to "🚶").

4.2 Warning signals

4.2.1 Battery charge level

INFORMATION

The electronics will turn off if the battery is discharged completely. While the lock can still be opened manually, it is a slightly cumbersome task.

Should the battery charge level fall below the pre-adjusted capacity value, the electronics will emit the warning signals listed here:

Signal	Visual feedback (LEDs)	Acoustic feedback	Tone / length
Battery low	Alternate flashing (white and red)		High-low series / individual signals
Battery empty			Falling / long

The warning signal is repeated every 15 minutes. Replacing the battery immediately is highly recommended. The number may vary. How quickly the battery is discharged depends on the ambient temperature: The colder it is, the faster the battery will be drained. The electronics will turn off if the battery power is insufficient. If this occurs, the knee joint can only be unlocked manually (see Page 12).

4.2.2 Failure of the electronics

The E-MAG Active is not vulnerable to failure. However, should disruptions occur, the electronics recognise the failure and send a warning signal.

Signal	Visual feedback (LEDs)	Audible feedback	Tone / length
Electronics malfunction	Continuously lit (all LEDs except green)		High/continuous tone

In the case of such a failure, the E-MAG Active knee joint can only be opened manually (see Page 12).

5 Maintenance and care

5.1 Maintenance

Have your E-MAG Active checked for wear and for proper function of the electronic components by your orthopaedic technician every six months.

5.2 Care

- 1) Keep your orthosis away from environments with high humidity. Use a slightly moistened cloth to clean the orthosis, E-MAG Active or any of the components from outside.
- 2) Since the E-MAG Control is a closed locking system, clothing fibres cannot get into the mechanism. However, please observe the maintenance intervals of six months, because particulate matter can settle in the electronics and mechanism and can disturb the E-MAG Active.

6 Technical data & accessories

Storage temperature range	-20 °C to +70 °C (-4 °F to 176 °F)
Operating temperature range	-15 °C to +50 °C (-5 °F to 122 °F)
Relative humidity for both ranges	0 % to 100 %
Power supply for joint	NiMh battery, 4.8 V nominal voltage
Range of a battery pack	Approx. 5,000 steps, corresponds to about 5 km

NOTICE

Only use accessories that have been approved for the E-MAG Active system.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

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Otto Bock HealthCare GmbH