8E33=*, 8E34=*

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	14
FR Instructions d'utilisation (Personnel spécialisé)	25
IT Istruzioni per l'uso (Personale tecnico specializzato)	36
ES Instrucciones de uso (Personal técnico especializado)	47
PT Manual de utilização (Pessoal técnico)	58
NL Gebruiksaanwijzing (Vakmensen)	69
SV Bruksanvisning (Fackpersonal)	80
JA 取扱説明書 (有資格担当者)	91

INFORMATION

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INFORMATION

Last update: 2016-11-03

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The 8E33=7-1, 8E34=7-1, 8E33=9-1, 8E34=9-1 System Electric Greifers are referred to as the product/prosthesis/terminal device in the following.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Design

The product consists of the following components:

1. **Fingertips**
The position can be rotated by $\pm 60^\circ$ to adapt to various tasks.
For other types of fingertips, see the section "Accessories" (see Page 18)
2. **Safety lever**
To unlock the prehensile arm drive (see Page 22)
3. **Dial**
For manual opening/closing
4. **On/Off switch**
To turn the electronics and flexion joint on/off
Position of the rocker switch towards the terminal device: electronics switched on and flexion brake released
Position of the rocker switch towards the prosthetic socket: electronics switched off and flexion joint blocked
5. **Quick-disconnect wrist unit** (only with 8E33*)
Mechanical connection to the lamination ring

2.2 Function

The product is a myoelectrically controlled terminal device of the MyoBock prosthesis system.

Various control programmes make optimal adaptation to the individual needs and abilities possible.

Only with 8E33=9* and 8E34=9*

The terminal device is based on the Ottobock DMC system (DMC=Dynamic Mode Control): This system, which has been developed by Ottobock, uses two independent measurement and control systems in order to optimally control the gripping speed and gripping force in accordance with the muscle signal.

2.3 Combination possibilities

This product can be combined with the following Ottobock components:

Power supply (rechargeable battery)

- 757B35=0 MyoEnergy Integral
- 757B35=1/757B35=3 MyoEnergy Integral
- 757B15 X-ChangePack
- 757B20/757B21 EnergyPack

Battery chargers

The following battery chargers can be used depending on the rechargeable battery:

- 757L14 battery charger
- 757L20 battery charger
- 757L35 battery charger

Lamination rings

- 10S1=50 lamination ring
- 9S110=50 lamination ring

Rotation

- 10S17 electric wrist rotator
- 13E205 MyoRotronic

Elbow

- ErgoArm Hybrid plus: 12K44
- ErgoArm Electronic plus: 12K50
- DynamicArm: 12K100N/12K110N

3 Application

3.1 Indications for use

The System Electric Greifers are intended **exclusively** for exoprosthetic upper limb fittings.

3.2 Use/Field of Application

The System Electric Greifers can be used for unilateral or bilateral amputations, and in case of dysmelia.

3.3 Environmental conditions

The prosthesis system was developed for everyday activities and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the prosthesis system should not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see Page 23).

3.4 Qualification

Fitting a patient with the product may only be carried out by orthopaedic technicians who have been authorised by Ottobock after completion of a corresponding training course.

4 Safety

4.1 Explanation of warning symbols

WARNING

Warning regarding possible serious risks of accident or injury.

CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

4.2 Structure of the safety instructions

⚠ CAUTION**The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

⚠ WARNING**Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ CAUTION**Independent manipulation of the product**

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION**Signs of wear on the product**

Injury due to faulty control or malfunction of the product.

- ▶ In the interest of the patient's safety and in order to maintain operating reliability, the product should be serviced at regular intervals.

4.4 Information on Alignment/Adjustment

⚠ CAUTION**Incorrect electrode settings**

Injury due to unexpected product behaviour.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.

- ▶ Set the electrode gain as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

⚠ CAUTION

Use of unsuitable prosthetic components

Injury due to unexpected product behaviour.

- ▶ Use the product only in combination with components listed in the section "Combination Possibilities" (see Page 15).

4.5 Information on Use

⚠ CAUTION

Risk of pinching between the prehensile arms/fingertips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the prehensile arms/fingertips when using the product.

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Operation outside of the permissible temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid operation in areas outside the permissible temperature range (see Page 23).

⚠ CAUTION

Penetration of dirt and humidity into the product

Injury due to unexpected behaviour of the prosthesis system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Use only the following products for cleaning/disinfecting the inner socket:

Cleaning: 453H10=1 Ottobock DermaClean

Disinfection: 453H16 Ottobock DermaDesinfect.

4.6 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the product with soiled or damaged contacts

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

NOTICE

Contact of the charging plug with magnetic data carriers

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

INFORMATION

The prosthesis system cannot be used as long as the charging plug is connected to the charging receptacle. The prosthesis system is deactivated for the duration of charging.

4.7 Notice for using the torch

⚠ WARNING

Use of the incorrect battery type and/or failure to note the polarity

Explosion due to escaping gases and their self-ignition caused by short circuit and excessive temperature.

- ▶ Only use type CR2032 batteries.
- ▶ When inserting the two batteries, note the polarity information on the cover of the battery compartment.

⚠ CAUTION

Looking directly into the beam of light from the torch

Eye damage due to looking directly into the beam of light.

- ▶ Do not use the torch for examinations, for instance of the eyes (pupil test).
- ▶ Do not aim the beam of light directly at the eyes of people or animals.
- ▶ Do not look directly into the beam of light.
- ▶ If the beam of light is aimed at the eyes, be sure to close your eyes and immediately move the head out of the beam of light.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 2 pc. CR2032 lithium battery

- 1 pc. 8E33=9-1 DMC VariPlus System Electric Greifer or
- 1 pc. 8E34=9-1 DMC VariPlus System Electric Greifer or
- 1 pc. 8E33=7-1 Digital Twin System Electric Greifer or
- 1 pc. 8E34=7-1 Digital Twin System Electric Greifer
- 1 pc. 647G1274 instructions for use (qualified personnel)
- 1 pc. 647G1275 instructions for use (user)

5.2 Accessories

Fingertips

- 9S138 tip set
- 9S145 rubber gripping pad set
- 9S149 fingertip blank set
- 9S234=PAA tip pair
- 9S278=PAA tip pair


6 Preparation for use

6.1 Control programmes for 8E33=7-1, 8E34=7-1

The Digital Twin control unit features a digital and dual-channel control mode. The desired control mode is selected with the integrated 13E189 function plug.

6.1.1 Digital mode programme (with function plug)


Control with 2 electrodes.

The terminal device **opens and closes at constant speed** after the electrode signal exceeds a threshold

6.1.2 Dual-channel mode programme (without function plug)

Control with one electrode.

Closing:

Slow, weak muscle signal

Opening:

Fast, strong muscle signal

6.2 Control programmes for 8E33=9-1, 8E34=9-1

The "757T13 MyoSelect" programming device is used to select the desired control programmes.

The following control programmes are recommended:

- **Control programme 1 DMC plus:** for patients with two strong electrode signals
- **Control programme 3 VarioControl:** for patients with one strong electrode signal

6.2.1 Programme 1: DMC plus

Proportional control with 2 electrodes for patients with two strong electrode signals
Depending on the strength of the electrode signal, the opening and closing speed can be varied proportionally between 8 mm/s and 200 mm/s. The gripping force can also be varied proportionally between 10 N and 160 N.

Control	2 electrodes, proportional
Closing	Proportional to the signal strength on the "CLOSE" electrode
Opening	Proportional to the signal strength on the "OPEN" electrode
Gripping force	Proportional to the signal strength on the "CLOSE" electrode

6.2.2 Programme 2: AutoControl - LowInput

Proportional control with 2 electrodes or a switch.
The speed is increased proportionally depending on how fast the opening signal on the electrode rises. When a switch is used, opening is at a constant speed. Closing is at the maximum speed of 200 mm/s. The gripping force increases proportionally to the time. The longer the signal on the "CLOSE" electrode is applied, the higher the gripping force.

Control	2 electrodes or switch, proportional
Closing	Constant speed with application of the signal on the "CLOSE" electrode or activation of the switch
Opening (electrode)	Proportional to the signal strength on the "OPEN" electrode
Opening (switch)	Constant speed with activation of the switch
Gripping force	Proportional to the signal duration on the "CLOSE" electrode
Indication	For patients with two weak electrode signals

6.2.3 Programme 3: VarioControl

Proportional control with 1 electrode or 1 linear control element.
The terminal device is opened or closed depending on how fast and with what intensity the signal rises or falls. By dropping the signal very slowly (very slow muscle relaxation), the terminal device stops in the open position.

- **Control with one electrode:** for patients with one strong electrode signal or co-contraction
- **Control with one linear control element:** for patients with an electrode signal that is too weak, or no electrode signal

Control	1 electrode or linear control element, proportional
Closing	Proportional to the increase in the strength and speed of the signal
Opening	Proportional to the decrease in the strength and speed of the signal
Gripping force	Proportional to the decrease in the strength and speed of the signal
Stopping	Through very slow drop in the signal

6.2.4 Programme 4: VarioDual

Proportional control with 2 electrodes for patients with two strong electrode signals.
The terminal device is opened or closed depending on how fast and with what intensity the signal on the first electrode rises or falls. The gripping force is determined by the subsequent or simultaneous electrode signal on the second electrode.

Control	2 electrodes, proportional
Closing	Proportional to the decrease in the strength and speed of the signal on electrode 1.
Opening	Proportional to the increase in the strength and speed of the signal on electrode 1.
Gripping force	Proportional to the signal strength on electrode 2.

6.2.5 Programme 5: Digital Control

Proportional control with 2 electrodes or a switch.

The terminal device opens and closes at the fastest speed. The gripping force increases proportionally to the duration of gripping.

Control	2 electrodes or switch, digital
Closing	Constant speed with detection of the "OPEN" signal on electrode 1 or the switch
Opening	Constant speed with detection of the "CLOSE" signal on electrode 2 or the switch
Gripping force	Proportional to the signal duration on the "CLOSE" electrode or switch. The longer the signal is applied, the higher the gripping force.
Speed	Constant speed of 200 mm/s for opening and closing

6.2.6 Programme 6: Double Channel Control

Control with 1 electrode or 1 linear control element.


The terminal device closes with a slow, weak signal and opens with a fast, strong signal. When one electrode is used, the upper switching threshold can be adjusted with the 757T13 MyoSelect.

Control	1 electrode
Closing	Slow, weak electrode signal that does not reach the upper threshold within 80 ms.
Opening	Quick, strong electrode signal that reaches the upper threshold within 80 ms.
Gripping force	Proportional to the signal duration. The longer the signal is applied, the higher the gripping force
Speed	Constant speed of 200 mm/s for opening and closing

7 Use

7.1 Torch


(for terminal devices with built-in torch)

1. **2x powerful LEDs**
to light the working environment
2. **On/Off switch**
Pressing this button turns the torch on or off. If this button is pressed immediately after switching on, the brightness can be changed in 3 levels by pressing it again. After about 15 minutes, the torch is switched off automatically.
3. **Battery compartment**

Inserting the batteries (2x CR2032)

To insert the batteries, complete the following steps:

- 1) Turn the cover on the battery compartment counter-clockwise and remove it.
- 2) Insert 2 CR2032 button cells with the + pole facing the cover.
- 3) Install the cover. In doing so, make sure the rubber seal is not damaged.
- 4) Close the cover by turning it clockwise.

7.2 Opening the terminal device in an emergency

This safety function allows the prehensile arms to be opened regardless of the control signals.

- Swing the unlock lever away from the dial (arrow) and open the prehensile arms manually.

7.3 Changing the 8E33=7*, 8E33=9* terminal device

Separating the terminal device from the socket

- 1) Rotate the terminal device once around its own axis until a slight resistance is felt (about 360°).
- 2) Overcome this resistance and pull the terminal device off the socket.

Attaching the terminal device to the socket

- 1) Insert the quick-disconnect wrist unit into the lamination ring and press it in firmly.
- 2) Turn the terminal device slightly to the left or right.
- 3) Verify proper attachment by pulling on the terminal device.

8 Cleaning

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.
Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

9 Maintenance

Since all moving mechanical parts are subject to wear, regular service inspections are recommended. These service inspections include an inspection of the product and, when necessary, the replacement of wear parts.

10 Legal information

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

11 Technical data

Environmental conditions	
Storage and transport in original packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing
Storage and transport without packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing
Operation	-10 °C/+14 °F to +60 °C/+140 °F max. 93 % relative humidity, non-condensing

Reference number	8E33=7-1	8E34=7-1	8E33=9-1	8E34=9-1
Operating voltage [V]	6/7.2			
Opening width [mm]	95			
Gripping force [N]	160		0 - 160	
Weight [g]	540	520	540	520
Service life [years]	5			
Fingertip torque value [Ncm]	40 - 50			
Proportional speed [mm/sec]	180		8 - 200	

Torch	
Illuminants [number, type]	2, LED
Total maximum luminosity of the two LEDs [mcd]	16400-24000
Battery [number, type]	2x CR2032
Maximum operating time with new batteries at highest brightness and at room temperature [h]	8
Brightness levels [number]	3
Automatic switch-off	after 15 minutes

12 Appendices

12.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

SN YYYY WW NNN

Serial number
 YYYY – year of manufacture
 WW – week of manufacture

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Otto Bock Healthcare Products GmbH