

MyoSelect 757T13

DE Gebrauchsanweisung	5
EN Instructions for Use	12
FR Instructions d'utilisation	19
IT Istruzioni per l'uso	26
ES Instrucciones de uso	34
PT Manual de utilização	41
NL Gebruiksaanwijzing	48
SV Bruksanvisning	55
PL Instrukcja użytkowania	62
CS Návod k použití	70
TR Kullanma talimatı	77
EL Οδηγίες χρήσης	84


7

8

9

3

10

4

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English

Date of the last update: 2014-02-26

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

1 Application

The 757T13 MyoSelect is to be applied for the identification and adjustment of Ottobock Myobock components such as System Electric Hands, System Electric Greifers, 4-Channel Processor II, ErgoArm Electronic plus and BionicArm. Do not use the 757T13 MyoSelect for testing and measuring purposes.

2 Safety Information

- Only certified persons may identify and adjust Ottobock Myobock Components with the 757T13 MyoSelect.
- Modifications to components are explicitly forbidden.
- Do not let foreign particles or fluids enter the MyoSelect.
- Do not subject the MyoSelect to intense smoke or dust, or to extreme mechanical vibrations, shocks or great heat.
- Avoid the vicinity of strong magnetic or electric fields such as transformer stations or transmitters.
- You may only open the 757T13 MyoSelect if the batteries need to be changed.
- Repairs of the 757T13 MyoSelect may only be done by an authorised Ottobock Service facility.

- Only take the 757T13 MyoSelect into operation within the confines of an orthopaedic workshop.
- Before using the 757T13 MyoSelect, check for any visible damages.
- If one or several components are damaged or malfunction, please contact Otto Bock Healthcare GmbH.

3 Scope of Delivery

- 1 pce. MyoSelect
- 1 pce. Instructions for Use
- 1 pce. Instructions - Control Types
- 1 pce. Coding Plug black

4 Technical Characteristics

Attention: In order to be able to make program modifications to the 8E38=8 / 8E39=8 SensorHand Speed, 8E38=5/8E39=5/8E41=5 SensorHand, 13E195 4-Channel Processor II and 12K50 ErgoArm Electronic plus with the 757T13 MyoSelect, all these components must be equipped with a black Coding Plug. Please read the Manual accompanying the individual Myo Component carefully.

4.1 Multi-Function Button (Fig. 1, Fig. 3)

The Multi-Function Button (Fig. 1, Pos. 3) controls the functions of the MyoSelect. It can be operated in three different ways:

- Pressing inwards (Fig. 3, Pos. 8)
- Pressing downwards (Fig. 3, Pos. 7)
- Pressing upwards (Fig. 3, Pos. 9)

4.2 LCD Display (Fig. 1, Pos. 2)

Shows information relating to the connected Myobock Component as well as to the MyoSelect itself (connection display, battery voltage alarm, etc.).

4.3 Connections (Fig. 1, Fig. 2)

Communication between prosthesis and MyoSelect is achieved by connecting a:

- coaxial plug (Fig. 1, Pos. 1) e.g. System Electric Hand and Electric Greifer with Quick Disconnect Wrist.
- coaxial socket (Fig. 1, Pos. 4) e.g. Four-Channel Processor II, ErgoArm Electronic plus, BionicArm.
- 4-pole socket (Fig. 2, Pos. 5) e.g. System Electric Hand and Electric Greifer without Quick Disconnect Wrist

4.4 Acoustic Signaller

The Acoustic Signaller gives additional acoustic feedback after the successful saving of the adjustments.

5 Description

The 757T13 MyoSelect replaces mechanical coding plugs and potentiometers through electronics. External adjustment is easier and quicker done without having to disassemble the Myobock Component.


The MyoSelect displays following functions for Myobock Component:

- Product number of the Myobock Component
- Software Version of the Myobock Component

- Serial number of the Myobock Component from a given production date on
- Read-off and adjustment of Myobock Component control type
- Read-off and adjustment of Myobock Component adjustment bar
- Whether adjustments to the Myobock Component are "customized"
- Automatic reset to NiCd/NiMH-Battery Technology
- Signals low battery charge for the MyoSelect

For further information concerning control type and function of the adjustment bar of the individual Myobock Component please refer to the Instruction Manual of the individual Myobock Component.

5.1 Display

8Exx=8

Product Number

Is used for the identification of the Myobock Component. Non-identifiable elements are characterised by a small "x".

U2.2

Software Version

Is used for the identification of Software.

S/N= ...

Serial Number (optional)

If available, the number given by the manufacturer is displayed.

02

Program Number of Control Type

The big-lettered 2-digit number shows the current active Control Type. Control type is determined by Function or Coding Plugs (display in grey) or electronically selected via MyoSelect (display blinking in black).

! 01

Customizing Symbol

If the Myobock Component has been "**Customized**" to the individual patient a warning triangle will be displayed in addition to the program number. No adjustments can be made (grey display).

Switch (optional)

Gives information whether the switching function (if available - see Addendum) specified by the Myobock component is switched ON or OFF.

Symbol for:

Switch "ON"

Switch "OFF"

Adjustment Bar (optional)

Displays the current adjustment of the Adjustment Bar. These are determined by potentiometer adjustments (display in grey) or electronically selected by MyoSelect.


Actual state of the Myobock component is displayed with stable grey display of **Program Number** or **Adjustment Bar**. No adjustments can be made.

5.2 MyoSelect Specific Displays

Batteries low

If the Battery symbol on the left side of the display starts to blink then the battery capacity is low. Functionality is not yet affected. Exchange batteries promptly as described in Chapter 7.

Batteries empty

A big blinking battery symbol will appear on your display when the batteries are empty. You cannot continue working with the MyoSelect. Change batteries immediately, as described in Chapter 7. "V1.0" on the screen refers to the version of the MyoSelect.

5.3 Coding / Function Plug Allocation

 Coding Plug	Program
13E184=1 white	1
13E184=2 red	2
13E184=3 green	3
13E184=4 blue	4
13E184=5 yellow	5
13E184=6 purple	6
13E184=7 orange	7

 Function Plug	Program
plugged in	1
unplugged	2

6 Operating MyoSelect

6.1 Connecting the Myobock Components (Fig. 5, Fig. 6)

Connect Myobock component to the corresponding plug as described in Chapter 4.3. Hold the MyoSelect in one hand and support the Myobock Component with your other hand (Fig. 5, Fig. 6)


Only **one** Myobock component at a time may be connected to the MyoSelect


6.2 Establishing Connection (Fig. 3)

Activate MyoSelect by a short push on the Multi-Function Button (Fig. 3, Pos. 8).

While MyoSelect tries to build-up a connection, the following will appear on the display:

If no MyoComponent is connected or no connection can be established, the following symbol will show up on the LCD display:

Each new push on the Multi-Function Button will restart the connection build-up attempt.

After a successful connection build-up the Myobock component will automatically be reset to NiCd/NiMH Battery Technology.

6.3 Making Adjustments (Fig. 3)

If the Myobock component is adjustable, then the display begins to blink.

1. CHOOSE

A **short press** on the Multi-Function Button (Fig. 3, Pos. 8) will choose the adjustable parameter *Program Number, Adjustment Bar 1, Adjustment Bar 2* or *Switch* (if available and displayed in black). The current choice is displayed by the blinking parameter.

2. ADJUST

By a **short flick upward or downward** on the Multi-Function Button (Fig. 3, Pos. 7, Pos. 9) you can increase or decrease the adjustment value of the parameter by one step at a time.

Keep flicking and hold the button in position for a longer period (Fig. 3, Pos. 7, Pos. 9): this will result in a faster increase or decrease of the parameter values.

3. SAVE

Hold the Multi-Function Button (Fig. 1, Pos. 3) **until the blinking stops**. You will hear a beep. The values are saved into the Myobock component. After 2 seconds the program number will begin to blink again and MyoSelect is back in adjustment mode.

Successful saving is confirmed by a beep. If saving was unsuccessful then the symbol for "no connection" will appear on the display.

4. FINISH

After **unplugging** the Myobock component the symbol for "no connection" will re-appear.

If the Myobock component was disconnected before saving (Item 3) then the values originally adjusted are maintained in the Myobock component.

In order to conserve the batteries, the MyoSelect will automatically switch-off if the device is not operated over a certain period of time. You can re-activate the device by pressing the Multi-Function Button (Fig. 3, Pos. 8).

6.4 Examples


Ex. 1: 13E195 Four-Channel Processor II with attached 13E196 Adjustment Cap

If Program Number and Adjustment Bars are **displayed in grey** then the electronics are equipped with traditional potentiometers and coding plugs that **only allow adjustment control**.


Program Number "09" refers to the **Adjustment Cap** of the Four-Channel Processor II. Potentiometer "A" is in (approx.) position 4,3, potentiometer "B" in position 2,8.

Ex. 2: SensorHand Speed

A **black coding plug** has been inserted into the Hand thus allowing for **electronic adjustment** of control type (**display in black**). The display blinks. Program Number "**02**" corresponds to Program "red coding plug". Independently of chosen control type you can adjust the Adjustment Bar "A" (display in black).

Ex. 3: SensorHand Speed (adjustment changed by Customizing)

This example shows the patient-specific adjusted Program "**01**" done with "**Customizing**". No adjustments with MyoSelect are possible.

7 Battery Exchange (Fig. 2, Fig. 4)

1. For inserting batteries open the housing by unscrewing the two torx screws counter-clockwise with a T10 screw driver (Fig. 2, Pos. 6).
2. Lift-off housing cover (Fig. 1, Pos. 11).
3. Remove empty batteries, insert new batteries. Ensure correct polarity. (Fig. 4, Pos. 10)!
4. Return housing lid and fasten with torx screws.

Do not overwind torx screws.

8 Technical Data

Power Supply	2 Batteries Type Micro (AAA)
Operating Voltage 3V/DC	
Operating Temperature	-10 °C to +50°C (14°F - 122°F)

9 Disposal

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

10 Trademarks

All product names mentioned in this accompanying document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners. Should trademarks used in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11 Service

The 757T13 MyoSelect does not necessitate a service during its life cycle. If a malfunction is suspected please send the device to your Ottobock Service Facility.

12 Warranty

Warranty only applies if the 757T13 MyoSelect is used according to the specified conditions and for the intended purposes. Ottobock warranty applies only to tested component combinations.

13 Declaration of Conformity

This device meets the requirements of the European Directive 2004/108/EC. The declaration of conformity was created by the manufacturer according to Annex II of the Directive.

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Otto Bock Healthcare Products GmbH