
1C70 Evanto

DE Gebrauchsanweisung	6
EN Instructions for use	12
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	23
ES Instrucciones de uso	29
PT Manual de utilização	35
NL Gebruiksaanwijzing	41
SV Bruksanvisning	47
DA Brugsanvisning	52
NO Bruksanvisning	58
FI Käyttöohje	63
PL Instrukcja użytkowania	69
HU Használati utasítás	75
CS Návod k použití	81
RO Instrucțiuni de utilizare	86
HR Upute za uporabu	92
SL Navodila za uporabo	98
SK Návod na používanie	104
BG Инструкция за употреба	109
TR Kullanma talimatı	115
EL Οδηγίες χρήσης	121
RU Руководство по применению	127
JA 取扱説明書	133
ZH 使用说明书	139
KO 사용 설명서	144

1

Bench Alignment, transtibial

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INFORMATION

Date of last update: 2023-05-26

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 1C70 Evanto prosthetic foot has two exchangeable heel wedges. The heel wedges influence the behaviour of the prosthetic foot in different phases of the gait cycle. This allows the prosthetic foot to be individually adapted.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

Combination with prosthetic knee joints	
From stiffness 4 and foot size 28	Do not combine with 3R85 Dynion
From stiffness 5 and foot size 26	Do not combine with 3R55

The prosthetic foot generates high moments in the ankle area. Use structural components with higher weight limits:

Body weight [kg]	87 to 100	101 to 125
Foot size [cm]	From 24	From 26
Structural component weight limit [kg]	≥125	≥150

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

TIP: Also try the next higher stiffness so the user can experience the difference.

Body weight [kg]	Spring stiffness
Up to 58	1

Body weight [kg]	Spring stiffness
59 to 72	2
73 to 86	3
87 to 100	4
101 to 125	5

2.3 Environmental conditions

Storage and transport
Temperature range -20 °C to +60 °C, relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions
Temperature range: -10 °C to +45 °C
Chemicals/liquids: fresh water, soap solution, chlorine water, salt water ≤ 3.5 % salt content
Moisture: Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions
Solids: Dust, occasional contact with sand
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see page 16).

Prohibited environmental conditions
Solids: highly liquid-binding particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand
Chemicals/liquids: perspiration, urine, acids, continuous use in liquid media

2.4 Lifetime

Depending on the patient's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The products are wear and tear parts that are subject to normal degeneration.

3 General safety instructions

CAUTION! Risk of injury and risk of product damage

- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 12).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

NOTICE! Danger of product damage and restricted functionality

- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Quick start guide	647G1809=*
1	Prosthetic foot	–
1	Protective sock	Size 22 to 25: 2U3=1-7 Size 26 to 30: 2U3=2-7
Additional heel wedges		
1	Upper heel wedge	2F70=*
1	Lower heel wedge	2F71=*

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- Do not alter the prosthetic foot or footshell.

5.1 Applying/removing the footshell

INFORMATION

- Pull the protective sock over the prosthetic foot to prevent noises in the footshell.

► **CAUTION! Always use the prosthetic foot with the footshell.**

Apply or remove the footshell as described in the footshell instructions for use.

5.2 Bench Alignment

Bench alignment, transtibial prosthesis

> **Required materials:** 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment apparatus (e.g. 743A200 PROS.A. Assembly)

- 1) Position the prosthetic foot in the alignment apparatus.
→ **Heel height: effective heel height of everyday shoe + 5 mm**
- 2) **Sagittal plane:** Align the prosthetic foot with the alignment reference line.
→ Position the alignment reference line on the markings of the footshell.
- 3) **Frontal plane:** Align the prosthetic foot with the alignment reference line.
→ Position the alignment reference line on the markings of the footshell (**external rotation approx. 5°**).
- 4) Install the structural components between the prosthetic foot and prosthetic socket.
- 5) **Sagittal plane:** Align the prosthetic socket with the alignment reference line.

- Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line.
 - Mark the socket reference point on the centre line: **For transtibial prostheses, at the height of the centre of the patella.** The alignment reference line runs vertically through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion + 5°**
- 6) **Frontal plane:** Align the prosthetic socket with the alignment reference line.
- Take the abduction or adduction of the residual limb into account.

Bench alignment, thigh/knee disarticulation prosthesis

- ▶ Observe the instructions for use for the prosthetic knee joint.

5.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- **Check the length of the prosthesis:** Compression of the prosthetic foot and other components (e.g. a vacuum pump) may require correction of the length.

The alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) can be requested from Ottobock as needed.

5.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

Walking and standing with Evanto looks and feels different. The user typically immediately feels increased flexibility when walking and standing.

Recommended action

- Let the user walk for a few minutes to get used to the new feel and flexibility of the Evanto. In doing so, observe the process of heel strike, rollover and toe-off.
- Have the user walk up ramps. Compared to the previously used prosthetic foot, this should be perceived as equally heavy or lighter.
- ▶ **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.
 - If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
 - If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.

5.4.1 Adjusting stiffness

The stiffness of the prosthetic foot can be adjusted by replacing the two heel wedges:

- **Upper heel wedge:** Determines the heel stiffness.
- **Lower heel wedge:** Determines the forefoot stiffness and heel stiffness.

The heel stiffness influences the knee flexion effect (soft heel wedge = less knee flexion effect. Hard heel wedge=more knee flexion effect). Forefoot stiffness influences the rollover behaviour of the prosthetic foot.

The table shows some examples:

Situation	Measures (several measures individually or in combination)
The user feels a gradual rollover (second half of the stance phase).	First insert a harder lower heel wedge. To make the heel softer again, insert a softer upper heel wedge.

Situation	Measures (several measures individually or in combination)
The user feels that the resistance is too high when walking up a slope.	Move the prosthetic foot gradually in the posterior direction.
The user feels a double rocking during load transfer.	Insert a harder upper heel wedge.
The user feels a dead point during the rollover (first half of the stance phase).	Insert a harder upper heel wedge.
Some users (mobility grade 3 and especially mobility grade 2) prefer a foot position slightly further anterior. This reduces the knee flexing effect, gives a secure feeling at heel strike and a stable rollover.	

Replacing heel wedges

Heel wedge degrees of hardness: The heel wedges are marked with numbers. The smaller the number, the softer the heel wedge.

To replace the heel wedges, the strap looped around the heel must be relaxed. The seat of the strap is then checked and the strap is re-tightened.

- 1) Lift off the cover of the prosthetic foot towards the heel (see fig. 3).
- 2) Fully unscrew the strap tensioner; the screw cannot fall out (see fig. 5).
- 3) Remove the lower heel wedge (see fig. 6).
- 4) **Optional:** Replace the upper heel wedge, aligning the recess in the heel wedge upwards (see fig. 7, see fig. 8).

INFORMATION: As soon as the upper heel wedge has been removed, the seat of the strap tensioners can be easily checked from the bottom.

- 5) Insert the lower heel wedge (see fig. 9).
- 6) Check whether the belt lies correctly in the guide:
 - The belt is positioned all around the guide (see fig. 10).
 - The Ottobock logo is centred on the side of the foot (see fig. 11).
 - The strap tensioner is not at an angle (see fig. 12).
- 7) Tighten the strap tensioner (tightening torque: **4 Nm**) (see fig. 13).

6 Cleaning

> **Allowable cleaning agent:** pH-neutral soap (e.g. 453H10 Derma Clean)

- 1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**
Clean the prosthetic foot and footshell with clear fresh water and pH-neutral soap.
- 2) **If present:** Use a toothpick to remove dirt from the contours that drain off the adapters or footshell and rinse.
- 3) Rinse the soap away with clear fresh water. In doing so, rinse the footshell until all dirt has been removed.
- 4) **In case of heavy soiling and noises:** Pull the footshell and protective sock off the prosthetic foot. Clean all products individually with water or compressed air.
- 5) Dry the product with a soft cloth.
- 6) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Direct the patient to observe the following maintenance instructions:
Check the entire prosthesis for damage after each use.
In case of changes in or loss of functionality and damage, do not continue using the prosthesis and have it checked by qualified personnel.

8 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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