

## 4R140

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for use .....	5
<b>FR</b>	Instructions d'utilisation .....	7
<b>IT</b>	Istruzioni per l'uso .....	9
<b>ES</b>	Instrucciones de uso .....	11
<b>NL</b>	Gebruiksaanwijzing .....	13
<b>DA</b>	Brugsanvisning .....	15
<b>NO</b>	Bruksanvisning .....	17
<b>PL</b>	Instrukcja użytkowania .....	19
<b>HU</b>	Használati utasítás .....	21
<b>CS</b>	Návod k použití .....	23
<b>HR</b>	Upute za uporabu .....	25
<b>EL</b>	Οδηγίες χρήσης .....	27
<b>RU</b>	Руководство по применению .....	29

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Last update: 2020-04-02

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

## 1 Indications for use

The 4R140 Ottobock One-Way Valve may **only be used** for transtibial prosthetic fittings. It is not suitable for the use in swimming prostheses.

## 2 Indication

In combination with a PUR- or TPE-coated knee sleeve the 4R140 One-Way Valve serves to create a vacuum in a transtibial socket, which helps to minimize shearing and friction forces.

## 3 Description (Fig. 1)

The 4R140 One-Way Valve ejects excessive air out of the socket during stance phase and closes during the swing phase to prevent further inflow of air. Thanks to this system, friction and shearing forces are minimized.

The 4R140 One-Way Valve consists of:

- |                     |            |
|---------------------|------------|
| (1) One-Way Valve   | 5X115      |
| (2) Securing Washer | 507S37=M10 |
| (3) Knurled Nut     | 502Z24=M10 |
| (4) Dummy Set       | 5X116      |
| (5) Dust Filter     | 5X118      |

## 4 Fabrication

Use the attached nails to affix the 5X116 Dummy to the plaster positive (Fig. 2). Afterwards laminate the socket in the usual way. Drill a hole ( $\varnothing$  10 mm ) in the center of the Dummy area. Apply silicone bonding agent (617H46) to the threaded part of the 5X115 One-Way Valve (Fig. 3) and press into the bore hole. Seal from outside with 617H46 silicone bonding agent (Fig. 4).

Next, screw the 607S37=M10 securing washer and the 502Z24=M10 knurled nut onto the threaded part of the valve (Fig. 5 and 6).

Finally, let the silicone bonding agent dry for about 1 hour. Insert, from inside, a 5X118 dust filter into the hexagon socket of the valve (Fig.7).

In case of Botta cosmetic covers or cosmetic covers coated with Superskin the 502Z24=M10 knurled nut can be sealed in place on the socket.

**Attention: Make sure that the 5X115 One-Way Valve does not come in contact with resin.**

If it is necessary to change the valve, unscrew it from the socket using a 4 mm allen wrench.

## 5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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