


4R182 Harmony P4 modular

4R183 Harmony P4 modular HD

DE Gebrauchsanweisung	4
EN Instructions for use	11
FR Instructions d'utilisation	18
IT Istruzioni per l'uso	25
ES Instrucciones de uso	32
PT Manual de utilização	39
NL Gebruiksaanwijzing	46
SV Bruksanvisning	53
DA Brugsanvisning	60
NO Bruksanvisning	66
PL Instrukcja użytkowania	73
HU Használati utasítás	80
CS Návod k použití	87
SK Návod na používanie	94
TR Kullanma talimatı	100
RU Руководство по применению	107
JA 取扱説明書	115

INFORMATION

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INFORMATION

Date of last update: 2022-04-01

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 4R182 Harmony P4 modular and 4R183 Harmony P4 HD modular mechanical vacuum pumps generate a vacuum in the prosthetic socket that improves the connection with the residual limb. The vacuum pumps have a four-hole connector for mounting, for example with a socket adapter.

While walking, the vacuum pump is alternately compressed and released again. This generates the vacuum and simultaneously dampens vertical impact loads. The design also permits slight torsion of the prosthetic socket relative to the foot.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

- The maximum approved body weight is specified in the technical data (see page 17). The product must be used in TT prostheses only.

2.3 Environmental conditions

Storage and transport
Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
Temperature range: -10 °C to +45 °C (14 °F to 113 °F)
Moisture: relative humidity: 20% to 90%, non-condensing
Chemicals/liquids: fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)
Solids: dust
Prohibited environmental conditions
Chemicals/moisture: salt water, perspiration, urine, acids, soapsuds, chlorine water
Solids: dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

2.5 Qualification

The product may only be integrated into the prosthesis by O&P professionals who have completed training as certified Harmony experts.

- **Fabricating the total surface weight-bearing socket:** Fabricating a total surface weight-bearing socket is required for the use of the product.

2.6 Absolute Contraindications

- Interim prosthesis
- Effectively non-weight-bearing residual limb end
- Neuromas in the area of the residual limb
- User limitations which prevent the safe use of the product

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 11).
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.

- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested life-time.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by reduced prosthetic suspension.

4 Scope of delivery

Fig. 1, item	Quantity	Designation	Reference number
–	1	Instructions for use	–
1	1	Harmony P4 modular	–
2	1	Filter unit	4X632
3	4	Countersunk head screw	501S41=M6x12-12.9
4	4	Countersunk head screw	501S41=M6x16
5	1	For 4R183 only: elastomer rod (red)	4Y347
6	1	Socket connector	2R117
7	1	Connection kit (one hose, one hose connector 1/8")	–

Spare parts	
Designation	Reference number
Set screw	506G3=M8x16
Elastomer rod (red)	4Y347
Elastomer rod (yellow)	4Y348
Harmony P4 housing screw connection (one threaded ring, four set screws M4x6)	4X904
Harmony P4 HD housing screw connection (one HD threaded ring, four set screws M4x6)	4X446
Cover for 4R182/4R183	4X356

5 Preparing the product for use

<p>⚠ CAUTION</p> <p>Incorrect alignment or assembly</p> <p>Risk of injury due to damaged prosthetic components</p> <ul style="list-style-type: none"> ▶ Observe the alignment and assembly instructions.

<p>⚠ CAUTION</p> <p>Improper assembly of the screw connections</p> <p>Risk of injury due to breakage or loosening of the screw connections</p> <ul style="list-style-type: none"> ▶ Clean the threads before every installation. ▶ Apply the specified torque values. ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

5.1 Adjusting the vacuum pump

Vacuum generation and damping of the vacuum pump can be adjusted using the elastomer rod. Plantar flexion and dorsiflexion as well as pronation and supination can be adjusted via the pyramid receiver.

When a socket adapter with rotatable pyramid receiver is installed, this can be used for rotation.

Adjusting the damping

> **Required materials:** Allen wrench 4 mm

- 1) Disassemble existing distal prosthetic components.
- 2) Screw in the adjustment screw in the distal end of the vacuum pump fully to the stop (see fig. 3).
- 3) Loosen the adjustment screw as described in the adjustment table below.
- 4) **Optional:** Fully unscrew the adjustment screw to replace the elastomer rod (see fig. 4). Unscrew the set screws of the pyramid receiver sufficiently in order to do so. After inserting a new elastomer rod, complete the adjustment procedure again.
- 5) Install the distal prosthetic components as described in the respective instructions for use.

Adjustment table (unscrewing the elastomer rod)			
Body weight [kg]	Product	Red rod	Yellow rod
50	4R182, 4R183	3	-
60 to 70		2.5	-
80 to 90		2	-
100		1.5	3
110	4R183	-	3
120 to 130		-	2.5
140 to 150		-	2

The optimal setting depends on the user's level of activity and preferences. **Loosen the adjustment screw of the elastomer rod by a maximum of 4 turns.**

5.2 Installing the vacuum pump

Countersunk head screws that must be used for installing a socket adapter are included with the vacuum pump. The countersunk head screws are chosen according to the socket adapter.

Socket adapter	Compatible countersunk head screw
Pyramid/pyramid receiver, fixed	501S41=M6x12-12.9
Pyramid/pyramid receiver, rotatable	501S41=M6x16

> **Required materials:** torque wrench (e.g. 710D20), 636K13 Loctite 241

- 1) Select the countersunk head screws that are compatible with the socket adapter (see selection table).
- 2) **INFORMATION: If the pyramid receiver or the pyramid is rotatable: Place the pressure plate on the adapter.**
Position the socket adapter.
- 3) Secure the screws with Loctite.
- 4) Tighten the posterior screws first (**12 Nm**), then tighten the anterior screws (**12 Nm**).
- 5) Install the other prosthetic components as described in the respective instructions for use.
- 6) Carry out the static alignment of the prosthesis.

INFORMATION: The length of the vacuum hose can be determined once all prosthetic components are correctly aligned.


5.3 Installing the socket connector

The vacuum pump is connected to the prosthetic socket by a socket connector. A filter is installed on the vacuum pump to protect it against dirt.

- 1) Install the 2R117 socket connector according to the corresponding instructions for use.

- 2) Shorten the hose of the socket connector to the optimal length and connect it to the filter (see fig. 5, see fig. 6, see fig. 7).

5.4 Checking vacuum generation

To check for leaks, the prosthesis is divided into areas.

- **Area 1:** Prosthetic socket and sealing sleeve
- **Area 2:** Socket connector and connection to the vacuum pump
- **Area 3:** Vacuum pump

Checking area 1

First check the donning process and the sealing sleeve. No pressure gauge is needed. The vacuum suctions the sealing sleeve against the socket brim. The socket brim is clearly visible under the sealing sleeve during the stance phase. If no vacuum is generated or the vacuum decreases, this can be identified on the sealing sleeve.

- Ask the user to put on the prosthesis. Check the process to exclude errors during donning.
- Get the user to walk in order to build up a vacuum. Get the user to stand and check whether the vacuum is maintained.
- **If the vacuum is not maintained:** Replace the sealing sleeve with a new one.
- Get the user to walk in order to build up a vacuum. Get the user to stand and check whether the vacuum is maintained.
- **If the vacuum is still not maintained:** The problem lies in area 2 or 3.

Checking area 3

- Loosen the hose of the socket connector from the filter.
- Connect the pressure gauge to the filter.
- Activate the vacuum pump two to four times until a vacuum of **500 hPa to 850 hPa** is reached.
- **If the vacuum is maintained:** The vacuum pump and filter are OK. The fault is in area 2.
- **If the vacuum is not maintained:** Either the filter is plugged or the vacuum pump is defective.
- Remove the filter from the vacuum pump and connect the pressure gauge directly to the vacuum pump.
- Activate the vacuum pump two to four times until a vacuum of **500 hPa to 850 hPa** is reached.
- **If the vacuum is maintained:** Replace the filter. Then check whether this measure was successful.
- **If the vacuum is not maintained:** Send the vacuum pump to the manufacturer for inspection.

Checking area 2

- Connect the pressure gauge with a T-piece between the socket connector and the vacuum pump.
- Seal the opening to the socket connector on the inside of the prosthetic socket with air-tight adhesive tape (e.g. PVC tape).
- Activate the vacuum pump until a vacuum of **500 hPa to 850 hPa** is reached.
- **If the vacuum is not maintained:** The socket connector is defective. Replace the socket connector. Then check whether this measure was successful.
- **If the vacuum is maintained:** The glued joint of the socket connector is likely defective.
- Seal the opening to the socket connector on the inside of the prosthetic socket with a small piece of paper and air-tight adhesive tape (e.g. PVC tape).
- Activate the vacuum pump until a vacuum of **500 hPa to 850 hPa** is reached.
- **If the vacuum is not maintained:** Glue the socket connector again. Then check whether this measure was successful.

6 Use

⚠ CAUTION

Use of powder (e.g. baby powder or talcum)

Skin irritation on the residual limb and loss of function of prosthetic components due to clogging with particles or deprivation of lubricant

- ▶ Keep the product away from powders.

6.1 Donning the prosthesis

- 1) Roll the liner over the residual limb. Avoid wrinkles, air pockets and soft tissue displacement as you do so.
- 2) **For liners without a textile coating:** Put on a suitable stocking (e.g. 451F21).
- 3) Slide into the prosthetic socket and unroll the sealing sleeve over the prosthetic socket up to the thigh.

6.2 Doffing the prosthesis

- 1) Unroll the sealing sleeve from the thigh onto the prosthetic socket.
- 2) Doff the residual limb and liner out of the prosthetic socket.

7 Cleaning

7.1 Cleaning the prosthetic socket

- ▶ Clean the product with a damp, soft cloth and then let it dry.

7.2 Cleaning the liner

- ▶ Clean the liner as described in the liner instructions for use.

7.3 Cleaning the vacuum pump

INFORMATION

Replace the filter when cleaning the vacuum pump.

INFORMATION

If you cannot operate the vacuum pump:

- ▶ Determine the current setting (number of turns of the adjustment screw to the stop).
- ▶ Release the elastomer rod by five turns before rinsing.
- ▶ Restore the original state after rinsing.

For qualified personnel: cleaning with compressed air

- 1) Loosen the hose of the socket connector from the filter.
- 2) With compressed air, carefully blow through the hole in the prosthetic socket to remove dirt from the socket connector.
- 3) Remove the filter from the the vacuum pump.
- 4) With compressed air, carefully blow through the suction connection of the vacuum pump to remove dirt from the vacuum pump and one-way valve.

For qualified personnel: rinsing with distilled water

- 1) Remove all distal prosthetic components.
- 2) Pull the filter off the vacuum pump.
- 3) Place a towel on the floor or a table and set the vacuum pump onto it.
- 4) Fill a container with approx. **60 ml** distilled water.
- 5) Connect a hose to the vacuum pump and immerse it in the container.

- 6) Operate the vacuum pump. The water is suctioned from the container and comes out of the outlet opening.

INFORMATION: Operate the vacuum pump until no more water comes out of the outlet opening.

- 7) Reconnect the filter to the vacuum pump.
8) Reinstall the distal prosthetic components.

8 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ Perform a visual inspection of the elastomer rod annually. Replace the elastomer rod in case of limited functionality or pronounced signs of wear (see page 14).

9 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

11 Technical data

Reference number	4R182	4R183
Weight [g]	550	675
System height [mm]	125	
Build height [mm]	107	
Body weight [kg]	50 to 100	50 to 150
Mobility grade	2 to 4	

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