1P9

DE	Gebrauchsanweisung	2
EN	Instructions for use	6
FR	Instructions d'utilisation	10
IT	Istruzioni per l'uso	14
ES	Instrucciones de uso	18
PT	Manual de utilização.....	22
NL	Gebruiksaanwijzing.....	26
SV	Bruksanvisning.....	30
DA	Brugsanvisning.....	34
PL	Instrukcja użytkowania.....	38
HU	Használati utasítás	42
CS	Návod k použití.....	46
RO	Instrucțiuni de utilizare.....	49
TR	Kullanma talimatı.....	53
EL	Οδηγίες χρήσης	57
RU	Руководство по применению	61

INFORMATION

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INFORMATION

Date of last update: 2020-03-31

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

1 Application

The 1P9 Pirogoff Foot is to be used **solely** for the prosthetic fitting of the lower limb.

1.1 Safety

1.1.1 Explanation of warning symbols

CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

1.1.2 General safety instructions

CAUTION

Use of a worn product and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- Inspect the product regularly for wear. Do not continue using a worn product.
- Only use the product for a single patient.

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- Only combine the product with prosthetic components that are approved for that purpose.
- Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- Do not expose the product to unallowable environmental conditions.
- If the product has been exposed to unallowable environmental conditions, check it for damage.
- If damage is apparent or in case of doubt, do not continue using the product.
- Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage

- Use caution when working with the product.
- If the product is damaged, check it for proper function and readiness for use.
- In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover characteristics) is an indication of loss of functionality.

2 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

Other environmental conditions

UV light

Not resistant

Hydrolysis

Not resistant

3 Handling

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

NOTICE

Grinding the prosthetic foot

Premature wear resulting from damage to the prosthetic foot

- Do not grind the prosthetic foot.

We recommend using 636W17 Adhesive and 636W26 Hardener to bond the Pirogoff Foot to the ankle block.

First, mix the adhesive with 4 % hardener. Apply this mixture to both bonding surfaces and let the surfaces air dry for 5–8 minutes. Firmly press the ankle and foot together for at least one hour. The highest degree of bonding strength will be reached after 24 hours of curing with the parts pressed together.

Tools used in this process may be cleaned using 634A23 Thinner.

4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

4.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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