


4R1

DE Gebrauchsanweisung	3
EN Instructions for use	6
FR Instructions d'utilisation	10
IT Istruzioni per l'uso	13
ES Instrucciones de uso	17
PT Manual de utilização	21
NL Gebruiksaanwijzing	24
SV Bruksanvisning	28

INFORMATION

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INFORMATION

Date of last update: 2021-12-10

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 4R1 sliding adapter is used to optimise the alignment of leg prostheses. The adapter can be repositioned in the frontal plane and the sagittal plane via the dovetail guides while the patient's weight is on the prosthesis. Adjustments are reproducible using scales. The adapter is removed again before finishing the prosthesis.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Approved for a body weight of **up to 100 kg**.

The product may be used **only** for fitting in a workshop.

2.3 Environmental conditions

Storage and transport
Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
Temperature range: -10 °C to +45 °C (14 °F to 113 °F)
Moisture: relative humidity: 20% to 90%, non-condensing
Unacceptable environmental conditions
Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water
Solids: dust, sand, highly hygroscopic particles (e. g. talcum)

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 7).
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Adapter	–

5 Use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

5.1 Mounting the adapter

Socket adapters are pre-assembled on the sliding adapter. It has a proximal pyramid receiver and a distal pyramid. The socket adapters can be replaced. The sliding adapter is installed in the prosthesis during the bench alignment (e.g. between the prosthetic knee joint and prosthetic socket or between the prosthetic foot and prosthetic socket).

Optional: Replacing socket adapters

- 1) Remove the existing socket adapter.
- 2) Install the socket adapter according to its instructions for use.

Connecting the pyramid adapter and pyramid receiver

> **Required materials:** Torque wrench (e.g. 710D20)

- 1) Install the set screws and tighten them with the torque wrench (**10 Nm**).
- 2) Replace any set screws that are protruding or are recessed too much with appropriate set screws (see selection chart).

Selection table for set screws	
Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

5.2 Adjusting adapter

The sliding adapter is used during the static and dynamic fitting. The position is adjusted with the adjusting screws and secured with cap screws.

- 1) Loosen the two cap screws (see fig. 1).
- 2) **NOTICE! To avoid increased wear and tear, loosen the cap screws before adjusting the adapter.**
Use the adjusting screws to adjust the position in the frontal plane and in the sagittal plane (see fig. 2).
- 3) Tighten the cap screws **slightly** to secure the position of the adapter (see fig. 1).

5.3 Removing the adapter

After finishing the static alignment and the dynamic fitting, the adapter is removed from the prosthesis and replaced with the definitive prosthetic components. Ottobock recommends the 743A16 transfer device for transferring the alignment parameters.

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10 Technical data

Reference number	4R1
Weight (without socket adapter) [g]	465
Weight (with socket adapter) [g]	610
Height without adapter [mm]	46

Reference number	4R1
Height with adapter [mm]	82
Build height [mm]	68
System height [mm]	68
Max. displacement medial/lateral [mm]	±15
Max. displacement anterior/posterior [mm]	±25
Max. body weight [kg]	100

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