


## 17B95, 17B96

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	10
<b>FR</b> Instructions d'utilisation .....	16
<b>IT</b> Istruzioni per l'uso .....	23
<b>ES</b> Instrucciones de uso .....	30
<b>PT</b> Manual de utilização .....	37
<b>NL</b> Gebruiksaanwijzing .....	43
<b>SV</b> Bruksanvisning .....	49
<b>DA</b> Brugsanvisning .....	55
<b>NO</b> Bruksanvisning .....	61
<b>PL</b> Instrukcja użytkowania .....	68
<b>CS</b> Návod k použití .....	74
<b>TR</b> Kullanma talimatı .....	80
<b>EL</b> Οδηγίες χρήσης .....	87
<b>RU</b> Руководство по применению .....	94

1

## INFORMATION

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### INFORMATION

Date of last update: 2021-05-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17B95\* and 17B96\* orthotic knee joints.

## 2 Product description

### 2.1 Components/design

Scope of delivery (see fig. 1)			
Item	Quantity [Piece (s)]	Designation	Reference number
1	1	Upper joint section	
2	4	Screws	501T7*
3	1	Joint screw	501S32*
4	1	Lower joint section	
6	1	Bearing nut	17Y93*

Scope of delivery (see fig. 2)			
Item	Quantity [Piece (s)]	Designation	Reference number
1	1	Perlon cable	21A18=2
2	1	Locking rocker	30Z22*
3	1	Notch pin	506K2*
4	1	Countersunk screw	501S41*
5	1	Threaded sleeve	21A12
6	1	Bushing	17Y89
7	1	Lock wedge	17Y88*
8	1	Compression spring	513D19*

Not included in scope of delivery			
Item	Quantity [Piece (s)]	Designation	Reference number
(see fig. 1), item 5	1	Side bar	17B38*, 17B39*, 17B40*, 17B41*, 17B84*
Not illustrated	1	Lamination bar	17Y128*, 17Y129*

## 3 Intended use

### 3.1 Indications for use

The orthotic knee joints are to be used **exclusively** in pairs to provide orthotic treatment for a patient's lower limbs.

Only use the lamination bars for the lamination and prepreg technique.

Only use the side bars for the thermoplastic and bar/clamp technique.

#### INFORMATION

##### Increased wear and tear

The orthotic knee joints are intended only for use in combination with free-motion or foot lifting orthotic ankle joints. Wear occurs more rapidly in combinations involving orthotic ankle joints with a dorsal stop.

### 3.2 Indications

- Partial or total paralysis of the leg muscles
  - Orthopaedic diseases of the lower limbs
- Indications must be determined by the physician.

### 3.3 Restrictions for use

Only use the 17B95=L/R16 and 17B96=16 orthotic knee joints up to a body weight of **max. 50 kg**.

### 3.4 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

## 4 Safety

### 4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

### 4.2 Safety instructions

#### Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

**⚠ CAUTION**

**Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**⚠ CAUTION**

**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

**⚠ CAUTION**

**Improper use of the lock**

Risk of injury due to unexpected deactivation of the orthotic joint lock

- ▶ Fully extend the orthotic joint to activate the lock. Check whether the lock is engaged.
- ▶ Only deactivate the lock when you are able to control the flexion of the orthotic joint.

**⚠ CAUTION**

**Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

**NOTICE**

**Thermal overloading of the product**

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

## NOTICE

### Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

## 5 Preparing the product for use

### INFORMATION

#### Parallel alignment of the orthotic joints

Use the 743R6 alignment fixture to position the orthotic joints on the plaster positive.

### INFORMATION

#### Corrosion protection

The side bars and orthotic joints are resistant to corrosion. Fine grinding and polishing increase corrosion protection. Use sintering powder (e. g. 618T40) for coating.

### Gluing in bars

**Required materials:** 636W18 special glue and 636W19 hardener or 636W28 special glue, degreasing cleaner

The bars must be glued in the insertion zones to ensure stability.

- 1) Clean the connecting surfaces with a degreasing cleaner.
- 2) Mix the special glue and the hardener.
- 3) Apply the mixture to the insertion zones.
- 4) Insert the bars.
- 5) Secure with the enclosed screws.
- 6) Allow to cure for at least **4 hours**.

**INFORMATION: Final bonding strength will be reached after 16 hours.**

## Connecting the perlon cable

- 1) Cut the perlon cable to length according to the patient.
- 2) Screw the threaded sleeve onto the first perlon cable (see fig. 3, item 4).
- 3) Thread the second perlon cable through the coupling piece.
- 4) Then slide the spring over the perlon cable and screw on the threaded sleeve.
- 5) Unscrew the preassembled threaded sleeve from the perlon cable.
- 6) Pass the perlon cable through the bore hole in the locking rocker and screw on the threaded sleeve (see fig. 2, item 1).

## Using protective devices

- ▶ Protect the product against dirt: Use suitable protective devices (e.g. joint protectors).

## 6 Cleaning

Promptly clean the product after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse system joints with clean fresh water.
- 2) Dry system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

## 7 Maintenance

### CAUTION

#### **Wear and tear of the locking mechanism**

Injuries due to changes in or loss of functionality

- ▶ In case of wear and tear of the locking mechanism, replace the lock wedge or the lock lever and/or the lower joint section.

### INFORMATION

#### **The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.
- ▶ **Compression spring:** Replace the compression spring at least once a year or if there are signs of wear and tear.

## Replacing the compression spring

Disassemble the locking rocker to replace the compression spring (see fig. 2).

- 1) Remove the notch pin and unscrew the countersunk screw.
- 2) Remove the bushing and the locking rocker from the upper joint section.
- 3) Replace the compression spring.
- 4) Position the bushing and the locking rocker on the upper joint section.
- 5) Press in the notch pin.
- 6) Install the countersunk screw and secure it with Loctite 241.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Only use 633F7 special lubricant.

Spare parts are listed in the section "Components/design".

## 8 Disposal

Dispose of the product in accordance with national regulations.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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