

17B97=*


DE	Gebrauchsanweisung	5
EN	Instructions for use	12
FR	Instructions d'utilisation	19
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	34
PT	Manual de utilização	42
NL	Gebruiksaanwijzing	49
SV	Bruksanvisning	57
DA	Brugsanvisning	63
NO	Bruksanvisning	70
FI	Käyttöohje	77
PL	Instrukcja użytkowania	83
HU	Használati utasítás	91
CS	Návod k použití	98
SK	Návod na používanie	105
TR	Kullanma talimatı	112
EL	Οδηγίες χρήσης	119
RU	Руководство по применению	127
JA	取扱説明書	135


1

2

INFORMATION

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INFORMATION

Date of last update: 2020-04-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17B97* orthotic knee joints.

2 Product description

2.1 Available sizes

Reference number	System width
17B97=L16 17B97=R16	16 mm
17B97=L20 17B97=R20	20 mm

2.2 Components

Scope of delivery (see fig. 1)			
Item	Quantity [Piece (s)]	Description	Reference number
1	4	Oval head countersunk screws	501T7=7.5x9xM5
2	2	Orthotic knee joints	17B23K=*
3	4	Lamination bars	17Y128=*
4	2	Compression spring	513D19=3.8x16
5	2	Lamination dummy	17Y103
6	6	Double hollow rivet	504H1=7-100
7	1	Perlon cable	21A18=2
8	1	Perlon cord	21A18=3
9	1	PVC profile material	17Y106=*

Scope of delivery (see fig. 2)			
Item	Quantity [Piece (s)]	Description	Reference number
1	1	Coupling piece	21A7
2	1	Perlon cable	21A18=2
3	1	Spring	21A25
4	2	Threaded sleeve	21A12

Not included in scope of delivery			
Item	Quantity [Piece (s)]	Description	Reference number
not illustrated	1	Special adhesive	636W18
not illustrated	1	Hardener	636W19
not illustrated	1	Orthotic joint alignment fixture	743R6
not illustrated	1	Plastaband	363K8=20x2x10

Not included in scope of delivery			
Item	Quantity [Piece (s)]	Description	Reference number
not illustrated	1	Lamination bar, 100 mm	17Y128=16x100/=20x10-0
not illustrated	1	Lamination bar with contoured calf	17Y129=L16/=R16/=L2-0/=R20

2.3 Design

The 17B97* components facilitate the fabrication of a leg orthosis using laminating technology and prepreg technology. The orthotic knee joints with a concealed lock and short lock lever are designed especially for unlocking with a pull-release cable. The perlon cable is positioned with the threaded sleeve in the blind hole of the lock lever.

The orthotic joints have a detachable connection with the lamination bars.

The 17Y106 PVC profile material is used to establish a connection between the two 17Y128 lamination bars.

3 Intended use

3.1 Indications for use

The orthotic knee joints are to be used **exclusively** in pairs to provide orthotic treatment for a patient's lower limbs and are intended for processing using lamination technology.

3.2 Indications

Partial or total paralysis of the leg muscles.

Indications must be determined by the physician.

3.3 Combination possibilities

The 17B97* orthotic knee joints can be combined with the 17B98* orthotic ankle joints (see 647H167 instructions for use).

4 Safety

4.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 General safety instructions

CAUTION

Incorrect alignment, assembly or adjustment

Injury and damage to product components

- ▶ Observe the alignment, assembly and adjustment instructions.
- ▶ Work on the product may only be carried out by qualified personnel.

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

CAUTION

Improper handling

Risk of falling due to accidentally opening the lock

- ▶ Prior to handing over the orthosis, ensure that the locks are engaged securely and in parallel.
- ▶ Avoid accidentally opening the locks.
- ▶ Instruct the patient in the correct usage of the system joints.

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

5 Preparing the product for use

⚠ CAUTION

Repeated bending at the same position

Risk of injury due to breakage of the bar, change in or loss of functionality

- ▶ Avoid repeated bending at the same position.

NOTICE

Orthotic joints not aligned in parallel

Malfunction of the locking mechanism

- ▶ Use the 743R6 orthotic joint alignment fixture.
- ▶ Make sure the orthotic joints lock simultaneously in the extension stop.

Preparing the orthotic knee joint

- 1) Secure the lamination bars in the orthotic knee joint using the included screws.
- 2) Remove the lock lever and insert the lamination dummy (see fig. 3).
- 3) **INFORMATION: The lamination dummy fixes the orthotic knee joint in the extended position and prevents resin from collecting inside the joint.**
- 4) Install the orthotic knee joint on the size 3 alignment insert (see fig. 4).

- 5) Use 636K8 Plastaband to seal the milled area of the bearing nut, the screw heads and the spring channel of the upper joint section.

Shaping the bars

The 17B97* kit contains four lamination bars with a straight design and a length of 80 mm. The lamination bars are optionally available with a length of 100 mm and a contoured calf. The lamination bars are profiled at the ends for better anchoring in the laminate (see fig. 3). The sharpened part with the bore hole (4.2 mm) supports the attachment to the 17Y106 PVC profile material.

- 1) Secure the bars on the orthotic knee joint using the screws and check the alignment on the model (see fig. 5).
- 2) Mark the insertion zone on the bar.
- 3) Loosen the screws of the bars and adapt the bars with bending irons.
- 4) Fasten the shaped bars to the orthotic knee joint using the screws.
- 5) Repeat the steps until the desired result is obtained.
- 6) **INFORMATION: When fabricating the orthosis, use 636W18 special adhesive and 636W19 hardener (not included in the scope of delivery) to bond and screw together the lamination bars in the insertion zone of the orthotic knee joint.**

Processing the PVC profile material

Use the PVC profile material as core material for the reinforcement and the reinforcement clamps. It serves as the connecting element between the orthotic knee joint and orthotic ankle joint.

- 1) Mark the position of the PVC profile material on the plaster model and cut it to length.
- 2) Heat the PVC profile material with a hot air gun (approx. 80 °C) and mould it into shape.
- 3) **INFORMATION: The material can also be shaped slightly while cold (see fig. 6).**
- 4) Use the included hollow rivets to connect the shaped PVC profile material to the lamination bars (see fig. 7). **Optional:** Use 636K11 quick-drying adhesive.

Using the perlon cord

The perlon cord (3 mm) is integrated in the reinforcement of the orthosis and removed after the laminate has hardened. The perlon cable (2 mm) is routed in the fabricated cable channel for the lock lever.

- 1) Determine the length of the perlon cord.
- 2) Heat the perlon cord with the hot air gun (180 °C) and mould it to the model (see fig. 8).

- 3) **INFORMATION: The perlon cord does not need to be isolated from the laminating resin.**
- 4) Cut the beginning and end of the perlon cord free and remove it from the laminate (see fig. 9).

Connecting the perlon cable

The coupling piece is used to give the ends of the perlon cables a detachable connection so the lock can be released.

- 1) Cut the perlon cable to length according to the patient.
- 2) Screw the threaded sleeve onto the first perlon cable (see fig. 2, item 4).
- 3) Thread the second perlon cable through the coupling piece.
- 4) Then slide the spring over the perlon cable and screw on the threaded sleeve.
- 5) Secure the perlon cable and preinstalled threaded sleeve in the blind hole of the lock lever.
- 6) Check that the threaded sleeve is securely positioned in the blind hole of the lock lever and use 636K14 Loctite to prevent it from being pushed out if desired.

6 Cleaning

Clean the orthosis joint promptly after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the orthosis joint with clean, fresh water.
- 2) Dry the orthosis joint with a cloth or allow it to air-dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer recommends a semi-annual inspection to verify functionality and check for wear and tear to the orthotic joints.

Spare parts are listed in the section "Components/design".

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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