
















## 8E44=6, 8E44=7

<b>DE</b> Gebrauchsanweisung .....	6
<b>EN</b> Instructions for use .....	14
<b>FR</b> Instructions d'utilisation .....	21
<b>IT</b> Istruzioni per l'uso .....	29
<b>ES</b> Instrucciones de uso .....	37
<b>PT</b> Manual de utilização .....	45
<b>SV</b> Bruksanvisning .....	53
<b>NL</b> Gebruiksaanwijzing .....	61
<b>NO</b> Bruksanvisning .....	69

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Last update: 2021-04-01

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- Please keep this document in a safe place.

## 1 Application

The 8E44=6 resp. 8E44=7 Transcarpal Hand is to be used **solely** for the fitting of upper limb prostheses. It allows for the fitting of very long residual limbs, including transcarpal amputations.

The 8E44=\* Transcarpal Hand is designed **exclusively** for use on one patient. Use of the product by other persons is not approved by the manufacturer.

## 2 Safety precautions

Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

- Do not use silicone spray when putting on the prosthetic glove, because silicone spray may compromise secure fit of the glove. We recommend using 633S2 Ottobock Procomfort Gel as a mounting aid.
- In order to reduce the malfunctioning risk of the Transcarpal Hand and the connected components due to strong electromagnetic radiation (e.g. when near high-tension power lines, transmitters, transformers or security systems for goods in department stores), the electrode gains should be set as low as possible.
- Before interrupting or establishing electrical connections (such as when removing the hand from the prosthesis), the system must be disconnected from the power supply. For this purpose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.
- Modifying component parts without the written consent of the manufacturer is strictly prohibited.
- The Transcarpal Hand has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Myo-Service.
- If you do not achieve the expected results by adjustment and by selection of the appropriate control mode, please contact an Ottobock Myo-Service facility.
- Use the "Patient Information" section of this pamphlet (Section 3) to inform your patient about the function and use of the Transcarpal Hand.
- Please provide your patient with the 646D44 Patient Information.

### 3 Patient information

- Do not let foreign particles or liquids enter the Transcarpal Hand.
- The Transcarpal Hand should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Avoid staying near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores), as this can lead to malfunction of the Transcarpal Hand.
- If the hand is not in use, it should be left in the open position to protect the mechanism.
- The Transcarpal Hand and Transcarpal Hand components may be opened or repaired only by certified Ottobock Myo-Service technicians.
- Keep the product away from the immediate vicinity of active implantable systems (e.g. pacemakers, defibrillators etc.) since electromagnetic radiation may cause temporary interference.
- If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour. Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.99 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points,...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

#### 3.1 Driving motor vehicles

**Attention:**

An upper limb amputee's ability to drive a motor vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb condition, design of the prosthesis) and the amputee's abilities. Drivers must observe driving laws applicable in their country. For insurance purposes, we strongly recommend having your driving ability examined and confirmed by an authorized agency.

For maximum safety and convenience, Ottobock recommends that a specialist evaluate the need to adapt the vehicle with specialized equipment. It is indispensable that the user can drive the motor vehicle safely with the Transcarpal Hand switched OFF. Driving with the Transcarpal Hand switched ON could pose a risk to public traffic in case of improper control of the Hand.

### 4 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product. The grace period is no more than one month before or three months after maintenance is due. In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover. Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty. The following components must always be sent in for maintenance and repairs: The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

## 5 Components (Fig. 1):

- (1) 1 pc. 8E44=6 Transcarpal Hand DMC plus/8E44=7 Transcarpal Hand Digital Twin
- (2) 1 pc. lamination plate
- (3) 1 pc. lamination cover
- (4) 1 pc. system inner hand
- (5) 1 pc. cable lamination cover
- (6) 1 pc. 7 mm special ring key
- (7) 1 pc. connecting cable
- (8) 1 pc. cable insulation
- (9) 1 pc. special thread for system inner hand

## 6 Basics

### **Selection of Control Mode** (Fig. 7; Fig. 8; Fig. 9)

In order to achieve optimal customization, two program variations are available for the DMC plus Control as well as for the Digital Twin control. These programs can be selected by inserting or removing the 13E185 Coding Plug. Before choosing the Control Mode partially pull down the Inner Hand in order to access the Coding Plug.

Afterwards pull the Inner Hand back on. To reset the system so that it will recognize the new control mode, you must remove the battery, then re-insert it.

Descriptions of the operating mode functions can be found with their respective controls.

### **Motor Connection - Coding Plug Connection**

The Motor is always connected to the electronics at the connection marked 'M', the Coding Plug at the connection marked 'S'. This is independent of left or right hand direction (Fig. 5).

### **ON/OFF Switch**

The ON/OFF Switch (Fig. 6) is operated by the amputee through the Inner Hand and Prosthetic Glove. To avoid accidental release the ON/OFF Switch must be held for approx. 1 second to activate/deactivate the prosthesis.

## 7 Description and function 8E44=6, DMC plus Control

The 8E44=6 Transcarpal Hand is equipped with a proportional DMC plus Control. The DMC plus Control is contained in a **black housing** with a **black adjustment ring**. The hand cable of the hand is marked with a **brown ring**.

**The transcarpal hand with the DMC plus control unit is delivered with the function plug connected in DMC mode.**

For ideal adaptation to the patient, there are two programming variations that can be selected with the 13E185 Function Plug:

### **7.1 Function plug connected** (Factory Settings Fig. 7)

The function is identical to the already familiar DMC Control. The grip speed and grip force are precisely determined by the height of the muscle signals. If the strength of the muscle signals changes, the grip speed and grip force adapt immediately and proportionately to the changed muscle signal.

## 7.2 Coding plug removed (Fig. 9):

Function as described above except that the DMC plus mode includes a "virtual hand switch". After gripping once with maximum force, the signal required to open the hand will be slightly higher. This improves grip security and facilitates eating, among other activities.

## 7.3 Adjusting the 13E200=\*/13E202=\* Electrodes

The best possible function of the DMC plus control is achieved through optimal adjustment of the electrodes. This is easy with the MyoBoy 757M11=X-Change.

**Procedure:** for detailed information see 647G265 MyoBoy Operating Manual.

In the PAULA software, select the Transcarpal Hand DMC plus. Adjust each electrode so that the patient can maintain the corresponding muscle signal for approximately two seconds at a point slightly above the value HIGH. (Fig. 2)

### **Attention:**

Allow the patient to rest during the adjustment of the electrodes to avoid fatigue, and thus incorrect electrode settings (i.e. excess sensitivity).

The electrodes should be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interferences from electric devices, the position of the electrodes should be verified and, if need be, repositioned. If the interference cannot be eliminated, please contact the Ottobock Myo-Service.

## 8 Description and function 8E44=7, Digital Twin Control

The Digital Twin Control includes a Digital and a Double-Channel Control Mode. The chosen control mode is selected with the 13E185 Coding Plug. The Digital Mode is controlled with two electrodes and the Double-Channel Mode with only one electrode (13E200=\* or 13E202=\*).

The Digital Twin circuit is placed in a **black housing with white adjustment Controls**.

The hand cable of the hand is marked with a **white ring**.

**The Transcarpal Hand Digital Twin is delivered with a Digital Mode Coding Plug installed.**

### 8.1 Control modes: Digital Twin Control

**Coding Plug inserted: Digital Mode** (as delivered, Fig.7)

The opening and closing speed of the Transcarpal Hand is constant after a threshold value has been exceeded by the corresponding electrode signal.

**Coding Plug removed: Double-Channel Mode** (Fig. 9)

The hand is opened / closed with only one electrode.

A long, weak muscle signal → Transcarpal Hand closes.

A short, strong muscle signal → Transcarpal Hand opens.

### 8.2 Adjusting the 13E202=\*/13E200=\* Electrodes

The best possible function of the Digital Twin Control is achieved through optimal adjustment of the electrodes. This is easy using the MyoBoy 757M11=X-Change.

**Procedure for use with Digital Mode (with Coding Plug)**

(Detailed information: see 647G461 instructions for use)

In the PAULA software, select the Transcarpal Hand Digital Twin.

Adjust each electrode so that the patient can maintain the corresponding muscle signal for approximately two seconds at a point slightly above the value ON. (Fig. 3)

## Procedure for use with Double-Channel Mode (= without Coding Plug)

(Detailed information: see 647G461 instructions for use)

In the PAULA software, select the Transcarpal Hand with dual-channel control. Position the indicator to 3.5 on the scale. The switch's Adjustment Control 2 is already pre-set to the corresponding value (Fig. 5).

Now the electrode has to be adjusted so that the second threshold – identified in the PAULA software by the black line at 3.5 – is barely reached with the muscle signal (Fig. 4). The Transcarpal Hand closes.

With a quick-and-strong muscle signal, a point above the threshold must be reached in order for the Transcarpal Hand to open.

If the muscle signal is too slow or too weak, use the adjustment control 2 to change the threshold (Fig. 5).

## 9 Battery management

Use only fully-charged batteries to operate the Transcarpal Hand. Always keep a charged, back-up Ottobock battery available as a replacement. An intelligent "battery manager" alerts the patient when the voltage is low: the Hand will become slower and the gripping force will decrease. At the same time, the battery is protected against harmful total discharge.

We recommend powering the Transcarpal Hand with the 757B20/21 EnergyPack or 757B35=\* MyoEnergy Integral to benefit from the Hand's full capacity over an extended period of time. Subject to functional limitations, the transcarpal hand can be operated with the 757B15 X-ChangePack. For detailed instructions on the use of the batteries see the information leaflet included in the batteries.

### Attention:

As soon as the Transcarpal Hand recognizes a **fully charged 757B20/757B21 EnergyPack or 757B35=\* MyoEnergy Integral**, it automatically switches to Li-Ion battery mode. Subsequently, only this technology should be used to operate the system. If an Ottobock 757B15 X-ChangePack is used anyway, its capacity can no longer be fully utilised. The 757T13 Ottobock MyoSelect can be used, however, to reset the Transcarpal Hand control for use with the X-ChangePack.

**For reasons of operational safety and reliability, use only 757B20, 757B21, 757B35=\* or 757B15 Ottobock Batteries.**

Upon delivery, the Transcarpal Hand is set for use with NiMH batteries.

## 10 Cleaning and Care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Make sure the product is not damaged and that no liquid can penetrate into system components.

Then, dry the prosthesis component with a soft cloth.

## 11 Disposal

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

## 12 Technical data

Average current consumption	250	mA
Electronic switch-off	720	mA
Static current	~1	mA
Operating temperature	0–70	°C
Opening width	100	mm
Proportional speed	15–130	mm/s
Proportional grip force	0–90	N
Weight (with System Inner Hand)	320	g (size 71/4)
Service life of terminal device	5	years
Battery service life	2	years

### Power supply

EnergyPack	757B20	– 7.2 V Li-Ion technology
EnergyPack	757B21	– 7.2 V Li-Ion technology
MyoEnergy Integral	757B35=*	– 7.4 V Li-Ion technology
X-ChangePack	757B15	– 6 V NiMH technology

### Ambient conditions

Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing

## 13 Symbols used

**MD** Medical device

## 14 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

## 15 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices. The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices. The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>.

# Lamination Instructions

## Transcarpal Hand 8E44=6/8E44=7

Please read the corresponding safety instructions prior to handling and processing 617H12 Pedilen Rigid Foam 200 and 617H19 Orthocryl lamination resin.

1. Manufacture the inner socket (Fig. 10).
2. Pull powdered PVC bag over inner socket and apply vacuum. Glue a 2 cm wide foam cutting strip along trim line.
3. Apply and form 617H12 Pedilen rigid foam 200 to inner socket (Fig. 11).
4. Working with the patient, establish the desired position of the hand using the Pedilen rigid foam. The position and rotation of the hand should have a natural appearance. Only minor adjustments can be made once the final position has been established.
5. Mark the perimeter of the lamination plate on the Pedilen rigid foam.
6. With the help of the included special ring key, separate the hand from the lamination plate.
7. Reduce the marked perimeter by 3 mm to provide enough clearance for the wire loops to fit over the Pedilen rigid foam (Fig. 12). Then shape the distal end of Pedilen rigid foam to this reduced perimeter. Make sure the RA-radius (thumb side) and UL-ulna markings are placed correctly on the model.

The lamination plate must sit level on the foam (Fig. 13a and 13b)!

8. Pull the PVA bag with distally sealed-seam over the model.
9. Cut a double-length of Perlon stockinette, sew it in the middle, and reflect it back over the rigid foam, making sure there are no wrinkles in the fabric (Fig. 14).
10. Screw the cable lamination cover onto the lamination plate. The opening in the middle of the plate comes with a piece of tape covering it (Fig. 15a, Pos. 1). Evenly snap the lamination cover into position (15b).
11. Using pigmented 617H21 Orthocryl Sealing Resin, glue the lamination plate under a load of 2 kg (Fig. 16) to the rigid foam over which 2 layers of perlon stockinettes have been pulled.
12. Cut off a double-length of Perlon stockinette, pull half over the model and tie-off beneath the lamination ring with thread (Fig. 17a). Reflect remaining half of the stockinette as shown in Fig. 17b.
13. If heavier loading of the prosthesis is expected, add additional layers of perlon stockinette or reinforce with carbon.
14. Carefully pull on the moistened PVA bag and begin the lamination process (617H19 Orthocryl Lamination Resin).
15. Once the resin has set, remove the lamination cover from the lamination plate (Fig. 18), and remove the tape. Separate the inner socket from the outer shell.

16. Remove the cable lamination cover. In order to maintain the cable relief, use a hot needle to completely burn out any resin in the opening (Fig. 19).
17. Sand the distal end of the shell smooth (Fig. 20) in order to avoid any damage to the inner shell.
18. Thread the 4-pole hand cable through the lamination plate (Fig. 21) and evenly press in the cable insulation (see cable installation, Fig. 7-Fig. 9).
19. Attach hand to lamination plate. To make adjustments, tighten the nut with the included special ring key.
20. Plug the 4-pole hand cable into the connection block. Secure the nuts with 636W60 Loctite 243.
21. Pull the inner shell over the hand and lamination plate and tie-off in the groove with the included special thread. Seal the slit between the System Inner Hand and the prosthetic socket with 633F11 Silicone Grease.
22. Pull on the prosthetic glove. If it is difficult to pull on, heat the glove and stretch it out.

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