



## 28R14

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	7
<b>FR</b> Instructions d'utilisation .....	11
<b>IT</b> Istruzioni per l'uso .....	15
<b>ES</b> Instrucciones de uso .....	19
<b>PT</b> Manual de utilização .....	23
<b>NL</b> Gebruiksaanwijzing .....	27
<b>SV</b> Bruksanvisning .....	31
<b>DA</b> Brugsanvisning .....	35
<b>PL</b> Instrukcja użytkowania .....	39
<b>HU</b> Használati utasítás .....	43
<b>CS</b> Návod k použití .....	46
<b>TR</b> Kullanma talimatı .....	50
<b>EL</b> Οδηγίες χρήσης .....	54
<b>RU</b> Руководство по применению .....	58

	<b>Größe Size</b>	<b>Taillenumfang Waist circumference</b>
		<b>cm</b>
	28R14=S	60–75
	28R14=M	75–90
	28R14=L	90–105
	28R14=XL	105–120

**INFORMATION**

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#### **INFORMATION**

Date of last update: 2021-10-21

► Please read this document carefully before using the product and observe the safety notices.

- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on adjusting and fitting-on the 28R14 hyperextension brace according to André Bähler.

## 2 Product description

### 2.1 Components

Scope of delivery (see fig. 1)					
Item	Component	Article number	Item	Component	Article number
1	Sternal support	29P13	6	Protective clip	29R145
2	Chest bar	29R138	7	Closure	29R137
3	Base plate	29R135	8	Back support	29P57
4	Stomach bar	29R139	9	Waistband	29R140
5	Symphysis support	29P14	10	Closing aid	29R151

Accessories	
Chest bar, short	29R138=200
Chest bar, long	29R138=350

### 2.2 Design

The product is ready for fitting when delivered, and can be adjusted to the patient by qualified personnel.

The ventral part consists of plastic-coated aluminum and has special ratchet closures. These accommodate the flexible waistband that is attached to the back support.

The sternal support and the symphysis support have an aluminum core that can be reshaped. The length of the sternal support can be adjusted by maximum 20 mm over the chest bar. It can be pushed freely along the chest bar.

## 3 Intended use

### 3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 3.2 Indications

- Stable traumatic fractures of the vertebral bodies in the lower thoracic and lumbar spine, without neurological deficits
- Osteoporosis

The indication must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

None known.

### 3.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the device is fitted, tumours, lymphatic congestion – including unclear soft tissue swelling distal to the device, and dysaesthesia in the trunk.

### 3.4 Lifetime

The product is designed for a maximum lifetime of **2 years**.

### 3.5 Mechanism of Action

The brace straightens the spine in the sagittal plane and stabilizes it. Through the back support, the sternal support, and the symphysis support, the three-point principle causes lordosing, relieving in this way the anterior vertebral bodies.

### 3.6 Restrictions for use

The product is designed exclusively for fitting on one patient.

Parts that are shaped, and parts that come directly into contact with the skin, could cause functional or hygienic risks if the product is used on another person:

- when cracks in the material form that cannot be visually detected if paint is undamaged
- in the case of material fatigue due to alternating loads
- through a decrease in elongation at break (cold fixing) in the case of multiple cold forming

Reuse of the product on other persons is not permitted.

The daily duration of use and application is based on the medical indication.

## 4 Safety

### 4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

### 4.2 General safety instructions

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

#### Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

## 5 Selecting the Size

- 1) Measure the waist circumference approx. **2 cm** above the iliac crest.
- 2) Determine the brace size (see size chart).

## 6 Fitting and Application

### CAUTION

#### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

### CAUTION

#### **Loosening the ratchet closures**

Loss of the support function, injuries due to wear or unintentional release of the ratchet closures

- ▶ At each application, ensure that the ratchet closures audibly engage in the waistband. The waistband must not slip out under load.
- ▶ Practice the application and removal of the product with the patient.
- ▶ In case of damage or noticeable changes in the product, the patient should consult a specialist workshop.

### INFORMATION

Continuous, direct skin contact with the product can lead to skin irritation. Only wear the product over clothing (e.g. a T-shirt or undershirt) in daily use. For limited times (e.g. water therapy), the product can also be worn directly on intact skin.

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

Ottobock recommends arranging the metal parts by hand, or using bending irons with rounded edges. Do not bend the adjustment area of the slotted parts or the guiding element in the base plate.

- 1) To adapt the length, release the screws inside the supports and the base plate (see fig. 2).
- 2) Adapt the sternal support and symphysis support to the height of the patient while observing the extension length. The length of the sternal support can be adjusted by around 20 mm over the chest bar. It can be pushed freely along the chest bar.
- 3) **Option:** A short chest bar 29R138=200 and a long chest bar 29R138=350 are available as accessories. To exchange the chest bar, remove the 3 screws of the base plate and reassemble the sternal support (see fig. 2, see fig. 3).
- 4) If pushing the slot is not enough for a short torso, the chest bar can be shortened (see fig. 4).
- 5) After adapting, check and tighten the screws.
- 6) Shorten the pad of the back support on the embossed cut edges to the desired circumference.
- 7) The brace can be fitted on without any problem via the ratchet closures. To open the brace, actuate one of the pushbuttons on the ratchet closure. The closing aid helps orientation when finding the correct setting of the waistband (see fig. 5).
- 8) Check the correct fit of the brace and adjust if necessary.

## **When the brace is fitted-on, but has not yet been handed over to the patient, check**

- Pressure, placement and fit of the 3 supports
- Tension and fitting of waistband and base plate
- Fit and hyperextension effect in the sitting position

### **INFORMATION**

During the first fitting, pay special attention that the waistband latches into the ratchet closure, as this helps prevent the waistband from slipping out under pressure.

## **7 Cleaning**

### **NOTICE**

#### **Contact with salt water or water containing chlorine as well as abrasive substances (e.g. sand)**

Damage and premature product wear and tear

- ▶ Promptly clean the product following any contact with the substances identified above.

The brace is not damaged by water, which means that personal hygiene is not restricted. Showering is possible when wearing the brace. Clean the brace afterwards:

- 1) Clean the brace with a damp cloth and neutral soap.
- 2) Rinse off the brace with water and dry with a cloth.

## **8 Disposal**

Dispose of the product in accordance with national regulations.

## **9 Legal Information**

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### **9.1 Liability**

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

### **9.2 CE Conformity**

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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