50A4 Rhizo Forsa

| | |
|--|----|
| DE Gebrauchsanweisung..... | 4 |
| EN Instructions for use..... | 9 |
| FR Instructions d'utilisation | 14 |
| IT Istruzioni per l'uso..... | 20 |
| ES Instrucciones de uso | 25 |
| PT Manual de utilização..... | 31 |
| NL Gebruiksaanwijzing | 36 |
| SV Bruksanvisning | 41 |
| DA Brugsanvisning..... | 46 |
| PL Instrukcja użytkowania | 51 |
| HU Használati utasítás | 57 |
| CS Návod k použití | 62 |
| TR Kullanma talimatı..... | 67 |
| EL Οδηγίες χρήσης | 72 |
| RU Руководство по применению | 78 |
| ZH 使用说明书 | 84 |

1

50A4

2

3

4

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INFORMATION

Last update: 2016-01-29

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 50A4 Rhizo Forsa carpometacarpal joint orthosis.

2 Application

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the upper limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Collateral ligament injuries of the metacarpophalangeal joint (e.g. skier's thumb)
- Sprains
- Irritation and wear of the carpometacarpal and metacarpophalangeal joints (e.g. osteoarthritis of the base of the thumb)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent, swollen scars, reddening and hyperthermia in the relevant area; pronounced varicose veins, especially with impaired blood return, lymph drainage disorders and unclear soft tissue swelling distal to the support; sensory and circulatory disorders in the area of the forearm and hand, e.g. associated with diabetic neuropathy.

2.4 Effects

The orthosis stabilises the carpometacarpal and metacarpophalangeal joints.

3 Safety

3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

3.2 General Safety Instructions

⚠ CAUTION**Contact with heat, embers or fire**

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.

⚠ CAUTION**Reuse on other patients and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.

NOTICE**Contact with salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.

NOTICE**Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- 1) Determine the wrist circumference.
- 2) Determine the orthosis size (see size chart).

4.2 Adaptation

The orthosis is delivered ready for trial fitting, but can be adapted by qualified personnel if required.

- 1) Cold mould the integrated reinforcing element to suit the individual user.
- 2) Undo the pull loop on the wrist strap and reposition it to suit the individual user.
- 3) Shorten the thumb closure along the embossing to suit the individual user.

4.3 Application

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.

⚠ CAUTION

Use on skin sensitive to heat

Skin irritation due to overheating

- ▶ Do not use the product in case of a known heat allergy.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- > The wrist strap is open.
- > The thumb closure is loosened.
- 1) Pull the orthosis over the thumb.
- 2) Guide the wrist strap through the guide loop (see fig. 2) and fasten it (see fig. 3).
- 3) Close the thumb closure (see fig. 4).

4.4 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

- 1) Fasten all hook-and-loop closures.
- 2) Hand wash the orthosis in warm water at **30 °C** with standard mild detergent.
- 3) Rinse well.

- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH