




ottobock.

CE

9X14, 9X18, 9X37, 9X25

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	14
FR Instructions d'utilisation (Personnel spécialisé)	26
IT Istruzioni per l'uso (Personale tecnico specializzato)	38
ES Instrucciones de uso (Personal técnico especializado)	49
PT Manual de utilização (Pessoal técnico)	61
NL Gebruiksaanwijzing (Vakmensen)	73
SV Bruksanvisning (Fackpersonal)	85
DA Brugsanvisning (Faguddannet personale)	96
NO Bruksanvisning (Fagpersonell)	107
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	118
RU Руководство по применению (Квалифицированный персонал) ...	130
JA 取扱説明書 (有資格担当者)	142

INFORMATION

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INFORMATION

Date of last update: 2021-03-03

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Function

This product was created for additional signal generation and/or signal support. The product is for use with patients who cannot generate their own muscular signals.

9X14 harness pull switch function

The harness pull switch is fastened to the stirrups between the definitive socket and shoulder harness. The harness pull switch features 2 control points that are controlled one after the other and independently. The control points are reached by pulling on the stirrup of the cable pull switch (see fig. 1).

9X18 cable pull switch function

The cable pull switch is installed in the definitive socket. The cable pull switch features 2 control points that are controlled one after the other and independently. The control points are reached by pulling on the cable of the cable pull switch (see fig. 2).

9X25 rocker switch function

The rocker switch is installed in the definitive socket. The rocker switch features 2 control points that are independently controlled. The control points are reached by pressing one side of the switch rocker (see fig. 3).

9X37 pressure switch function

The pressure switch is installed in the definitive socket. The pressure switch features 2 control points that are controlled one after another and independently. The control points are reached by pressing the switch rocker (see fig. 4).

Explanation of the control signals

If control point S1 is reached on the switch, control signal A is emitted. If control point S2 is reached, control signal B is emitted. Functions can be assigned for every control point/control signal depending on the prosthetic component used.

2.2 Combination possibilities

This product can be combined with the following Ottobock components:

- DynamicArm: 12K100N=* / 12K110N=*
- ErgoArm Hybrid plus: 12K44=*
- ErgoArm Electronic plus: 12K50=*
- AxonMaster 13E500
- System Electric Hand DMC plus: 8E38=6, 8E39=6
- System Electric Hand Digital Twin: 8E38=7, 8E39=7
- Sensor Hand Speed: 8E38=8, 8E39=8
- MyoHand VariPlus Speed: 8E38=9, 8E39=9, 8E41=9
- Transcarpal Hand DMC plus: 8E44=6
- System Electric Greifer DMC VariPlus: 8E33=9, 8E34=9
- 8E70=* bebionic hand EQD
- 8E71=* bebionic hand Short Wrist
- 8E72=* bebionic hand Flex
- 10S17 electric wrist rotator
- 13E205 MyoRotronic

3 Intended use

3.1 Indications for use

The product is intended **exclusively** for upper limb exoprosthetic fittings.

3.2 Conditions of use

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see page 24).

3.3 Indications

- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 CAUTION
Use of a damaged product
Injury due to loss of product functionality.
▶ Before use, perform a visual inspection to check if all parts of the product are undamaged.
▶ Replace the product immediately in case of damage.

 CAUTION
Non-observance of safety instructions
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety instructions and the stated precautions in this accompanying document.

⚠ CAUTION

Independent manipulation of the product

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ Opening and repairing the product or repairing damaged components is not permitted.

⚠ CAUTION

Signs of wear on the product

Injury due to faulty control or malfunction of the product.

- ▶ If a malfunction of the product is noted, it has to be replaced.
- ▶ Opening and repairing the product or repairing damaged components is not permitted. .

⚠ CAUTION

Faulty adjustment

Injury due to faulty control or malfunction of the product.

- ▶ When mounted into a provisional socket, the product may only be worn within an orthopaedic workshop or as an inpatient in a rehabilitation facility.

⚠ CAUTION

Penetration of liquids

Injury due to faulty control or malfunction of the product.

- ▶ Make sure that no liquids can penetrate into the product.
- ▶ Do not use the product for bathing prostheses.

⚠ CAUTION

Damage to the connection cable

Injury due to faulty control or malfunction of the product.

- ▶ When installing the connection cable, make sure there are no tight angles or kinks.
- ▶ Replace damaged connection cables immediately.

Additional safety instructions for cable pull switch/harness pull switch

CAUTION

Faulty adjustment

Injury due to faulty control or malfunction of the product.

- ▶ A load may only be applied to the product in the pulling direction.
- ▶ Please note that the cable travel of the product is limited. To prevent excessive strain on the product, make sure the harness is correctly fitted and accommodates the patient's movements.

CAUTION

Damage to the Bowden cable

Injury due to faulty control or malfunction of the product.

- ▶ When installing the Bowden cable, avoid sharp kinks and tight angles and check whether the cable runs smoothly.
- ▶ Replace damaged cables immediately.

4.4 Patient information

WARNING

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

WARNING

Operating vehicles and machinery

Injury due to faulty control or malfunction of the product.

- ▶ Driving vehicles of any kind or operating machinery is strictly prohibited while wearing a provisional socket and the product.

CAUTION

Penetration of dirt and humidity into the product

Injury due to faulty control or malfunction of the product.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

⚠ CAUTION

Independent user changes or modifications made to the product

- > Injury due to faulty control or malfunction of the product.
- > Destruction of the product.
- ▶ Only have any changes or modifications to the product carried out by a certified Ottobock Service Centre.

⚠ CAUTION

Mechanical loads

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product and its housing for visible damage before each use (e.g., cracking or breakage). Should the product be damaged or not function properly, please contact your prosthetist immediately.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see page 24).

NOTICE

Improper care

Damage or breakage through the use of acetone, petrol or similar solvents.

- ▶ Clean the product and the contacts with a damp cloth and mild soap only.

Additional safety instructions for cable pull switch/harness pull switch

⚠ CAUTION

Mechanical loads

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product, the cable stirrup, and the pull-cable and its housing for visible damage before each use (e.g., cracking or breakage). Should the product be damaged or not function properly, please contact your prosthetist immediately.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. Instructions for use (qualified personnel)
 - 1 pc. 9X14 harness pull switch
- or
- 1 pc. 9X18 cable pull switch
- or
- 1 pc. 9X25 rocker switch
- or
- 1 pc. 9X37 pressure switch

5.2 Accessories

The following components are not included in the scope of delivery and must be ordered separately:

- 13E50=* connection cable (only in combination with 10S17 electric wrist rotator)
- 13E97=* connection cable (only in combination with 10S17 electric wrist rotator)
- Connection cable 13E99=1200

6 Preparing the product for use

6.1 Preparing the switches

6.1.1 Required tools

- 706Z2 side cutters

6.1.2 Preparing the cable pull switch

Mounting the pull-cable to the stirrup

When mounting the cable pull switch to a definitive socket, you must connect the free end of the pull-cable to the stirrup.

1. Determine the necessary cable length.
2. Loop the pull-cable and feed it into the stirrup.
3. Place the wedge into the stirrup and tighten the pull-cable (see fig. 5).
4. Check the length of the cable and adjust if necessary.
5. After the correct cable length is set, cut the pull-cable to the correct length with an appropriate tool.

6.1.3 Adjusting the control points

Two adjustment screws are located on the back of the pressure switch. These adjustment screws can be used to determine the pressure at which the respective control point is activated.

1. Turn the right adjustment screw (see fig. 6, item 2) clockwise to increase the intensity of control point 1. Turn the adjustment screw anticlockwise to decrease the intensity of control point 1.
2. Turn the left adjustment screw (see fig. 6, item 1) clockwise to increase the intensity of control point 2. Turn the adjustment screw anticlockwise to decrease the intensity of control point 2.

6.2 Connecting the connection cable

6.2.1 Required tools

- Allen keys
- 633F11 silicone grease

6.2.2 Connect connection cable to harness pull switch

1. Unscrew the housing of the harness pull switch.
2. Carefully lift the switch unit out of the housing.
3. Grease the cable socket of the connection cable with sufficient silicone grease.
4. Connect the cable socket to the switch contact.
5. Place the switch unit into the housing.
6. Screw the housing of the harness pull switch together.

6.2.3 Connect connection cable to cable pull switch/pressure switch/rocker switch

1. Grease the cable socket of the connection cable with sufficient silicone grease.
2. Connect the cable socket to the switch contact.

7 Use

7.1 Control signals

The following events can be executed by the control signals generated by the switch:

- Controlling the system components
- Switching/pulse

Controlling the system components

The system component can be controlled by using control signals A and B. The following movement patterns can be generated in combination with the various system components:

- Pronation/supination
- Flexion/extension
- Opening/closing

Switching/pulse

Control signal B can be used to switch between the various system components (elbow, rotation, gripping components). Furthermore, the elbow joint lock can be activated by control signal B with certain system components.

8 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1-N Derma Clean) when needed.
Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

9 Legal information

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

9.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

10 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

9X14 harness pull switch

General	
Article number	9X14
Product service life	5 years
Switching voltage	12 V
Switching current	0.5 A

9X18 cable pull switch

General	
Article number	9X18
Product service life	5 years
Switching voltage	12 V
Switching current	0.5 A

9X25 rocker switch

General	
Article number	9X25
Product service life	5 years
Switching voltage	12 V
Switching current	0.5 A

9X37 pressure switch

General	
Article number	9X37
Product service life	5 years
Switching voltage	12 V
Switching current	0.5 A

11 Appendices

11.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Lot number (PPPP YYYY WW)
PPPP – plant
YYYY – year of manufacture
WW – week of manufacture

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