


ottobock.

21Y12, 21Y13

DE	Gebrauchsanweisung	6
EN	Instructions for use	10
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	22
PT	Manual de utilização	27
NL	Gebruiksaanwijzing	31
SV	Bruksanvisning	35
PL	Instrukcja użytkowania	39
HU	Használati utasítás	43
EL	Οδηγίες χρήσης	47

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English

Last update: 2020-04-07

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

Explanation of Symbols

NOTICE Warnings regarding possible technical damage.

1 Application

The **21Y12 Screw Valve Kit** is to be used **solely** for installation in **thermo-plastic** transfemoral [AK] sockets.

The **21Y13 Screw Valve Kit** is to be used **solely** for installation in **laminated** transfemoral [AK] sockets.

NOTICE

Risk of loss of function due to incorrect environmental conditions. Using the prosthetic components in incorrect environmental conditions can lead to malfunctions of the valve up to loss of functionality.

Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids. Penetration of liquids, dust, sand, etc. as well as the use at places with excessive ambient moisture must be avoided.

NOTICE

Damage due to using powder. Any sort of powder on the inside of the socket or on the skin can reduce or disrupt valve function during use. Avoid using powder in areas close to the valve.

2 Components / Service Parts

2.1 Screw Valve Kit 21Y12 (fig. 1):

- (1) Lower screw ring
- (2) Sealing ring 21Y103
- (3) Upper screw ring
- (4) Upper valve piece 21Y145
- (5) Vacuum forming template 21Y157
- (6) Screw 503S3
- (7) Lamination dummy 21Y160
- (8) Countersunk head screw 501T19=M5×18
- (9) O-ring 627F13=25×2
- (10) Pin wrench

2.2 Screw valve kit 21Y13 (fig. 2):

- (1) Small screw ring
- (2) Upper valve piece 21Y145
- (3) Lamination dummy
- (4) Screw 503S3
- (5) O-ring 627F13=25×2

3 Installation instructions

The screw valves can be installed medially or laterally at the distal end.

The template (dummy) provided is inserted before vacuum forming or laminating, making it possible to accurately install the valves into various socket materials.

3.1 Screw valve kit 21Y12

Saturate a 99B25 Nylon Sock with N-Wax=1 and pull it on the plaster positive. Attach the 21Y157 Vacuum forming template to the plaster positive using the 503S3 Screw provided (fig. 3).

After vacuum forming the test socket or the ThermoLyn definitive inner socket (fig. 4), smooth the vacuum forming template by grinding it until it is flush with the ThermoLyn material (fig. 5).

Attach the lamination dummy to the vacuum forming template using the 501T19=M5×18 Countersunk Head Screw provided (fig. 6).

For frame reinforcement, the 5Z1 Lamination Kit with the aid of the Technical Information Booklet 2.3.5 or the 5Z4 Lamination Kit for CAT-CAM can be used.

After lamination, evenly grind off the excess lamination around the lamination dummy and remove the countersunk head screw.

Remove the socket frame and inner socket from the plaster positive.

Disassemble the screw valve and apply **617H46 Bonding Agent** to the contact surface of the lower screw ring. Position the lower screw ring. Slightly lubricate the sealing ring of the upper screw ring with 633F11 Silicone Grease. Tighten upper screw ring with attached sealing ring (fig. 8). For tightening the lower and upper screw ring the Pin Wrench can be used.

NOTICE

Risk of damage to the valve. The valve is not intended for the fixation of the inner socket and frame socket. It has no holding function. Only attach the valve to the test socket or to the flexible inner socket.

Grind the laminate of the socket frame in the valve area only until the groove for the cosmetic sock on the upper screw ring lies free with the inner socket inserted.

Put the socket frame and inner socket together and screw in the upper valve piece.

If the valve has become dirty, remove the 21Y145 Upper valve piece using a pin wrench. Clean the individual components and reassemble.

3.2 Screw valve set 21Y13

Flatten the plaster positive in the valve area, apply talc and pull on a moistened 99B81 PVA Bag. Reinforce the PVA Bag in the valve area with 672B4 Polyethylene Adhesive Tape.

Remove the upper valve piece from the screw valve. Lubricate the inside of the small screw ring with 633F11 Silicone Grease and mount the lamina-

tion dummy. Lubricate the distal surface on the plaster positive in the area of the reinforced PVA bag with silicone grease (fig. 9).

Attach the small screw ring with the lamination dummy to the plaster positive using the 503S3 Screw provided (fig. 10). Seal the screw head area with 636K8 Plasta Tape.

Cut 616G6 Dacron Felt as reinforcement material and pull it on the plaster positive. Weld the edges to each other and remove the excess material covering the screw valve (fig.11).

Pull on a double layer of 623T3 Perlon Stockinette. Tie off the first layer in the screw valve's groove (fig.12). Reinforce the area around the socket brim with carbon. Next, pull on another double layer of Perlon stockinette and tie it off in the screw valve's groove.

Pull on a 99B81 PVA Bag and laminate.

Before the lamination resin hardens, apply polyethylene adhesive tape to the small screw ring above the groove in order to squeeze out any excessive lamination resin (fig.13).

After the lamination resin has hardened, remove the socket from the plaster positive and remove the lamination dummy.

Screw in the upper valve piece (fig.14).

If the valve has become dirty, remove the 21Y145 Upper valve piece using a pin wrench. Clean the individual components and reassemble.

4 Cleaning

The valve should be cleaned regularly with clear, warm water.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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