


4R136

DE	Gebrauchsanweisung	3
EN	Instructions for use	6
FR	Instructions d'utilisation	9
IT	Istruzioni per l'uso	12
ES	Instrucciones de uso	15
PT	Manual de utilização.....	18
NL	Gebruiksaanwijzing.....	22
SV	Bruksanvisning.....	25
DA	Brugsanvisning.....	31
NO	Bruksanvisning.....	31
FI	Käyttöohje.....	34
PL	Instrukcja użytkowania.....	37
HU	Használati utasítás	40
CS	Návod k použití.....	43
RO	Instrucțiuni de utilizare.....	46
HR	Upute za uporabu	49
BG	Инструкция за употреба	52
TR	Kullanma talimatı.....	55
EL	Οδηγίες χρήσης	58
RU	Руководство по применению	61
JA	取扱説明書.....	65
ZH	使用说明书.....	68

INFORMATION

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INFORMATION

Last update: 2020-04-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

Explanation of Symbols

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

The scope of delivery is shown on the front page.

1 Description

1.1 Indications for use

The Ottobock 4R136 exhaust valve is intended exclusively for use in prosthetic fittings of the lower limbs.

1.2 Field of application

In combination with a sealing sleeve, the 4R136 exhaust valve serves to create a vacuum in a prosthetic socket. Thanks to this technology, shearing and friction forces are effectively reduced.

1.3 Safety instructions

Please make your patients aware of the following safety instructions:

NOTICE

Risk of corrosion. Prosthetic components must not be exposed to substances which could cause metal parts to corrode, such as fresh water, salt water and acids.

Using this medical device in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

Failure to observe the installation instructions can result in complete loss of vacuum and therefore significantly compromise the adhesion.

1.4 Function

The 4R136 exhaust valve ejects excessive air out of the socket during the stance phase and closes during the swing phase to prevent further inflow of air. Thanks to this system, shearing and friction forces are effectively reduced.

2 Components

- V4 valve, 1 pc.
- Socket connector, angled, 1 pc.
- Socket connector, straight, 1 pc.
- Tube, 1 pc.
- Washer, 2 pcs.

3 Handling

3.1 Installation preparation

The following tools are needed in addition for installing the socket connector:

- Pliers

- 5, 6, 10 mm drill bit
- 2-component epoxy adhesive (5 min. curing time) or sealing resin-talcum mixture
- Insulating or polyethylene tape
- 1/4 UNF tap

3.2 Completion

Laminate the socket in the usual way. Drill a 10 mm bore hole into the posterior side of the socket, as close to the distal end as possible. Remove any burrs caused by drilling. Secure from inside with polyethylene adhesive tape. Fill the hole with lamination resin to achieve air tightness (Fig. 1). Drill a 5 mm hole in the centre of the sealing resin fill and cut the thread with a 1/4 UNF tap, or drill a 6 mm hole in the centre of the sealing resin fill and cut the thread in the hole using the connector. Verify that the thread of the elbow connector is not longer than the laminate's wall thickness. If need be, use a washer (Fig. 2).

Mix the 2-component epoxy adhesive or the sealing resin-talcum mixture. Apply a thin layer of the mass to the thread of the connector (Fig. 3). When applying the adhesive, be careful not to let any adhesive get into the open end of the connector.

Screw the connector fully into the prepared hole using the pliers (Fig. 4). Slide the tube onto the connector.

Install the exhaust valve on the end of the hose. The arrow on the valve shows the direction of air flow (Fig. 5).

4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

4.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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