
17B26*, 17B3*, 17B43*, 17B47*

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	20
ES	Instrucciones de uso	26
PT	Manual de utilização	33
NL	Gebruiksaanwijzing	39
SV	Bruksanvisning	45
DA	Brugsanvisning	50
NO	Bruksanvisning	56
PL	Instrukcja użytkowania	62
HU	Használati utasítás	68
CS	Návod k použití	74
TR	Kullanma talimatı	80
EL	Οδηγίες χρήσης	85
RU	Руководство по применению	92

1

INFORMATION

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INFORMATION

Date of last update: 2021-05-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17B3*, 17B26*, 17B43* and 17B47* orthotic knee joints.

2 Product description

2.1 Available sizes

Article number	Medial and lateral joint	System width	Posterior placement	Material
17B26=- L/R16	Outside straight, inside bottom contoured	16 mm	16 mm	Stainless steel
17B26=- L/R20		20 mm		
17B47=- 16	Outside and inside straight	16 mm		
17B47=- 20		20 mm		
17B3=L/ R16	Outside straight, inside bottom contoured	16 mm	22 mm	
17B3=L/ R20		20 mm		
17B43=- 16	Outside and inside straight	16 mm	16 mm	
17B43=- 20		20 mm		

2.2 Components/design

Scope of delivery (see fig. 1)			
Item	Quantity [Piece (s)]	Designation	Article number
1	1	Upper joint section	
2	4	Screws	501T7*
3	1	Joint screw	501A6=4x5xM5
4	1	Ball bearing	509K11=5x16x4
5	1	Lower joint section	

Not included in scope of delivery			
Item	Quantity [Piece (s)]	Designation	Article number
(see fig. 1), item 6	1	Side bar	17B38*, 17B39*, 17B40*, 17B41*, 17B84*
not illustrated	1	Lamination bar	17Y128*, 17Y129*

3 Intended use

3.1 Indications for use

The orthotic joint is intended **exclusively** for use as a pair for orthotic fittings of the lower limbs on one patient.

3.2 Indications

- Partial or total paralysis of the leg muscles
- Orthopaedic diseases of the lower limbs

Indications must be determined by the physician.

3.3 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

3.4 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

3.5 Combination possibilities

The 17B3*, 17B26*, 17B43* and 17B47* orthotic knee joints can be combined with the 17Y32 or 17Y35 system extension assist.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

 **CAUTION**

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

 **CAUTION**

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **CAUTION**

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.

- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below **-10 °C (14 °F)** or above **+60 °C (140 °F)** (e.g. sauna, excessive sunlight, drying on a radiator).

INFORMATION

Increased wear and tear

The orthotic knee joints are intended only for use in combination with free-motion or foot lifting orthotic ankle joints. Wear occurs more rapidly in combinations involving orthotic ankle joints with a dorsal stop.

5 Preparing the product for use

INFORMATION

Parallel alignment of the orthotic joints

Use the 743R6 alignment fixture to position the orthotic joints on the plaster positive.

Gluing in bars

Required materials: 636W28* special adhesive, degreasing cleaner

The bars must be glued in the insertion zones to ensure stability.

- 1) Clean the connecting surfaces with a degreasing cleaner.
- 2) Apply the special adhesive to the connecting surfaces.
- 3) Insert the bars.
- 4) Secure with the enclosed screws.

- 5) Allow to cure for at least **4 hours**.

INFORMATION: Final bonding strength will be reached after 16 hours.

Improving corrosion protection

Polish or sinter the product for improved corrosion protection. Ottobock recommends the 618T40* sintering powder for surface coating.

- ▶ Do not heat the material longer than **5 minutes** at max. **150 °C** for surface coating.

Using protective devices

- ▶ Use joint protectors to protect the product against soiling.

6 Cleaning

Promptly clean the product after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e. g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Spare parts are listed in the section "Components/design".

Only use 633F7 special lubricant.

Replacing the ball bearing

NOTICE

Canting the ball bearing

Mechanical damage

- ▶ Do not cant or damage the ball bearing while inserting it into the orthotic joint.

- 1) Loosen the screws and disassemble the orthotic joint.
- 2) Remove the ball bearing from the orthotic joint and replace it with a new one.
- 3) Assemble the orthotic joint and tighten the screws.

- 4) For the final assembly, secure the screw connection with Loctite 241.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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