




17KF10, 17KL40, 17KL20

DE	Gebrauchsanweisung	3
EN	Instructions for use	9
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	19
ES	Instrucciones de uso	24
PT	Manual de utilização	29
NL	Gebruiksaanwijzing	35
SV	Bruksanvisning	40
DA	Brugsanvisning	45
NO	Bruksanvisning	50
PL	Instrukcja użytkowania	55
HU	Használati utasítás	61
TR	Kullanma talimatı.....	66
EL	Οδηγίες χρήσης	71
RU	Руководство по применению	77

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Explanation of symbols

English

⚠ WARNING Warnings regarding possible risks of severe accident or injury.

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Last update: 2016-10-14

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place

INFORMATION

The Knee Joints can also be used in conjunction with the Ottobock Modular System Ankle Joints. For this purpose special adapters are at your disposal (17M60/61) allowing the connection of knee joints with the modular system.

1 Application

The 17KF10; 17KL40 and 17KL20 Knee Joints are to be used **exclusively** for the orthotic fitting of the lower extremities.

2 Indication

Partial paralysis or complete paralysis of leg muscles. Indication must be diagnosed by a physician.

3 Processing

INFORMATION

Please note that the Ottobock Joint Bar Systems, due to different alignment parameters such as the created forefoot load or flexion, are subjected to increased load.

Please observe the following processing instructions:

3.1 Bending the joint bars

The bending radius should be at least 20mm.

While bending, ensure that no bending load is applied on joint heads so as to prevent the knee joints from premature wear or damage.

To enable you to identify this area we have laser etched dotted lines on both upper and lower joint bars indicating not to bend within this area.

3.2 Heating the joint bars

All joint bars are heat treated. The aluminum joint bars should not be exposed to temperatures exceeding 150 °C (318 °F), the steel joint bars should not be exposed to temperatures exceeding 300 °C (572 °F).

3.3 Drill holes

It is permitted to drill holes into the joint bars up to dia. 5mm and the hole must be positioned centrally of the joint bar. **Never** repair or fill the drill holes.

3.4 Fastenings

The locking function of the joint bar system must never be obstructed by closures of the attachment straps or knee pads.

3.5 Setting-up the joints

Use 633F7 Special Lubricant to grease the metal moving parts. All screw heads must be secured with 636K13 Loctite.

4 Adjustments and application

4.1 Documentation

To ensure traceability in terms of the product liability as provided by law, the article number and batch number must be entered into patient's medical chart.

4.2 Activity level

The Knee Joints were designed exclusively as an anatomical support/walking aid. They are not suited for excessive strain such as sports activities, gardening or other similar activities (lifting and carrying heavy objects, among other things) which would involve heavy physical strain.

4.3 Locking knee joints

The knee joints with locks such as ring, bale and lever are delivered ready for use. No further adjustments are necessary on the locking mechanism or knee joint.

To ensure a safe function, the bale arms **should not** be rigidly connected to each other in bale lock knee joints.

INFORMATION

Before handing over an orthosis with bale lock knee joints, check that the locks are fully engaged and instruct the patient how to use the orthosis.

5 Maintenance and service

Use only 633F7 Special Lubricant as a lubricant.

The knee joint mechanism should be protected from dirt by a suitable protective device (e.g. joint protectors) for proper functioning of the orthosis.

Designation	Reference number	17KL40=*	17KL20=13	17KL20=16, 17KL20=20	17KF10=16	17KL19=2
Joint screw	501S32=M5X9X12					
	501S32=M5X9X10	•			•	
	501S32=M5X9X18					•
	501S32=M5X9X16		•			
	501S32=M6X11X16			•		
Red washer	507U5=11.4X7.4X1.5	•				
	507U5=10.6X7.4X1.5				•	•
	507U5=19X14		•			
	507U5=25X20			•		
Bushing	17Y17=10X4.5X7	•			•	
Bearing nut	17Y93=M5X7X6.4	•			•	
	17Y93=M5X9X8		•			
	17Y93=M6X9X10			•		
Ball cage	501D1	•				
Lock ring (optional ring lock)	17Y13=9	•				
Lock spring (optional ring lock)	513D83=0.5X3.7X22	•				
Auto screw	501S32=M3X12	•				
Lever	17Y204					•
	17Y20=1		•			
	17Y20=2			•		
Cam bushing	17Y17=7X11.5					•
Nylon washer	507U57=12X14					•
Lock ring	507S96=4		•	•		
Cap screw	507U57=16X5.4		•			
	507U57=22X7.6			•		
Ball pin	4A101=5X12.8		•	•		
	170D5=120			•		
Spring	513D83=1X4.1X22.8		•			
	513D83=1X4.1X26.2			•		

INFORMATION

The knee joints must be checked for wear and tear every six months. When the knee joints are subjected to extreme environmental influences (e.g. high air humidity, industrial dust, patient's incontinence etc.), the inspection must be carried out every three months at least.

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7 Patient information

Please make sure that the patient is informed about the following items:

CAUTION

Improper use by the patient

Injuries due to falling

- If your orthosis is provided with knee locks, be certain that both knee locks are activated before standing up!
- Do not make any adjustments or modifications on the knee or ankle joints.

Regular inspections

Check your orthosis/joint regularly for damages of adjustable straps/buckles (such as tears, loosening). If you ascertain any problems, please contact your orthopedic technician.

If the knee locks fail or any problems with the orthosis (joint bar) occur, please contact your orthopedic technician immediately.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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