



# ottobock.

## 12V10

|           |  |    |
|-----------|--|----|
| <b>DE</b> | Gebrauchsanweisung (Fachpersonal) .....                      | 3  |
| <b>EN</b> | Instructions for use (Qualified Personnel) .....             | 7  |
| <b>FR</b> | Instructions d'utilisation (Personnel spécialisé) .....      | 12 |
| <b>IT</b> | Istruzioni per l'uso (Personale tecnico specializzato) ..... | 16 |
| <b>ES</b> | Instrucciones de uso (Personal técnico especializado) .....  | 21 |
| <b>PT</b> | Manual de utilização (Pessoal técnico) .....                 | 26 |
| <b>NL</b> | Gebruiksaanwijzing (Vakmensen) .....                         | 30 |
| <b>SV</b> | Bruksanvisning (Fackpersonal) .....                          | 35 |
| <b>DA</b> | Brugsanvisning (Faguddannet personale) .....                 | 39 |
| <b>NO</b> | Bruksanvisning (Fagpersonell) .....                          | 43 |
| <b>PL</b> | Instrukcja użytkowania (Personel fachowy) .....              | 47 |
| <b>HU</b> | Használati utasítás (szakszemélyzet) .....                   | 52 |
| <b>CS</b> | Návod k použití (Odborný personál) .....                     | 56 |
| <b>RO</b> | Instrucțiuni de utilizare (Personal de specialitate) .....   | 61 |
| <b>HR</b> | Upute za uporabu (Stručno osoblje) .....                     | 66 |
| <b>TR</b> | Kullanma talimatı (Uzman personel) .....                     | 69 |
| <b>EL</b> | Οδηγίες χρήσης (Τεχνικό προσωπικό) .....                     | 74 |
| <b>RU</b> | Руководство по применению (Квалифицированный персонал) ..... | 79 |
| <b>JA</b> | 取扱説明書 (有資格担当者) .....   | 84 |

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English

Date of the last update: 2020-11-27

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

## 1 Intended use

### 1.1 Medical purpose

The 12V10 Tube Valve for the Suction Socket is to be used **exclusively** for assembly in standard and suction sockets of upper limb prostheses.

### 1.2 Application

The Tube Valve for the 12V10 Suction Socket can be used for unilateral or bilateral amputees.

### **1.3 Application conditions**

The 12V10 Tube Valve for the Suction Socket can be used in standard sockets.

The 12V10 Tube Valve for the Suction Socket can be used in suction sockets in combination with the 13E202=\* Suction Socket Electrode.

### **1.4 Prosthetist qualifications**

The fitting of a patient with the 12V10 Tube Valve for the Suction Socket must be carried out exclusively by orthopedic technicians.

## **2 Safety instructions**

Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

- Before you process the 617H21=\* Orthocryl Sealing Resin and 617P37=\* Hardening Powder, please first read the relevant safety data sheet.
- Use the Patient Information section of this pamphlet (Section 3) to inform your patient about the function and use of the tube valve.

## **3 Patient information**

- Regularly clean the socket, the connection tube and the tube valve using a cloth slightly moistened with a soap solution. Rinse thoroughly.
- To degrease the tube valve and connection tube, use 640F18 Donning Spray for Silicone Liners.
- Make sure that no foreign particles get into the tube valve.
- The tube valve should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- In case of noticeable leakage, the tube valve should be removed from the connection tube. Rinse the tube valve and connection tube with water and dry it. Reinsert the tube valve.
- Put the tube valve into the connection tube just deep enough to be in slight contact with the skin, then close the valve

## 4 Scope of delivery and spare parts

- 1 12V10 Tube valve for Suction Socket
- 1 99B13=21 Connection Tube
- 1 Instructions for use (qualified personnel)

## 5 Function

The 12V10 Tube Valve for the Suction Socket is a one-way valve, which allows air to escape from the socket but prevents air from entering. Due to the movement of the residual limb in the inner socket, air is constantly being pumped out through the tube valve. The resulting, minor vacuum helps create optimal linkage with the patient's residual limb.

## 6 Application

Pull the residual limb into the socket using the pull-in sock or Gliding Sock. Once the pull-in sock has been pulled completely out of the connection tube, reinsert and close the tube valve.

The function of the Suction Socket or re-equipped Standard Socket in myoelectric prostheses is only guaranteed when used in combination with the 12V10 Tube Valve and 13E202=\* Suction Socket Electrode(s).

## 7 Fabrication instructions

**Note:** For clear presentation the outer and inner socket have been shown cut open in the following sections.

Finish the Inner and Outer Socket according to the Instructions for Use 647G334 for the 13E202=\* Suction Socket Electrode. Insert the 99B13=21 connection tube into the inner socket and create an **airtight** seal using 636K11 Cyamet Quick Drying Adhesive.

**Note:** The connection tube must be flush with the bottom of the socket to prevent chafing of the residual limb (Arrow Fig. 1).

Mount the inner and outer socket. The connection tube must be sanded flush with the exterior side of the outer socket (Fig. 2).

### 7.1 Cutting the Tube Valve to length

The tube valve must be flush inside the connection tube to avoid chafing of the residual limb.

- 1.) Place the tube valve into the connection tube and close it. It must be flush with the outer side of the socket (Fig. 3)
- 2.) Mark the tube valve at the bottom of the socket (Arrow Fig. 3).

**Note:** When cutting off the tube valve, take care that at least 20 mm are left up to the valve head (Fig. 4).

- 3.) Cut off the tube valve as marked, trim off any burrs and sand smooth. Remove sticker.

## **7.2 Closing the Tube Valve**

The open end of the valve must be sealed with Orthocryl Sealing Resin to prevent the residual limb from being drawn in:

- 1.) Cut off an approx. 4 cm long strip of 627B2=19 Monofilament Adhesive Tape.
- 2.) Mix the 617H21=\* Orthocryl Sealing Resin and 617P37=\* Hardening Powder and fill approx. 10 ml into the tube valve you previously cut to length.
- 3.) Close the tube valve with Monofilament tape (Fig. 5).
- 4.) Tilt the tube valve and place the taped end on a flat surface. Hold it in this position until the Orthocryl-Sealing Resin has cured. Be careful that no Orthocryl-Sealing Resin leaks out! (Fig. 6).
- 5.) Remove the adhesive tape and any Orthocryl-Sealing Resin residuals.

**Note:** Figure 7 shows a sectional view of the suction socket with tube valve.

## **8 Technical information**

Service life: 5 years

## **9 Symbols used**

Medical device

## **10 Liability**

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manu-

facturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

## **11 Trademarks**

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## **12 CE conformity**

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European erganzen!European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

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