
ottobock.

1D10

DE	Gebrauchsanweisung	3
EN	Instructions for use	6
FR	Instructions d'utilisation	9
IT	Istruzioni per l'uso	12
ES	Instrucciones de uso	15
PT	Manual de utilização.....	18
NL	Gebruiksaanwijzing.....	21
SV	Bruksanvisning	24
DA	Brugsanvisning.....	27
PL	Instrukcja użytkowania	30
HU	Használati utasítás	33
CS	Návod k použití.....	36
RO	Instrucțiuni de utilizare.....	39
HR	Upute za uporabu	42
TR	Kullanma talimatı.....	45
EL	Οδηγίες χρήσης	48
RU	Руководство по применению	51
JA	取扱説明書.....	54
ZH	使用说明书.....	57

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Date of the last update: 2014-12-08

- Please read this document carefully.
- Follow the safety instructions.

1 Description

1.1 Intended use

The Ottobock 1D10 Dynamic prosthetic foot is to be used **exclusively** for lower limb prosthetic fittings.

1.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

Size 22 – 30

Recommended for amputees with **mobility levels 1 and 2** (indoor walkers and restricted outdoor walkers).

Permitted for **patients weighing up to 150 kg (330 lbs)**.

Note: The Ottobock MOBIS Mobility System is NOT intended to be used in the USA as a guide to obtaining reimbursement for prosthetic components.

2 Safety instructions

Attention!

Improper use or use that does not correspond to the intended use may result in wear of the components of the prosthetic foot. To avoid risk to the patient, the foot must not be used after a malfunction has taken place. A noticeable malfunction may include decreased forefoot resistance.

Measure

Have the prosthetic foot examined at a qualified prosthetic facility where it may be replaced or repaired depend on applicable warranty terms.

Avoid exposing prosthetic components to elements that can:

- cause corrosion on metal parts (such as freshwater, salt water and acids) or

- have an abrasive effect (such as sand).

Using this medical device under such conditions will render all claims against Otto Bock HealthCare null and void.

The material is not UV or hydrolysis resistant.

The surface of the foot must not be damaged by grinding as this leads to premature wear.

Please forward this information to your patients.

3 Description and functional wualities

A specially designed rollover element together with a sophisticated polyurethane foam provide the main functions of this foot. The 1D10 is designed for use with modular prostheses and is delivered with an integrated titanium modular adapter.

The rollover element and foam technology offer the following functional advantages:

- Comfortable heel strike
- Elastic return of the forefoot thanks to sophisticated foam technology
- Rollover smoother than with SACH foot due to flixible keel design

4 Handling

4.1 Alignment recommendations

4.1.1 Basic alignment recommendation

Foot size	Mid-foot anterior to alignment line	Heel height
22-30 cm	30 mm	10 mm +/- 5 mm

4.1.2 Static alignment recommendations

Align the prosthesis using the LASAR Posture, if available.

Follow the Ottobock alignment recommendations when aligning lower limb prostheses.

4.1.3 Dynamic alignment recommendations

Align the prosthesis in the coronal (ML) and sagittal (AP) plane through sliding or by changing the angle to achieve proper heel contact, smooth rollover

and optimal forefoot weight transfer onto the contralateral side. With trans-tibial amputees, observe the physiological knee flexion during stance phase.

4.2 Maintenance instructions

Note:

As a basic principle, all Ottobock modular foot components are subjected to tests involving two million load cycles. Depending on the amputee's activity this corresponds to a service life of two to three years.

We recommend carrying out regular safety checks once a year.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH